




OWNER'S MANUAL OF

MODEL ○ 800○800D○80-1○80-2



Maintenance

- 1、 Operating and maintaining correctly can prolong the life of centrifuge.
- 2、 Cut off the power and clean the centrifuge if you don't use the centrifuge for a long time.

AC 220V only, 1 year warranty.

Safety Information

- It **must** be connected to an electrical outlet that is properly grounded.
- The power cord must be connected to an electrical outlet that is neat the product and easily accessible.
- Refer service or repairs, other than those described in the operating instructions, to a professional service person.

Content

Instruction 1

- Usage & Intruduction
- Packing list & Specifications
- Features & Function

Operation Instruction 2

- Operation explanation
- Operation steps

Trouble shooting guide 3

Maintenance 4



➤ Usage and Introduction

Centrifuge is possessed of these special features following: It is beautiful in appearance. Large in capacity but small in bulk and has a complete function and a stable performance, you can chose the speed that you need at the setrange, and it can keep balance automaticly, in addition, it has a little rate of increase in temperature and can be used efficiently and widely. It is a dependable instrument for qualitative analysis of serum, plasma and immunity in hospitals or chemical and biochemical laboratories.

➤ Packing list

Make sure you have the following items:

Name	Centrifuge	User's manual	Package
Number	1	1	1

If items are broken or damaged, contact your point of purchase.

➤ Specifications

Model:	800	800D	80-1	80-2
Max Speed:	4000rpm			
Max RCF:	1,430(×g)			1,790(×g)
Rotors (included):	Six-place×20ml (fixed-angle)			Twelve-place×20ml
Motor Power:	40W	25 W		40 W
Time Range:	---	0~60min	0~30min	0~60min
Dimensions:	250×230×200	250dia×220h	275×245×235	321×298×239
Shpg wt:	6 kg		7 kg	12 kg
Power:	AC 220V±10% 50Hz±2%			

➤ Features & Function

- ▶ Compact size saves room on your benchtop
- ▶ Rubber feet keep unit stable
- ▶ Easy to use

Operation explanation

- For your safety, it **must** be connected to an electrical outlet that is properly grounded.
- Please read this user's manual carefully for first use..
- Operating temperature range:0℃~40℃,Humidity: below 80%.
- Cut off power for long time no use.

Operation steps

- Connect the centrifuge power cord to the centrifuge. Then plug the centrifuge power cord into an electrical outlet.
- Turn on the power.
- Adjust the speed and time.