Rev. 05

September 2015

PROBES - CLEANING, DISINFECTION, STERILIZATION
SONDE - PULIZIA, DISINFEZIONE, STERILIZZAZIONE
SONDES - NETTOYAGE, DÉSINFECTION, STÉRILISATION
SONDEN - REINIGUNG, DESINFEKTION, STERILISATION
SONDAS - LIMPIEZA, DESINFECCIÓN, ESTERILIZACIÓN
SONDAS - LIMPEZA, DESINFECÇÃO, ESTERILIZAÇÃO

141015300

Table of Contents

1	Agents and Systems	1-i
1	Agenti e sistemi	1-iii
1	Agents et Systèmes	1-v
1	Mittel und Systeme	1-vii
1	Agentes y Sistemas	1-ix
1	Agentes e sistemas.....	1-xi
	Agents/Agenti/Agents/Mittel/Agentes/ Agentes	1-1
	Agents/Agenti/Agents/Mittel/Agentes/ Agentes	1-3
	System/Sistemi/Systèmes/Systeme/Sistemas/Sistemas.....	1-5
	Agents/Agenti/Agents/Mittel/Agentes/ Agentes	1-7
	System/Sistemi/Systèmes/Systeme/Sistemas/Sistemas.....	1-8
	Needle Guides/Guide d'ago/Guides d'aiguille/Biopsie Kits/Guias de aguja/Guías de agulhas de biópsia	1-9
	Manufacturers of Recommended Agents and Systems	1-11
	Produttori degli agenti e dei sistemi raccomandati.....	1-23
	Fabricants des agents et des systèmes recommandés	1-35
	Hersteller der empfohlenen Mittel und Systeme	1-47
	Fabricantes de los agentes y los sistemas recomendados	1-59
	Fabricantes de agentes e de sistemas recomendados.....	1-71

Agents and Systems

The following tables list the **MyLab** probes and the recommended cleaning, disinfection and sterilization agents and systems.

The symbol “√√” means that the agent/system has been tested as both **chemical compatible and biological efficacy** with the probe.

The symbol “√” means that the agent/system has been tested as **chemical compatible** with the probe.

“N” means that the agent/system has been tested as **NOT compatible** with the probe.

The symbol “⊗” means **not tested**.

Cleaning, disinfection and sterilization are allowed up to the immersible part limit.

PC

Refer to the “Probes and Consumables” manual for detailed information on how to clean, disinfect and sterilize probes and needle guides.

CAUTION

Do not use not tested agents/ systems either to clean, or to disinfect or to sterilize the probe. The probe could be damaged if such agents/systems are used.

Note

Any damage caused by the use of not recommended agents/ systems is not covered by the warranty.

WARNING

The listed disinfection/sterilization agents and systems are recommended because of chemical compatibility with the probe materials.

The tables in the next pages indicate which agents/systems have been positively tested for both biological effectiveness and chemical compatibility and which for chemical compatibility only. In the latter case follow the guidelines and recommendations of the manufacturer for the biological effectiveness of agents/systems.

Use of solutions other than those referenced is not recommended. They may damage the probe housing or the acoustic lens. Esaote takes no responsibility for damage caused by using non-approved products.

Follow the instructions provided by the manufacturer of the agent for proper use. Observe specifically soak times and dilution rates.

Overexposure to the disinfection fluid can damage the probe.

Personnel should adopt all necessary protective measures during the probe cleaning, disinfection and sterilization processes (for example gloves, protective glasses).

Never attempt to clean or disinfect the probes while they are connected to the system.

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System/Sistemi/Systèmes/Systeme/Sistemas/Sistemas

Follow the guidelines and recommendations of the manufacturer for the correct use of the system.

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- a. The frequent use of Ultraviolet C (UVC) may cause a yellowing of the cable and of the plastic parts, though not compromising the safety and reliability of the probe.

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All compatible sterilization methods for surgical instruments can be applied to the needle guide kits ABS421, ABS523, ABS424, ABS33A, DBSE12X, WBAC54X and WBSL33X.

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Manufacturers of Recommended Agents and Systems

Table 1: Advanced Sterilization Products

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Cidex®	No	Liquid	Glutaraldehyde	USA	Advanced Sterilization Products (www.sterrad.com/about_asp/locations/index.asp)
Cidex Activated Dialdehyde Solution	No	Liquid	Alkalyne glutaraldehyde solution		
Cidex® OPA	Yes	Liquid	Ortho-phthalaldehyde (OPA)		
Cidexplus®	No	Liquid	Glutaraldehyde		
CidezimeXTRA	No	Enzymes	Cellulase <5% Lipase <5% Amylase <5% Protease <5% Sodium tetraborate decahydrate <8.5%		
DisOPA	No	Liquid	Ortho-phthalaldehyde (OPA)		
Enzol/Cydezime	No	Liquid	Enzymatic detergent		
NuCidex	No	Liquid	Peracetic acid		
Sterrad® NX	No	Liquid	Hydrogen Peroxyde		
Sterrad® 100NX	No	Liquid	Hydrogen Peroxyde		
Sterrad® 100S	No	Liquid	Hydrogen Peroxyde		

Table 2: Alkapharm UK ltd

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Alkazyme	No	Powder	IPA, Quaternary Ammonium Chloride	UK	Alkapharm UK. Ltd (www.alkapharm.co.uk)

Table 3: Amity International

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Virusolve	No	Wipes	Dodecylamine	UK, USA	Amity International (www.amityinternational.com)

Table 4: Biosteril

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Peraction	No	Liquid	5% Peracetic acid	Italy	Biosteril (www.biosteril.it)

Table 5: BODE Chemie

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Korsolex extra	No	Liquid	Glutaraldehyde	Germany	BODE Chemie GmbH (www.bode-chemie.de)
Mikrobac	No	Tissues	Quaternary ammonium		

Table 6: Borer Chemie

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Deconnex 3-Zyme	No	Liquid	Enzymes	Switzerland	Borer Chemie AG (www.borer.ch)

Table 7: Day-Impex

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Virkon	No	Liquid	Potassium peroximonosulphate	UK	Day-Impex Ltd (www.day-impex.co.uk)

Table 8: Dr.Schumacher

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Cleanisept Wipes	No	Wipes	0.25% didecydimethyl ammonium chloride 0.25% alkyl dimethyl benzyl ammonium chloride 0.25% alkyl dimethyl ethylbenzyl ammonium chloride	Germany	Dr.Schumacher Gmbh (www.schumacher-online.com)
Septi Wipes	No	Wipes	Quaternary ammonium		

Table 9: Dr.Weigert

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Mediclean	No	Liquid	Enzymes	Germany	Dr.Weigert Gmbh (www.drweigert.com)

Table 10: Ebiox

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Trionic D	No	Spray, Wipes	0.5% alkyl dimethyl benzyl ammonium chloride	UK	Ebiox (www.ebiox.co.uk)

Table 11: Ecolab

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Asepti-Wipes II	No	Wipes	14.85% Isopropyl Alcohol 0.125% Alkyl Dimethyl Benzyl Ammonium Chloride 0.125% Alkyl Dimethyl Ethylbenzyl Ammonium Chloride	USA	Ecolab (www.ecolab.com)
Sekusept Aktiv	No	Powder	Peracetic acid		
Ster-Bac blu	No	Liquid	Quaternary ammonium		

Table 12: Euro Trading

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Medister	No	Powder	In aqueous solution-hydrogen peroxide and peracetic acid	Italy	Euro Trading S.r.l. , Via Cassia per Siena 158, Loc.Bargino 50023 S.Casciano V.d.P. (Fi) Italy

Table 13: Germitec

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Antigermix	No	System	Ultraviolet C Rays (UVC)	France	Germitec (www.germitec.com)

Table 14: Laboratories Anios

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Anios Detergent Disinfectant	No	Liquid	Quaternary ammonium propionate, guanidinium acetate	France	Laboratoires Anios, (www.anios.com)
Anioxyde 1000	No	Liquid	Peracetic acid		
Hexanios G+R	No	Powder	Polyexanide, quaternary ammonium chloride		
Salvanios	No	Powder	Quaternary ammonium propionate, guanidinium acetate		
Steranios	No	Liquid	Glutaraldehyde		

Table 15: Metrex Research Corporation

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Cavicide	No	Spray	17.2% Isopropyl alcohol 1-5% 2-Butoxyethanol	USA	Metrex Research Corporation (www.metrex.com/company/contact/index.cfm)
Caviwipes	No	Wipes	Isopropyl alcohol, Benzethonium Chloride		
Compliance TM	No	Liquid	Hydrogen peroxide and peracetic acid		
Metricide®14	Yes	Liquid	Glutaraldehyde		
Metrizyme	No	Enzymes	Proteinase subtilisin <2%		
Omnicide 14NS	No	Liquid	2.6% Glutaraldehyde		
Omnicide 28NS	No	Liquid	2.5% Glutaraldehyde		

Table 16: Milton Pharamceutical

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Milton	No	Liquid	Sodium hypochlorite	France	Milton Pharmaceutical Ltd , (www.milton-tm.com/milton-baby-hygiene-products-contact.htm) or Laboratoires RIVADIS

Table 17: Nanosonic

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Trophon	Yes	System	-	Australia	Nanosonics Ltd. (www.nanosonics.com.au)

Table 18: Parker Laboratories

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Protex					
Transeptic	Yes	Spray	Isopropyl alcohol	USA	Parker Laboratories Inc. (www.parkerlabs.com)

Table 19: Pharmaceutical Innovation

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
T-Spray	No	Spray	IPA, quaternary ammonium chloride	USA	Pharmaceutical Innovations Inc (www.pharminnovations.com)
T-Spray II	No	Spray			

Table 20: Professional Disposables International

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Sani Cloth Active	No	Wipes	Quaternary ammonium compound, monoethanolamine	USA	Professional Disposables International, (www.nicepak.com/divisions/pdi/default.asp)
Sani Cloth HB	No	Wipes	0.07% alkyl dimethyl benzyl ammonium chloride 0.07% alkyl dimethyl ethylbenzyl ammonium chloride		
Sani Cloth Plus	No	Wipes	0.125% alkyl dimethyl benzyl ammonium chloride 0.125% alkyl dimethyl ethylbenzyl ammonium chloride		
Sani Cloth Super	No	Wipes	0.025% alkyl dimethyl benzyl ammonium chloride 0.25% alkyl dimethyl ethylbenzyl ammonium chloride		
Super Sani Cloth	No	Wipes	Quaternary ammonium, Isopropyl alcohol		

Table 21: Safetec

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
SaniZide Plus	No	Spray, Wipes	0.105% alkyl dimethyl benzyl ammonium chloride 0.105% alkyl dimethyl ethylbenzyl ammonium chloride	USA	Safetec (www.safetec.com)

Table 22: Schülke&Mayr

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Gigasept (10% vol.)	No	Liquid	4.5% Formaldehyde 4.5% Dimethoxytetrahydrofuran 6.8% Succindialdehyde 3-8% Methanol 25-45% Ethanol <5% Alcohols	Germany	Schülke&Mayr GmbH, (www.schuelkemayr.com/int/en/contact/smi044_adresses.htm)

Table 22: Schülke&Mayr

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Gigasept AF	No	Liquid	10% Dimethyldioctyl-ammonium chloride 15.6% Cocosalkylpropylendiaminbiguanidinium-diacetat 5-10% N-dodecylpropane-1,3-diamine 5-15% Tridecylpolyethylenglycoether <5% Ethanol <5% propan-2-ol		
Gigasept FF	No	Liquid	11.9% Succindialdehyde 3.2% Dimethoxytetrahydrofurane 5-15% Ethanol 5-10% Methanol 1-5% Alkyl-polyethylen-glycolpolypropylen-glycoether 1-5% 2-(2-hexyloxy-ethoxy)ethanol		
Mikrozid AF	No	Wipes	24% Ethanol 35% Propan-1-ol		
Mikrozid PAA	No	Wipes	Peracetic acid		
Mikrozid Sensitive	No	Wipes	0.26% benzyl-C12-16-alkyldimethylchlorides 0.26% Didecyldimethylammonium chloride 0.26% C12-14-alkyl [(ethylphenyl)methyl] dimethyl		

Table 23: Seton

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Medi-Prep	No	Wipes	1% Cetrimide BP	UK	Seton , Wildmere Industrial Estate, Oxfordshire, OX16 5JU, UK, E-mail: contacus@seton.co.uk

Table 24: Steris

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Klenzyme	No	Liquid	Enzymatic detergent	USA	STERIS Corporation (www.steris.com)
Prolystica	No	Liquid	Enzymes		
Resert HLD (ex Secure HLD)	No	Liquid	Hydrogen peroxide		
Resert XL HLD	No	Liquid	1.4-2.3% Hydrogen Peroxide <2.5% 2-Furoic acid		
Revital OX Resert XL HLD	No	Liquid	Accelerated hydrogen peroxide		
Steris V-Pro 1	No	System	Hydrogen peroxide		

Table 25: Tristel Solution Limited

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Tristel Duo	No	Spray	Chlorine dioxide	UK	Tristel Solution Limited (www.tristel.com)
Tristel Sporcidal Wipes	No	Wipes	Chlorine dioxide		

Table 25: Tristel Solution Limited

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Tristel Stella 5	No	System	Chlorine dioxide		

Table 26: Virox Technologies

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Accel RTU	No	Wipes	Hydrogen peroxide	Canada	Virox Technologies (www.virox.com)
PreCept RTU	No	Wipes	Hydrogen peroxide		

Table 27: Wave Energy Systems

Agent	FDA Approved	Presentation	Main ingredients	Country of manufacturer	Manufacturer
Wavicide-01®	Yes	Liquid	Glutaraldehyde	USA	Wave Energy Systems IN, 216 Little Falls Road, Ceder Grove NJ 07009 USA

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Tabella 1: Advanced Sterilization Products

Agente	Approv. FDA	Presentazione	Ingredienti principali	Paese di produzione	Produttore
Cidex®	No	Liquido	Glutaraldeide	USA	Advanced Sterilization Products (www.sterrad.com/about_asp/locations/index.asp)
Cidex Activated Dialdehyde Solution	No	Liquido	Soluzione alcalina di glutaraldeide		
Cidex® OPA	Sì	Liquido	Ortoftalaldeide (OPA)		
Cidexplus®	No	Liquido	Glutaraldeide		
CidezimeXTRA	No	Enzimi	Cellulasi <5% Lipasi <5% Amilasi <5% Proteasi <5% Borace decaidrato <8.5%		
DisOPA	No	Liquido	Ortho-phthalaldehyde (OPA)		
Enzol/Cydezime	No	Liquido	Detergente enzimatico		
NuCidex	No	Liquido	Acido peracetico		
Sterrad® NX	No	Liquido	Perossido d'idrogeno		
Sterrad® 100NX	No	Liquido	Perossido d'idrogeno		
Sterrad® 100S	No	Liquido	Perossido d'idrogeno		

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Tabela 5: BODE Chemie GmbH
Tabela 6: Borer Chemie AG
Tabela 12: Euro Trading S.r.l.
Tabela 15: Metrex research Corporation
Tabela 22: Schülke&Mayr GmbH
Tabela 24: Steris Corporation