



FY "M
May 88%

MyLab

58J5B798 CD9F5H=CBG

OPERATOR MANUAL

Doc # 81B16EN13

Warning Symbols


Warning symbols are used to identify hazardous materials and procedures. The symbols are used to identify hazardous materials and procedures. The symbols are used to identify hazardous materials and procedures.

Warning symbols are used to identify hazardous materials and procedures. The symbols are used to identify hazardous materials and procedures.

- Warning symbols are used to identify hazardous materials and procedures.
- Warning symbols are used to identify hazardous materials and procedures.
- Warning symbols are used to identify hazardous materials and procedures.
- Warning symbols are used to identify hazardous materials and procedures.
- Warning symbols are used to identify hazardous materials and procedures.
- Warning symbols are used to identify hazardous materials and procedures.

Warning symbols are used to identify hazardous materials and procedures. The symbols are used to identify hazardous materials and procedures.

Warning symbols are used to identify hazardous materials and procedures. The symbols are used to identify hazardous materials and procedures.



Warning symbols are used to identify hazardous materials and procedures. The symbols are used to identify hazardous materials and procedures.

Warning symbols are used to identify hazardous materials and procedures. The symbols are used to identify hazardous materials and procedures.

Be sure to understand and observe each of the cautions and warnings.

GC: HK 5F 9? 9MG G97 HCB

H\g`gWcb`lg`hYg`Zk`fY`_Yng`Uj`Uj`WY`b`h`Yj`Uf`ci`g`a`cXg`b`: fYnY`UbX`Z`f`f`Ya`ch`Wbfc`c`Zh`Yj`Xc`f`Wf`Xf`-h`Ug`[`j`Yg`ga`Yh`d`g`cb`ck`lc`i`g`Y`h`Y`Wbfc`g`Uj`Uj`WY`lc`cd`h`a`h`Y`ja`U`Yei`U`h`n`i`

H\Yg`Wcb`g`cf`[`Ub`h`X`U`g`Z`ck`g`:

€ 7\Udhf`%`6!`A`c`X`Y`

H\g`W`U`d`h`f`lg`h`Y`g`Z`k`f`Y`_`Y`n`g`U`j`U`j`W`Y`b`6!`A`c`X`Y`U`b`X`Y`d`U`b`g`ck`lc`i`g`Y`h`Y`a`lc`cd`h`a`h`Y`h`Y`ja`U`Y`

€ 7\Udhf`&`A`!`A`c`X`Y`

H\g`W`U`d`h`f`lg`h`Y`g`Z`k`f`Y`_`Y`n`g`U`j`U`j`W`Y`b`A`!`A`c`X`Y`

€ 7\Udhf`."`8`c`d`d`Y`f`

H\g`W`U`d`h`f`lg`h`Y`g`Z`k`f`Y`_`Y`n`g`U`j`U`j`W`Y`b`DK`8`c`d`d`Y`f`Z`c`f`a`U`h`U`b`X`Y`d`U`b`g`ck`lc`i`g`Y`h`Y`Wbfc`g`Uj`Uj`WY`lc`cd`h`a`h`Y`ja`U`Yei`U`h`n`i`

€ 7\Udhf`(.`7`c`c`f`:`ck`A`U`d`b`l`f`f`7`A`E`U`b`X`D`c`k`Y`f`8`c`d`d`Y`f`f`DK`F`8`L`

H\g`W`U`d`h`f`lg`h`Y`g`Z`k`f`Y`_`Y`n`g`U`j`U`j`W`Y`b`7`c`c`f`:`ck`A`U`d`b`l`U`b`X`Y`d`U`b`g`ck`lc`i`g`Y`h`Y`a`"

€ 7\Udhf`)`.:`f`Y`n`Y`

H\g`W`U`d`h`f`lg`h`Y`g`Z`k`f`Y`_`Y`n`g`U`j`U`j`W`Y`b`: fYnY`UbX`Y`d`U`b`g`ck`lc`i`g`Y`h`Y`a`"

€ 7\Udhf`*`6`c`X`m`A`U`f`_`g`

H\g`W`U`d`h`f`Y`d`U`b`g`ck`6`c`X`m`A`U`f`_`g`U`f`Y`c`f`[`Ub`h`X`U`b`X`ck`lc`i`g`Y`h`Y`a`"

€ 7\Udhf`+`5`b`b`c`h`U`j`c`b`g`

H\g`W`U`d`h`f`Y`d`U`b`g`ck`lc`i`g`Y`h`Y`U`b`c`h`U`j`c`b`Z`b`W`c`b`V`c`h`V`r`i`k`c`f`X`U`b`X`V`r`g`b`h`b`W`

€ 7\Udhf`,`I`g`f`D`f`Y`g`h`g`

H\g`W`U`d`h`f`lg`h`Y`d`U`f`a`Y`n`f`g`h`U`h`U`f`Y`U`j`U`j`W`Y`b`h`Y`I`g`f`D`f`Y`g`h`c`d`h`c`b`U`b`X`Y`d`U`b`g`ck`lc`W`b`Z`[`i`f`Y`h`Y`a`"

A m@UV ' ' 5 8 J 5 B 7 9 8 ' C D 9 F 5 H = C B G

Dqxlpq aa 8/4
 Qxct {"Crr rlec vqp aaa 8/4
 Ecplpg aaa 8/4
 Cdf qo kpciCrr rlec vqp aaa 8/4
 Ecplpg aaa 8/5
 Ectf lq"Crr rlec vqp aaa 8/5
 Ecplpg aaa 8/5
 Tgr tq"Crr rlec vqp aaa 8/5
 Gsvlpq aaa 8/6
 Cdf qo kpciCrr rlec vqp aaa 8/6
 Gsvlpq aaa 8/6
 Ectf lq"Crr rlec vqp aaa 8/6
 Gsvlpq aaa 8/6
 Qxct {"cpf "Tgr tq"Crr rlec vqp aaa 8/6
 Gsvlpq aaa 8/7
 Vgpf qp"Crr rlec vqp aaa 8/7
 Hgkpg aaa 8/7
 Cdf qo kpciCrr rlec vqp aaa 8/7
 Hgkpg aaa 8/8
 Tgr tq"Crr rlec vqp aaa 8/8
 Cplo cnHcto aaa 8/8
 Tgr tq"Crr rlec vqp aaa 8/8
 Uy lpg aaa 8/8
 Tgr tq"Crr rlec vqp aaa 8/8
 Qxlpq"cpf "Nico cu aaa 8/8
 Tgr tq"Crr rlec vqp aaa 8/8

7 - ANNOTATIONS _____ 7-1

Cppqv vqpu'lp'O {Ncd aaa 9/3
 I rquct {"d {"Y qtf aaa 9/4
 I rquct {"d {"Ugvpqpeg aaa 9/4

8 - USER PRESETS _____ 8-1

Wgt 'Rt gwg'v'O gpw'lp'O {Ncd aaa /3
 Et gc'g'c'Rt gwg'v'lp'O {Ncd aaa /4
 Rt gwg'v'K6 r qt vGzr qt v'lp'O {Ncd aaa /4
 Rt gwg'v'Rct co gvgf uaa /6
 Ugw'kpi "Qj gt'Rctco gvgf uaa /8
 Tgwg'v'Rt gwg'v'lp'O {Ncd aaa /8

1 - B-Mode


This chapter lists the software keys available in B-Mode format and how best to use them to optimize the image quality.

Activation of B-Mode Format

The system activates automatically the B-Mode format each time the system is switched on. The 2D format can be redisplayed when deactivating any other mode.

Control Panel Keys in B-Mode on MyLab

Upon activation, the following screen image will be visible:




In B-Mode, under the User environment, the following functions and commands are available:

GD97=9G`f Vtbt	7C@CF
5DD@=75H=CB	DK F`8
DF 9G9H	8CDD@9F
: F 9E I 9B7M	DF C; F 5A A 56@9 ? 9M
89DH<	DF C; F 5A A 56@9 ? 9M
H; 7	DF C; F 5A A 56@9 ? 9M
58>I GH	A CF 9I
: F 99N9	
=A 5; 9 CF `7@=D	

***MORE** menu for **IMAGE**, **REVERSE**, **B UPDATE**, **M MODE**, **ZOOM**, **B+B**, and **L/R**.

Control Keys in B-Mode on MyLab

Additional settings are present under the **SETTINGS** tab, the following command parameters are listed on the right of the screen:



The command parameters in Real Time 2D are organized in pages (on the right of the screen) that can be scrolled by use of the **SETTINGS** key at the bottom of the list. The following command keys are available:

H9=	: C75@DC=BHG
89BG=HM	GH99F
G=N9	DCK 9F
L!J =9K	LJ =9K `GA CCH<
; F 9MA 5D	LJ =9K `9B< 5B79
5A 6=9BH`@=: < H	LJ =9K 89H5=@
8MB5A =7 `F 5B; 9	A!J =9K
7C@CF =N9	HD!J =9K
G< 5F DB 9GG	CF =9BH5H=CB
D9F G=GH9B79	7@=D`8I F 5H=CB
	6F =; < HB 9GG

Tips for 2D Scanning with MyLab

2D Format Optimization

The **ORIENTATION** key can be used to change the sector orientation (normal/upside down), while the **REVERSE** key changes right/left or left/right orientation depending on the application.

- Use the **DEPTH** button to increase or decrease scanning depth.

The 2D can be enlarged; the enlargement factor is variable.

- Press **ZOOM** to display the ROI box
- Position the ROI box in the screen area by use of the sliders
- If necessary, the zoom level and therefore the ROI size can be modified.
- Press the **ZOOM** key to activate the selected Zoom zone
- Press **ZOOM** again to return to a normal 2D

To Enlarge 2D Area



HAY[!YUf`hYrca
 YYZ!YgUY`hY
 FC=k]`W


2D Display Optimization

First of all, the gain must be properly adjusted in order to clearly display the structures being examined; then, fine optimizations can be performed interacting with the display commands or with the acoustic parameters of the probe.


Acoustic Parameters

The **FREQUENCY** key allows the imaging frequency to quickly change (higher frequency to optimize resolution, lower frequency to increase penetration).

Use the **DEPTH** command to increase or decrease scanning depth.



The gain can be adjusted within the entire sector by the gain slider that will appear when the area beneath the presets is touched. Slide the button on the bar to the up/down direction to decrease/increase the gain.



H 7

To adjust the gain on individual areas of the sector use the **TGC** potentiometers.

The **TGC** potentiometers appear on the screen by touching the same area as for the gain when the **TGC** key is selected. After 5 seconds they are not touched they will disappear from the screen again. Moving the sliders to the left decreases and to the right increases the gain. With the **RESET** button all sliders can be reset in the centre position. The position of the TGC sliders can be saved for each probe under preset: set the TGC sliders position and then press **SAVE** followed by **STORE SETTING**.

In order to fully employ the dynamics of the system, ESAOTE recommends keeping the general gain at a relatively high level.



TEI mode can be activated on lower frequency probes by pressing the **TEI** key. This mode generally improves the brightness of the image by decreasing acoustic noise. Because of the non-linear response of tissues to ultrasound energy, **TEI** mode is not recommended for high frequency probes.

Note

Acoustic parameters and gain interact with one another; it may be necessary to review the adjustment of gain when an acoustic parameter changes.



DENSITY, available with linear and convex probes, optimizes image quality.




With the **SIZE** key the image sector size can be adapted. This function is available with convex-array probes only.

View Commands

These commands are mainly “subjective” and patient-dependent.



The **M-VIEW** key is displayed when either a Linear or Convex probe is used. To activate the MView mode, press the corresponding key: the B-Mode image will be the result of different bidimensional images acquired with different steering angles. Different MView values correspond both to different view lines and steering angles.



The **TP-VIEW** key is active when Linear Array probes are used. TP VIEW activates trapezoidal view, enlarging the probe's field of view.



The **X-VIEW** key allows a real-time image enhancement of tissue margins and tissue resolution to increase diagnostic confidence, eliminating speckle and noise artefacts. The available values are 1, 2 and C. X-view values 1 and 2 are fixed presets, while with X-view value C the X-view parameters can be adjusted using **XVIEW SMOOTH**, **XVIEW ENHANCE** and **XVIEW DETAIL**.



The **GREY MAP** key changes the post-processing map. The gray map can be replaced with one of the twelve chromatic scales to improve contrast, ranging from minimum to maximum contrast.



The **AMBIENT LIGHT** key allows tissue structures to be characterized, reacting to the compression of reflected echoes. Ten different levels can be selected: the lower the selected value, the greater the contrast.



The **DYN RANGE** key changes the dynamic range improving the image contrast



A chrominance scale can be selected using the **COLORIZE** key. The system has various different colorizations.




The **SHARPNESS** key accentuates the edges and the small differences in tissues. Five different levels can be selected.



The **PERSISTENCE** key allows the persistence level applied to the Real Time view to be controlled; increasing persistence levels increases the perception of the image, but decreases the discrimination of moving structures.



The focal point in transmission is positioned according to the initial depth, in order to optimize the entire analysis area. The focal point can be changed to increase the resolution and sensitivity of a specific area of the 2D; it is also possible to activate more than one focus in transmission and to increase resolution over a larger area.



Note

Several focuses in transmission can be activated; in this case, the relative distance between focuses is pre-established.


Pressing the **FOCAL POINTS** key will change the number of active focuses. 1 to 4 focal points can be chosen. The Focus Touch Technology gives control over the relocation of the focal points by simply touching the screen at any region of interest on the ultrasound image.



The **GH99F** key (active for linear probes only) modifies the slope of the B-Mode image.



The **POWER** key gives the possibility to adjust the acoustic output power between 10 and 100% of the maximum Acoustic Output allowed by the application.



HAY POWER Wbfc`' lg' i gX hc' WUb[Y hY hfUga]hX dckYf/ i gY hY
 a]b]a i a 'dck Y' Wa dUjVYk]h' UX]U bcg]WYj Y' cZh Y]a U Yg" =Zh YfY
 lg']bg' Z]M]bh gYb]h]j]hã a U_Y g fY hY [U]bž Z'W' dc]bh UbX dfcVY
 ZYei YbWã U YVYb gYhWffWãhVZcfY]bWUg]b] h Ydck Yf''



The **BRIGHTNESS** key gives the possibility to adjust the the brightness of the screen.




2

This chapter lists the software keys available in M-Mode format and explains how to best use the controls available to optimize the image quality.

M-Mode Analysis: cfa Uh

- Press **M-MODE** key to activate M-Mode analysis.
- Position the relative B-Mode line with the slider.
- Press **M-MODE** again to activate the M-Mode; the display will split showing the B-Mode and the M-Mode (rapid sequence of the B-Mode).



By pressing the **B UPDATE** key a real time update of the B-Mode image is displayed. Pressing the **M-MODE** key will return to the 2D screen.

7 cblfc`? Yng]b`A!A cXY`cb`An@UV

There are different command parameters in M-Mode. In the **SETTINGS** tab in M-Mode

In M-Mode the following functions and commands are available:

- X-VIEW**
- GREY MAP**
- AMBIENT LIGHT**
- DYNAMIC RANGE**
- COLORIZE**
- SHARPNESS**
- SWEEP TIME**
- POWER**

H]dgZcf`A!A cXY`GMbb]b[`k]h`An@UV

A!A cXY`8]gd`UriCdh]a]rU]cb

The adjustments for this mode are the same as those for B-Mode. In order to obtain a good M-Mode trace, it is fundamental to optimize the 2D image from which the trace will then be sampled. During the exam the **B UPDATE UPDATE** key freezes the trace acquisition and the reference 2D is temporarily re-activated. Normally, further interactions are not necessary.

View commands

The following keys: **X-VIEW**, **GREY MAP**, **AMBIENT LIGHT**, **DYNAMIC RANGE**, **COLORIZE**, **SHARPNESS**, **SWEEP TIME** and **POWER** are similar in function and operation to B-Mode.

The **SWEEP TIME** command can be used to change the sweeping speed of the M-Mode trace.

Refer to B-Mode for further explanation



3 - Doppler

This chapter lists the software keys available in PW Doppler and explains how to use the controls available to optimize image quality.

Activation of Doppler Format

DOPPLER

Doppler

- Press **DOPPLER** key to activate the sample volume on the 2D image.
- Position the sample volume with the sliders, by dragging it or by pointing on a preferred spot.
- Press **DOPPLER** key again to activate PW Doppler analysis.
- **IF** necessary, press **SETTINGS** tab to modify PW Doppler parameters.



Pressing the **DOPPLER** key will return to the 2D screen. During the exam the **B UPDATE** key freezes the trace acquisition and the reference 2D is temporarily re-activated. To switch back to the Doppler trace press the key once more.

Control Keys in Doppler on MyLab

In the **SETTINGS** tab, the Doppler command parameters can be adjusted:

The command parameters in Doppler are organized through pages (on the right of the screen) that can be scrolled by use of the **SETTINGS** key at the bottom of the list.

In Doppler mode, the following functions and commands are available:

PRF DOPPLER	GREY MAP
THETA ANGLE	COLORIZE
STEER	AUDIO
BASELINE	SWEEP TIME
WALL FILTER	TRACE TYPE
SAMPLE SIZE	TRACE METHOD
POWER	DYNAMIC RANGE
	PLEX

Tips for Doppler Scanning with MyLab

Doppler Format Optimization

The **BASELINE** key moves the baseline up or down. The **BASELINE** function of the Doppler frequency spectrum presentation can also be used directly from the **USER** page. By just touching the Doppler spectrum area the baseline can be set directly to the preferred level.

If the application defaults to the correction-angle, use the **ANGLE** control to align the angle vector with the flow direction. The **STEER** key allows to orient the Doppler line.

The **SWEEP TIME** command can be used to change the scanning speed.

Doppler Optimization

The gain must first be optimized until a clear envelope of the spectral analysis is obtained; the **WALL FILTER** must be set in order to eliminate wrong low-speed signals caused by moving structures. Interaction with other commands or the acoustic parameters further improves the spectrum quality.



Baseline



Theta Angle



Sweep Time



Wall Filter

Gain and Wall Filters

This knob affects the Doppler video and audio components.



Note

The Doppler volume can be changed independently by using the **AUDIO** key (under **SETTINGS** tab).



Display Commands

The **GREY MAP** key allows the post-processing map to be changed. The gray map can be replaced using one of the six chromatic scales (**COLORIZE**) to improve contrast.



In PW Doppler, the **SAMPLE SIZE** key allows the sizes of the sample volume to be changed.


Acoustic Parameters

The **FREQUENCY** key allows the **Doppler frequency** to be changed: lower frequency increases penetration.



Refer to the "Safety and Standards" manual for further information.

The **POWER** control is used to change transmitted power; use the minimum power compatible with a diagnostic level of the images. If there is insufficient sensitivity, the user should ensure that the gain, focus and probe frequencies have been set correctly before increasing the power.



The **PRF** key activates the Doppler PRF, allowing to increase the available maximum PRF value to measure higher velocities by using more sample volumes.

Real time Doppler measurements



TRACE TYPE enables the calculations of parameters during real time PW trace. It also defines how the trace is followed, by Maximum values, mean values or both.



TRACE METHOD defines the part of the spectrum that has to be traced for real-time measurements calculation. It is possible to trace the full spectrum, or just the positive or the negative part.

CAUTION

Automatic Doppler measurements represent just an immediate method to get a general idea of the importance of the pathology under examination. For a precise pathology evaluation we recommended to use the application calculation package. - For a correct diagnostic evaluation, it is recommended to use the angle correction factor, in order to obtain the right flow alignment. Make sure that the profile of the automatically detected Doppler flow (yellow line) corresponds to the real profile.

The **DYNAMIC RANGE** key changes the dynamic range improving the image contrast.



The **PLEX** key updates or freezes the reference 2D image, keeping the trace in real time.

4- Color Flow Mapping (CFM) and Power Doppler (PWR D)


This chapter lists the software keys available in CFM and explains how they are used to optimize the image.

Activation of Color Flow Mapping Format

COLOR

Color

- Press **COLOR**.
- Use the sliders to move and change the area of the color box. or change the position by dragging.



Pressing the **COLOR** key will return to the 2D screen.

Control Keys in CFM on MyLab

In the **SETTINGS** tab the command parameters for Real Time Color Flow Mapping mode can be adjusted:

he command parameters in Real Time Color Flow Mapping are organized through pages (on the right of the screen) that can be scrolled by use of the **SETTINGS** key at the bottom of the list.

The following command keys are available:

PRF COLOR	MAS THRESHOLD
SMOOTH	COLORMAP
STEER	C CORRELATION
BASELINE	AXIAL INTEGRATION
WALL FILTER	POWER
PERSISTENCE	ORIENTATION
SENSITIVITY	

WARNING

When the steering is set to the maximum step, color dots could be displayed because of artefacts. Should this happen, reduce the steering by one step.

Tips for Color Flow Mapping Scanning with MyLab



Baseline

WARNING

Optimization of the Color Flow Mapping Format

The **BASELINE** key is used to move the zero line up or down. The **STEER** key allows changing the steering of the color box.

When the steering is set to maximum step, color dots could be displayed (artifacts). Should this happen, reduce the steering by one step.

The width of the ROI CFM and the B-Mode angle must be as small as possible in order to maximize the CFM frame rate.

Note

The B-Mode angle (**SIZE** key) can be reduced to maximize the frame rate and to enlarge the analysis area using the **ZOOM** key.

Display Optimization

First of all the operator must adjust the gain, to optimize the display of the analysis area. The commands and the acoustic parameters controls allow then to further refine the signal.

Gain

The user is recommended to adjust the gain to obtain the most useful signal level.

Note

Color mode gain must be properly adjusted in order to obtain a good CFM signal; excessive gain may “mask” the flow.



Color Map

View Commands

The **COLOR MAP** key allows a different CFM map to be selected: ten different color maps are available.

The **PRF** and **DENSITY** keys affect Color “filling”.

Control	CFM filling
PRF ↑	↓
PRF ↓	↑
DENSITY ↑	↑
DENSITY ↓	↓

Maximum PRF is proportional to the B-Mode depth



Color sensitivity can be adjusted using the **SENSITIVITY** control. Three levels are available (Low, Medium and High). Flow continuity increases with high persistence levels and reduce the system’s frame rate; the **PERSISTENCE** option optimizes this parameter. **WALL FILTER** can be used to reduce the artefacts caused by acoustic decoupling or moving structures. The **SMOOTH** key makes the flow representation homogenous.

Acoustic Parameters

Only one focal point is active in transmission in CFM, despite the B-Mode settings and is automatically positioned at the centre of the ROI CFM. The CFM frequency can be changed using the **FREQUENCY** key; the higher frequency helps to show low speeds. With the **POWER** key the acoustic output can be adjusted.



Refer to the “Safety and Standards” manual for further information.


The **POWER** control is used to change transmitted power; use the minimum power compatible with a diagnostic level of the images. If there is insufficient sensitivity, the user should ensure that the gain, focus and probe frequencies have been set correctly before increasing the power.

5 - Freeze

This chapter lists the software keys available in freeze and explains how to use the keys.

Control Panel Keys in Freeze on MyLab

Pressing the **FREEZE** key freezes the image and displays the following screen:



Functions available when Freeze mode is active are: **MARK**, **CALIPER**, **IMAGE**, **CLIP**, **PRINT**, **ANNOT**, **ZOOM** and **CINE**. In certain applications, where dedicated measurements are available, the key **MEASURE** is also available.

*Memory scrolling
cursor*

After freezing the image, the system displays the memory-scrolling bar, where the images acquired immediately prior to putting the system into freeze mode are temporarily saved and can be recalled on screen.

Use of the Control Panel Keys with MyLab



Begin/End Key

The **BEGIN/END** key is used to automatically position the user at the start or end of the sequence.

Moving the bar scrolls through the images acquired before activating the freeze mode.



Stop and Play Keys

When the **PLAY** mode is active, the sequence of stored images is reviewed kinetically. The **STOP** key stops play mode and keeps the latest played image on the screen.

6 - Body Marks

This chapter explains how body marks are organized and how to use them.

Body Marks

Body marks are schematic drawings of anatomical sections. A vector overlays the mark to indicate the probe's position. The active body mark is displayed at the bottom left of the screen.


Body marks are organized in groups: each species and application has its specific set of body marks.

Body Marks Activation

The **MARK** key allows adding a bodymark to the current image. When selecting the **MARK** key a pattern list opens on the left of the screen with the bodymarks available for the current application. Body marks can be activated and displayed in **FREEZE** mode.

Control Panel Keys

The following figure gives an overview of the options:



Number	Icon
1	Change group of bodymarks
2	Select bodymark
3	Confirm bodymark
4	Define left and right bodymark for B+B mode in all other modes, only the left is shown.
5	Position arrow on bodymark
6	Change angle of arrow on bodymark
7	Cancel bodymark

The body marks list is display on the left of the screen and can be scrolled to review all the body marks available.



To choose a different group of marks, use the **CHANGE GROUP** key. Beneath the icon available groups will be listed. Select the preferred group.

Use the sliders for positioning the arrow in the horizontal and vertical direction and adjust the angle of the arrow. The arrow can also be positioned by touching the wanted position within the bodymark

The user can alter the orientation of the image by choosing the **LEFT** or **RIGHT** key, which will alter the color of the arrow respectively into green of yellow.

Once the icon has been selected and the arrow positioned, press the **OK** key to confirm the body mark and arrow position on screen.



Select the **MARK OFF** key or press **CANCEL** to exit without displaying any body mark.

Species and Application Body Marks

Bovine Ovary Application



Canine Abdominal Application





Uterus Ovaries



Stomach



Liver



Kidney



All



Pancreas



Gallbladder

**Canine
Cardio Application**



Heart general



Short axis



Short axis PV



5 chamber view



4 chamber view



Long axis



Short axis MV



Short axis AV



2 chamber view



Aorta

**Canine
Repro Application**



General



Right



Left



Back



Uterus Ovaries



Equine

**Equine
Abdominal Application**



General



Stomach



Liver



Kidney



Head



Pancreas



Gallbladder



All

**Equine
Cardio Application**



Heart general



Short axis



Short axis PV



5 chamber view



4 chamber view



Long axis



Short axis MV



Short axis AV



2 chamber view



Aorta

**Equine
Ovary and Repro Application**



General



Uterus Ovaries

**Equine
Tendon Application**



General



Front limb



Hind limb

**Feline
Abdominal Application**



Feline



General



Left



Right



Back



Uterus Ovaries



All



Stomach



Pancreas



Liver




Gallbladder



Kidney

**Feline
Repro Application**



**Animal Farm
Repro Application**



**Ovine and Llamas
Repro Applications**



7 - Annotations

This chapter explains how to use the free text annotation function and how to activate and use the glossary available in the annotation session.

Annotations on MyLab

ANNOT

When the **FREEZE** key is pressed during the exam, the **ANNOT** key will show up. The **ANNOT** key allows the user to work with a configurable glossary:

Below is the screen in annotation mode.



In annotation mode the word keys are displayed on the left of the screen. The key **BY WORD** respectively activates the by word glossary associated to the application. The term keys (**E1**, **E2**, **E3** and **E4**) are displayed at the bottom of the list when the glossary by sentence is active. Each **TERM** key contains a specific glossary, with which a sentence can be made.

The text home position and Annotation mode can be changed in the Annotation Setup menu (**MENU – MISCELLANEOUS – ANNOTATION SETUP**).

The words can be repositioned in any direction with the touch screen, preferably with the stylus.

Glossary by Word

On the left of the image the system displays the list of the available words. To select a word follow this procedure:

Procedure

- Scroll the list through and select the desired word (highlighted in blue).
- The selected word will be displayed on screen in yellow. The word can be placed on the screen by the sliders or by pressing the position on the screen where the word is desired.
- Press **OK** to confirm the selected word and close the annotation environment or press **NEW** to confirm the word and add a new annotation word on the screen.

Use the **KEYBOARD** key to activate the text entry. The keyboard will appear and it is possible to modify a selected word on the screen or edit a new word not included in the list. Pressing **ENTER** confirms the editing and closes the keyboard.

The **POINTER** key gives the possibility to add an annotation in the form of a pointer. This pointer can be an arrow or a circle. Select the region of interest on the touchscreen to place the pointer.

Use the **KEYBOARD** key to activate the text entry with the selected word in the text entry field; use the **FREE** key to enter text in an empty text entry field.

The **CLEAR** key cancels the text without exiting while the **CANCEL** key closes the annotation environment without confirmation.

The **UNDO** key cancels the last operation.

Pressing the **OK** key confirms and closes the annotation session.

According to the annotation setup, the text on the screen can be removed or remain visible when turning to live mode.

The **CLEAR** key deletes all the text on screen while the **CANCEL** key closes the annotation session without confirmation. If the “Hide with live image” field in the “**MENU, MISCELLANEOUS, ANNOTATION SETUP** and **GENERIC**” option (**SETTINGS** key) is activated, the text is automatically erased as soon as real time is resumed.

The **UNDO** key deletes the last entered word / sentence or action.

Glossary by Sentence

On the screen the system displays the sentence that the system automatically composes using the first words of the term lists. The list of the available words for

the first term of the sentence is displayed on the left of the image. The four terms composing the sentence are displayed on the bottom of the list, one key for each term: press the keys to select the desired term.

Procedure

- Scroll the lists using the **TERM** keys and select the desired words (highlighted in blue). The sentence will be displayed in yellow on the screen and, selecting different words through the terms, the sentence will be automatically updated.
- Press a point on the screen or use the sliders to change the position of the sentence.
- Pressing **NEW** confirms the last sentence and starts a new one.
- Pressing **OK** closes the annotation environment with confirmation.
- The sentence can be modified or completely edited by the **KEYBOARD** key that activates text editing.
- Press **ENTER** to close the keyboard and confirm.

Use the **KEYBOARD** key to activate the text entry with the selected word in the text entry field; use the **FREE** key to enter text in an empty text entry field.

The **CLEAR** key deletes all the text on screen while the **CANCEL** key closes the annotation session without confirmation. If the “Hide with live image” field in the “**MENU, MISCELLANEOUS, ANNOTATION SETUP** and **GENERIC**” option (**SETTINGS** tab) is activated, the text is automatically erased as soon as real time is resumed.

The **UNDO** key deletes the last entered word / sentence or action.

8 - User Presets

This chapter lists and explains the menu available for creating and saving presets in MyLab.

User Presets Menu on MyLab

The preset selection is made by the preset button positioned on top left of the screen, under the Species and Application buttons.

The system permits to adjust parameters in different modes and save those configurations under user-defined presets.



Create a Preset on MyLab

This procedure enables the creation of a preset on **MyLab**.

Pressing the **SETTINGS** tab displays the parameters available in the active operation mode.

The user can modify the parameters and save the preferred settings by pressing the **SAVE** key.

Pressing the **SAVE** key opens a menu for preset management on the right side of the screen.

- The **STORE SETTING** key will save the parameters configuration in the current preset.
- The **MAKE STARTUP** key defaults the current preset to be active as the system starts up.
- The **LOAD FACTORY** key reloads the factory presets erasing the eventual customizations made previously on the presets.
- The **RENAME SETTINGS** key will enable the changing of names of the presets



With the **LOAD FROM SETTING** key, the user can choose to copy the parameters from a preset (selected by the user in the preset list) into the current preset.

When the **LOAD FROM SETTING** key is pressed, the Copy Preset menu appears on the left side of the screen. The user can then select the preset whose settings have to be imported and copied into the current one.

To exit the procedure without performing any selection press the **CANCEL** key.

Preset Import/Export on MyLab


Presets can be imported/exported on/from **MyLab** units through USB memory media.

Procedure to save presets to USB

Presets can be saved to a USB medium by the following procedure:

- Connect a USB memory medium to the system.
- Press the **SAVE TO USB** key (or the **SAVE ALL TO USB** key to save all current presets at once).
- The system displays the list of the presets in the system. When preferred, the user can change the name of the preset by keyboard. The keyboard appears when touching the field of the name.
- Press the **OK** key to finalize the procedure.

The following screen will be visible when pressing the **SAVE TO USB** key:


**Procedure to load presets from USB**

Presets can be loaded from a USB memory medium by the following procedure:

- Connect a USB medium to the system.
- Press the **LOAD FROM USB** key (or the **LOAD ALL FROM USB** key to load all presets at once).
- The system displays a list of the presets already saved on a USB medium.
- Highlight the preset to be loaded (information about the file will be visible on the right side of the window) and press **OK** to load the preset.

Repeat the operation for all the desired presets to be loaded.

The following screen will be visible when pressing the **LOAD FROM USB** key:



CAUTION

Presets can only be saved to USB when the presets are stored first. First, press the **STORE SETTING** key to save the presets prior to using the **SAVE TO USB** key command.

The **CLOSE** key exits the menu without saving any modifications and returns to the Preset menu.

Note

The system allows saving presets with the same name: Check therefore that the assigned name is not already available in the list.

Preset Parameters

Different parameters can be set for each mode, by activating the **SETTINGS** tab.

The following tables list the value or initial status of parameters that can be set mode by mode.

B-Mode

Field	Action
XVIEW	Sets the initial value of the XView algorithm.
GREY MAP	Sets the grey scale active in 2D.
AMBIENT LIGHT	Sets the ambient light settings
DYN RANGE	Sets the dynamic range.
COLORIZE	Sets the chromatic scale active in 2D.
SHARPNESS	Sets the image's sharpness value (low, average, high).

Field	Action
PERSISTENCE	Sets the persistence level.
FOCAL POINTS	Number of focuses active in transmission.
STEER	Number of angels for the steering (Linear probes).
POWER	Sets the acoustic output power.
XVIEW SMOOTH	Finetuning of XView algorithm (only available when XVIEW is set to c).
XVIEW ENHANCE	Finetuning of XView algorithm (only available when XVIEW is set to c).
XVIEW DETAIL	Finetuning of XView algorithm (only available when XVIEW is set to c).
M VIEW	Sets the initial value of the Mview mode. Different Mview values correspond both to different view lines and steering angles.
TP VIEW	TP View sets the status (active/inactive) of the trapezoidal view on linear probes.
ORIENTATION	Image orientation (up/down).
CLIPS DURATION	Sets the length of clips for live procedures. Clip duration can be to: 1 s, 3 s, 5 s, 10 s, 20 s or 30 seconds.

COLOR / PWR D

Field	Action
PRF COLOR	Sets the PRF.
SMOOTH	Sets the smooth level.
STEER	Number of angels for the steering (Linear probes).
BASELINE	Sets the position of the Baseline.
WALL FITLER	Sets the wall filter value.
PERSISTENCE	Sets the Persistence level.
SENSITIVITY	Sets the sensitivity level.
MAS THRESHOLD	Moving artefact suppression (movement of transducer, breathing of patient, etc.)
COLORMAP	Sets the colors of the colorscale (maps Doppler frequency to color)
C CORRELATION	Averaging length of the color correlation
AXIAL INTEGRATION	Averaging in depth over color pixels
POWER	Sets the acoustic output power.
ORIENTATION	Image orientation (up/down).

M-Mode

Doppler

Field	Action
PRF DOPPLER	Sets the PRF.
THETA ANGLE	Sets whether the angle correction factor is oriented with a freely configurable angle correction factor.
STEER	Number of angels for the steering (Linear probes).
POWER	Sets the acoustic output power.
WALL FILTER	Sets the wall filter value in Pulsed Wave Doppler.
SAMPLE SIZE	Sets the default value of the sample volume.
SWEEP TIME	Sets the sweeping speed of the M-Mode spectrum
GREY MAP	Sets the scale of the greys active in M-Mode.
COLORIZE	Sets the chromatic scale active in M-Mode.
AUDIO	Sets the initial sound value.
BASELINE	Sets the position of the Base Line.
TRACE TYPE	Sets the method for Automatic Doppler Measurements (ADM)
TRACE METHOD	Sets the trace direction for ADM
DYNAMIC RANGE	Sets the dynamic range.
PLEX	Updates or freezes the reference 2D image, keeping the trace in real time.

Setting Other Parameters

The initial values of the parameters listed in the following table can be set for each application.

Field	Action
DEPTH (cm)	Sector (image) depth.
SIZE	Image sector angle (for convex probes)
FREQUENCY (MHz)	Imaging, CFM or Doppler frequency or TEI mode.
REVERSE	Sector orientation (right/left).
DENSITY	Image density (High/ Low)

Imaging Frequency

If available, the TEI (Tissue Enhancement Imaging) mode can also be selected.

Setting the Gains

For each mode (Imaging, CFM, Doppler) the initial gain value can be set.

Reset Preset on MyLab

Use the **LOAD FACTORY** key to set the factory presets with their default parameters.

CALCULATIONS SECTION

This section explains how to use the calculations packages offered by **MyLab**. The section is organized as follows:

- **Chapter 1: General Information and Generic Measurements**
This chapter provides general information about the proper manner in which to take a measurement and lists the generic measurements available in each application.
- **Chapter 2: Cardiology Calculations**
This chapter lists all the measurements, formulas and accuracies available in the Cardiac calculations packages.
- **Chapter 3: Abdominal Measurement Package**
This chapter lists measurements included in the Abdominal measurement package.
- **Chapter 4: Equine Tendon Package**
This chapter lists measurements included in the Equine Tendon package.
- **Chapter 5: Reproductive Calculation Packages**
This chapter lists measurements and formulas available in the reproductive packages.
- **Chapter 6: Animal Science Calculation Packages**
This chapter lists measurements and formulas available in the animal science packages.
- **Chapter 7: Measurements Accuracy**
This chapter provides the accuracy of measurements taken with **MyLab**.

Table of Contents

1 - GENERAL INFORMATION AND GENERIC MEASUREMENTS	1-1
General Information	1-1
Diagnosis Based on Measurements	1-2
Activating Generic Measurements	1-2
How to Take Measurements	1-3
Flow Measurements	1-4
Generic Measurements	1-4
2- CARDIAC CALCULATIONS	1
Cardiac Application	1
Cardiac Measurements	1
B-Mode Measurements	1
Ejection Fraction (Simpson-Biplane)	2
Ejection Fraction (Area-Length)	3
Fractional Area Change (FAC)	3
Left Ventricle	4
Left Ventricle Outflow Tract	5
Aorta	5
Left Atrium	5
Right Ventricle	5
Pulmonary Artery	5
Mitral	6
Cardiac Calculations in M-Mode	6
Left Ventricle	6
Left Ventricle Mass	6
Ejection Fraction	7
Aorta and Left Atrium	7
Mitral	8
Cardiac Calculations in Doppler	8
Mitral	8
Mitral Regurgitation	9
Aorta	9
Aortic Area	9
Aortic Regurgitation	10
Descending Aorta	10
LVOT Flow	10
Tricuspid	10
Tricuspid Regurgitation	11
Pulmonary Veins	11
Pulmonary Artery	11
Pulmonary Regurgitation	12
Cardiac Output - LVOT	12
Cardiac Output - Aorta	12
Cardiac Output - Pulmonary	13
Setup	14
Cardiac Report	15
B-Mode	17
M-Mode	17
Doppler	18
Additional	19
Bibliographic References and Formulas in B-Mode	20

Bibliographic References and Formulas in M-Mode	23
Bibliographic References and Formulas in Doppler	25
Formulas and Bibliographic References for Body Surface Area in Canine and Feline Cardio	31
3 - ABDOMINAL MEASUREMENTS PACKAGE	1
Abdominal Application	1
Abdominal Measurements	2
B-Mode Measurements	2
Doppler Measurements	6
Setup	7
Abdominal Report	8
B-Mode	10
Doppler	11
Additional	11
4 - EQUINE MSK / TENDON PACKAGE	1
Equine MSK/Tendon Application	1
Equine Tendon Measurements	1
B-Mode Measurements	1
Setup	7
DISTANCE METHOD	8
ZONE METHOD	8
MSK / Tendon Report	9
5 - REPRODUCTIVE CALCULATION PACKAGES	1
Reproductive Calculations in B-Mode	2
Gestational Age in Reproductive Canine Application	2
Gestational Age in Reproductive Feline Application	3
Gestational Age in Reproductive Equine Application	5
Gestational Age in Reproductive Bovine Application	6
Gestational Age in Reproductive Porcine Application	7
Gestational Age in Reproductive Ovine Application	8
Gestational Age in Reproductive Llamas Application	9
Setup	9
Gestational Age	10
General settings	12
Reproductive Report	13
Graphs Display	15
6 - ANIMAL SCIENCE PACKAGE	1
Animal Science Calculations in B-Mode	3
Yield Grade in Bovine Animal Science Application	3
Yield Grade in Porcine Animal Science Application	4
Yield Grade in Ovine Animal Science Application	5
QUIP index in Bovine Animal Science Application	6
Setup	7
Animal Science Report	8
Animal Science measurements, explanation	11
BEEF MEASUREMENTS	12

SWINE MEASUREMENTS	14
SHEEP MEASUREMENTS	15
7- MEASUREMENTS ACCURACY	1
Derived Data	2

1 - General Information and Generic Measurements

Refer to "Getting Started" manual for applications available on MyLab model.

This chapter provides general information about performing measurements with MyLab.

This chapter also lists the generic measurements available in the applications and how to use them. Generic measurements enable the user to quickly take measurements such as distance, area, volume and angle.

General Information

Measurements can be taken on frozen, stored and archived images. The available measurements are displayed at the right of the image. Messages displayed on the screen provide information guiding the operator in the procedure for making the measurement. The results are displayed on the left of the image.

Clips are compressed for digital storage. Compressed files involve a minimal loss of information (see specifications). Image features, if compared to the original, may not be optimal for the reporting functions.

To select the views and the positioning of the cursors, Esaote urges the operator to act according to current medical practice and to the instructions of specialists in this subject.

Note

Always enlarge the format to maximize the structure/signal to be measured.

The system cannot be used to measure images with ambiguous calibrations. An error message is shown on such images when the measurement is made.

CAUTION

If measurements are taken in dual mode (B+B) make sure that the images are aligned before taking measurements.

Diagnosis Based on Measurements

MyLab calculation packages have to be used by qualified personnel as a diagnostic tool. The diagnosis does not have to be based on the measurements only, but these are to be integrated with other clinical data.

All formulas of **MyLab** advanced calculation packages refer to a number of clinical bibliographic references that are listed for each application in the corresponding chapters of this section. Users are kindly encouraged to consult the original references to draw their conclusions on clinical consistency of the measurements.

Note


The user is responsible for customizations on measurements and calculations.

Activating Generic Measurements



When in freeze mode, the **CALIPER** key activates the generic measurements menu. The system displays the list of available measurements on the right of the screen, which are automatically identified according to the active mode and application.

The following screen will be available when Generic Measurements are enabled:



bottom of the list.

Some measurements include sub parameters (identified by the symbol **+**). To display the measurements included in a group, press the symbol and the measurements will be visible.

The user can modify the list of measurements (displayed when the **CALIPER** key is pressed) by accessing the **SYSTEM SETUP – GENERIC MEASUREMENTS** in the **MENU**.

How to Take Measurements

The **CALIPER** key is used to quickly select the measurement required. The measurement highlighted is then operative.

Following the instructions on the screen, position the touch screen stylus and press on the screen to confirm the position. The measurement cross is displayed just above the point pressed on the screen, in this way the measurement cross is always visible, hence it can be positioned more accurately.

The **UNDO** key can be used to re-start a measurement before it has been confirmed. This key is also used to erase a dotted line, point by point.

Some measurements (Length by vertex for example) require the **FINISH** key to end and confirm the measurement.

The **DELETE** key deletes a selected measurement.

The **ADD TO REPORT** key adds the generic measurement to the report.

To move a measurement, select it with the stylus. A box will appear around the measurement point, select any point within the box and then drag the measurement point to a preferred position.

The value being measured is displayed in real time on the left of the image.

How to Select a Measurement

- Freeze the image and press **CALIPER**.
- Select the required measurement of a group by pressing the key or using the touch screen stylus (the selected group is displayed on top of the list).
- To follow the entire measurement procedure, press the selected group.

To measure a specific parameter from a measurement group, press the **+** symbol to expand, and then select the required item.

Selective Clearing of a Measurement

- Position the stylus on the measurement to be cleared (the measurement is displayed in white).
- Press the **DELETE** key to clear the selected measurement.

The **CLEAR** key cancels all the measurements and the values displayed in the measurements field from the screen.

Flow Measurements

A flow measurement is structured into two stages.

Performing a Measurement

- Acquire a Doppler trace and press **FREEZE**.
- Press **CALIPER** to activate the calculations menu.
- Select the desired flow measure.
- Follow the instructions displayed on the screen to trace the velocity profile.

Generic Measurements

The tables below list the measurements available in each mode.

2D Measurements

Parameter	Calculation	Measurement	Displayed results
Distance	Distance	Distance	Distance
Distance Ratio	Distance ratio	Two distances	Two distances, Distance1/Distance2
% Diam	Diameter reduction	Two distances	Two distances, Δ Distance/Distance1
Vx-Length	Length (approximately straight)	More distances	Global distance
Tr-Length	Length (Profile)	Distance	Distance
A-Area	Area (Ellipse)	Distance, Area	Area, Perimeter
Vx—Area	Area (approximately straight)	More distances	Area, Perimeter
Tr-Area	Area (Profile)	Profile	Area, Perimeter
A-Ratio	Areas ratio	Two	Two areas, Area1/Area2
% Area	Reduction area	Due areas (on Profile)	Two areas, Δ Area/Area1
El-Volume	Volume (Ellipse)	Distance, Area	Area, Volume
Tr-Volume	Volume (Profile)	Profile, Distance	Area, Distance, Volume
Bi-Volume	Volume (Axes)	Three distances	Three distances, Volume
Angle	Angle	Angle	α/β angle

TAM Measurements (M-Mode and Doppler)

Parameter	Calculation	Measurement	Displayed results
Distance	Distance	Distance	Distance
D-Ratio	Distance ratio	Two distances	Two distances, Distance1/Distance2
Time	Time	Time	Time
% Time	Time ratio	Two Times	Two times, Time1/Time2
HR	Heart rate	Distance	R-R Interval, Heart rate
Velocity	Velocity	Velocity	Distance, Time, Velocity
% Veloc	Velocity ratio	Two velocities	Two velocities, Velocity1/Velocity2
FVI	Vascular FVI	Spectral envelope	FVI, Minimum, mean and maximum velocity
SVel/DVel	Systolic and Diastolic Velocities ratio	Two velocities	Systolic velocity, Diastolic velocity, Systolic velocity/Diastolic velocity
RI	Resistive index	Two velocities	Two velocities, Resistive index
Slope	Slope	Distance	Acceleration, PHT

2- Cardiac Calculations



Canine



Feline




Equine

Cardiac application
icon

This chapter describes the Cardiac application package (measurements and report environments) and lists the available measurements.

The Cardiac application package applies for the following species: canine, feline equine.



Note

This application calculation package is only available when the corresponding license is activated. Please refer to the System Configuration manual for all the options/licenses of the **MyLab** and their status.

Cardiac Application

The Cardiac application can be activated selecting the **CARDIO** key from the application list. When the Cardiac application is selected, in freeze mode, press the **MEASURE** key to activate the cardiac measurements.

Cardiac Measurements

B-Mode Measurements

Pressing the **MEASURE** button on the keyboard of the system in B Mode makes the list of the cardiac measurements appear on the right side of the screen.

The measurements are listed below:

Ejection Fraction (Simpson-Biplane)

EF (SIMP SI-B)

The Simpson method requires that measurements be taken on two cardiac views (Apical four chambers, A4C; Apical two chambers, A2C). The parameters that can be measured are:

Parameter		Abbreviation	Measurement
4CAd	4C Diastolic area	4CAd	Profile, Distance
4CAs	4C Systolic area	4CAs	Profile, Distance
2CAd	2C Diastolic area	2CAd	Profile, Distance
2CAs	2C Systolic area	2CAs	Profile, Distance

Since the Simpson method requires measurements on two different cardiac views, it is recommended that you use a Dual 2D so that both views are available. The group parameters are measured in four phases.

Performing the Measurements

- Acquire a B-Mode image and press **FREEZE**.
- If necessary, scroll memories to select the end-diastolic A4C image and press **MEASURE** to activate the menu.
- Select the “EF (Simpson)” group
- Follow the instructions on the screen to trace the diastolic profile and the LV long axis.
- Press the touchscreen to select the end-systolic image in A4C
- Follow the instructions on the screen to complete the measurements.
- Repeat the phases in A2C.

If an A2C image has not been stored, it is possible to return to Real Time with **B-MODE** to complete the acquisition. Freeze the image and press **MEASURE** to complete the measurement.

These parameters are automatically calculated once the measurements have been taken:

Parameter	
LVVd	LV diastolic volume
LVVs	LV systolic volume
EF	Ejection fraction
SV	Stroke volume
SI	Stroke index
HR	Heart rate
CO	Cardiac output
CI	Cardiac index

LV: Left ventricle

Ejection Fraction (Area-Length)

EF (A-L)

The ejection fraction measurement taken using the Area-Length method requires the measurement of diastolic and systolic parameters. The images must be from the same cardiac cycle. The following parameters can be measured:

Parameter		Abbreviation	Measurement
DIA AREA	LV diastolic area	LVAd	Profile
DIA AXISI	LV diastolic axis	LVLd	Distance
SYS AREA	LV systolic area	LVAs	Profile
SYS AXISI	LV systolic axis	LVLs	Distance

LV: Left ventricle

Ejection Fraction can be calculated on a single-frame 2D or on a Dual 2D image. The Dual 2D (particularly in simultaneous mode) is recommended. The group parameters are measured in two phases.

Performing the measurements

- Acquire a cardiac cycle and press **FREEZE**.
- If necessary, scroll memories to select the tele-diastolic image and press **MEASURE** to activate the menu.
- Select the “EF (Area-Length)” group
- Follow the instructions on the screen to trace the diastolic profile and the Ls long axis.
- Select the tele-systolic image and follow the instructions on the screen to complete the measurements.

If an entire cardiac cycle has not been stored, it is possible to return to Real Time with **B-MODE** to complete the acquisition. Freeze the image and press **MEASURE** to complete the measurement.

The following parameters are automatically calculated once the measurements have been taken:

Parameter	
LVVd	LV diastolic volume
LVVs	LV systolic volume
EF	Ejection fraction
SV	Stroke volume
SI	Stroke Index
HR	Heart rate
CO	Cardiac output
CI	Cardiac Index

LV: Left ventricle

Fractional Area Change (FAC)

FAC

The images must be from the same cardiac cycle. The following parameters can be measured:

Parameter		Abbreviation	Measurement
DIA AREA	LV diastolic area	LVAd	Profile
SYS AREA	LV systolic area	LVAs	Profile

LV: Left ventricle

Fractional Shortening can be measured on a single-frame or on a Dual 2D format, the latter display (in simultaneous mode) is recommended to display side-by-side a systolic and diastolic frame of the same cardiac cycle. The group parameters are measured in two phases.

Performing the Measurements

- Acquire a cardiac cycle and press **FREEZE**.
- If necessary, scroll memories to select the tele-diastolic image and press **MEASURE** to activate the menu.
- Select the “FAC” group and go to one of the calculations.
- Follow the instructions on the screen to trace the diastolic profile and the Ls long axis.
- Use the trackball to select the tele-systolic image and press **MEASURE**.
- Follow the instructions on the screen to complete the measurements.

If an entire cardiac cycle has not been stored, it is possible to return to Real Time with **B-MODE** to complete the acquisition. Freeze the image and press **MEASURE** to complete the measurement.

Fractional area change (FAC) is automatically obtained once the measurements have been taken:

Left Ventricle

LV

This measurement requires an entire cardiac cycle be acquired. The group includes the following:

Parameter		Abbreviation	Measurement
IVS dia	Intraventricular septum - Diastole	IVSd	Distance
LVD dia	LV diameter - Diastole	LVDd	Distance
PW dia	Posterior wall -Diastole	PW	Distance
LVD sys	LV Diameter-Systole	LVDs	Distance

LV: Left ventricle

The following parameter is automatically calculated once the measurements have been taken:

Parameter	
FS	Fractional shortening

Note*

The ejection fraction is calculated using the Teichholz formula for the M-Mode. The use of this formula in B-Mode is not supported by any bibliographic reference but it is based on widely clinical practice.

Left Ventricle Outflow Tract

OUTFLOW TRACT

This group includes the following measurements:

Parameter		Abbreviation	Measurement
LVOT DIAM	LV outflow tract diameter	LVOD	Distance
LVOT AREA	Outflow tract area	LVOA	Profile

LV: Left ventricle

The **area** (OTA) of the outflow tract is obtained once the measurement is completed.

Aorta

AORTA

This group includes the following measurements:

Parameter		Abbreviation	Measurement
AO DIAM	Aortic diameter	AOD	Distance
AO PLAN	Aortic planimetry	AVA	Profile
AO OPENIN	Aortic valve opening	AVO	Distance

The aortic **area** (AOA) is obtained once the measurement is completed.

LA

Left Atrium

This group includes the following measurement:

Parameter		Abbreviation	Measurement
LA DIAM	Left atrium diameter	LAD	Distance

RV

Right Ventricle

This group includes the following measurement:

Parameter		Abbreviation	Measurement
RVDd	RV diameter - Diastole	RVDd	Distance

RV: right ventricle

PA

Pulmonary Artery

This group includes the following measurement:

Parameter		Abbreviation	Measurement
PA DIAM	Pulmonary artery diameter	PAD	Distance

The following parameter is automatically calculated once the measurement has been taken:

Parameter	
PAA	Pulmonary artery area

Mitral

MITRAL

This group includes the following measurements:

Parameter		Abbreviation	Measurement
ANN DIAM	Mitral annulus diameter	MAN	Distance
ANN AREA	Mitral annulus area	MAA	Profile
MIT AREA	Mitral area	MVA	Profile

Cardiac Calculations in M-Mode

Left Ventricle

LV

The following parameters can be measured:

Parameter		Abbreviation	Measurement
RV DIAM d	RV Diameter-diastolic	RVDd	Distance
IVS dia	IV Septum-diastolic	IVSd	Distance
LV DIAM d	LV Diameter-diastolic	LVDd	Distance
PW dia	Post wall-diastolic	PWd	Distance
IVS sys	IV Septum-systolic	IVSs	Distance
LVD sys	LV Diameter.systolic	LVDs	Distance
PW sys	Post wall-systolic	PWs	Distance

LV: Left ventricle

RV: Right ventricle

The following parameters are automatically calculated once the measurements are complete:

Parameter	
EF	Ejection fraction
FS	Fractional shortening
LVVd	Diastolic volume
LVVs	Systolic volume
CO	Cardiac output
S%	Septum thickening
PW%	Post wall thickening
LVM	Left ventricle mass
LVMI	Left ventricle mass index

LV Mass

Left Ventricle Mass

The following parameters can be measured:

Parameter		Abbreviation	Measurement
IVS dia	IV Septum-diastolic	IVSd	Distance
LV DIAM d	LV Diameter-diastolic	LVDd	Distance
PW dia	Post wall-diastolic	PWd	Distance

LV: Left ventricle

The following parameters are automatically calculated once the measurements are complete:

Parameter	
LVM	Left ventricle mass
LVMI	Left ventricle mass index

EF

Ejection Fraction

The following parameters can be measured:

Parameter		Abbreviation	Measurement
LV DIAM d	LV Diameter-diastolic	LVDd	Distance
LVD sys	LV Diameter.systolic	LVDs	Distance

LV: Left ventricle

The following parameters are automatically calculated once the measurements are complete:

Parameter	
EF	Ejection fraction
FS	Fractional shortening

Aorta and Left Atrium

AORTA/LA

This group is composed of the following measurements:

Parameter		Abbreviation	Measurement
AO DIAM	Aortic diameter	AOD	Distance
LA	Left atrium	LA	Distance
AO OPENIN	Aortic valve opening	AVO	Distance
EJECT TIM	Ejection time	ET	Time
R-R INT	R-R Interval	R-R	Time
AO COAPT	AO Coaptation Line	AOC	Distance

Once the measurements are complete, the following parameters are calculated automatically:

Parameter

LA/A	Left atrium/aortic diameter
AOEI	AO Excentricity index

Mitral

MITRAL

The following parameters can be measured:

Parameter		Abbreviation	Measurement
E SEPTUM	E Septum	ESD	Distance
EF Slope	EF Slope	EFS	Slope

Cardiac Calculations in Doppler

Mitral

MITRAL

This group includes the following parameters:

Parameter		Abbreviation	Measurement
MIT FVI	Mitral flow profile	MFVI	Profile
PEAK V-E	Mitral Peak Vel E Wave	MEVp	Velocity
PEAK V-A	Mitral Peak Vel A Wave	MAVp	Velocity
PHT MIT	PHT	MPHT	Slope
E ACC TIM	Mitral E Wave Acc Time	MEAT	Time
E DEC TIM	Mitral E Wave Dec Time	MEDT	Time
ISOV REL	Mitral IsoV Relax Time	MIRT	Time

Once the measurements are complete the following parameters are automatically calculated:

Parameter	
MEGp	Mit Peak Grad (E)
MAGp	Mit Peak Grad (A)
MVmn	Mean Velocity
MGmn	Mean Gradient
MVA	Mitral Area
E/A	Mit peak Vel E Wave/Mit peak Vel A Wave

MITRAL AREA

This group includes the following parameter:

Parameter		Abbreviation	Measurement
MPHT	Mitral area	MPHT	Time

Once the measurements are complete the following parameter is automatically calculated:

Parameter	
MVA	Mitral Area

Mitral Regurgitation

MIT REGURGE

This group includes the following parameter:

Parameter	Abbreviation	Measurement
REG VEL	Regurgitation velocity	MVrg Velocity

Once the measurements are complete, the **regurgitation gradient** (MGrg) is calculated automatically.

Aorta

AORTA

This group includes the following parameters:

Parameter	Abbreviation	Measurement
AO FLOW	Aortic flow	AFVI Profile
AO VEL	Aortic velocity	AVp Velocity
DIA VEL	Aortic diastolic velocity	AVd
ACC TIME	Acceleration time	AAT Time
EJECT TIME	Ejection time	ET Time

Once the measurements are complete, the following parameters are calculated automatically:

Parameter	Measurement
AVmn	Mean velocity
AGmn	Mean gradient
AGp	Peak gradient

Aortic Area

AO EFF A

This measurement requires a CW trace for the aortic flow, a trace for the LVOT flow and B-Mode image for measuring the diameter. This group includes the following parameters:

Parameter	Abbreviation	Measurement
AO FLOW	Aortic Flow Profile	AFVI Profile
LVOT FLOW	LVOT flow profile	OFVI Profile

The group parameters are measured in three phases.

Performing the Measurements

Acquire the aortic flow and press **FREEZE**.

Press **MEASURE** to activate the menu, select the “Aortic Area Aortic Flow profile **ENTER**.

Follow the instructions on the screen for measuring the velocity integral.

Press **FREEZE** to return to Real Time and acquire the left ventricle outflow.

Press and then **MEASURE** to activate the menu.

Follow the instructions on the screen to measure the velocity integral of the left ventricle.

Press **FREEZE** and then **B-MODE** to return to Real Time and acquire the left ventricle image.

Press **FREEZE** and then **MEASURE** to reactivate the menu and to complete the measurements

Once the measurements are complete, the following parameters are calculated automatically:

Parameter	
OTA	Outflow tract area
AVA	Aortic effective valve area
LVOD	LVOT diameter

Aortic Regurgitation

AO REGURGE

This group includes the following parameters:

Parameter	Abbreviation	Measurement
REG PHT	AO regurge PHT	APHT Slope

Descending Aorta

DESC AORTA

This group includes the following parameters:

Parameter	Abbreviation	Measurement
SYS Vpeak	DA Systolic Peak Velocity	DAVp Velocity
PDA	Patent ductus artery	PDA Velocity

Once the measurements are complete, the **Peak gradient** (DAGp) is calculated automatically.

LVOT Flow

LVOT FLOW

This group includes the following parameters:

Parameter	Abbreviation	Measurement
LVOT FLOW	LVOT flow profile	OFVI Profile
LVOT VEL	LVOT peak velocity	OVp Velocity

Once the measurements are complete, the following parameters are calculated:

Parameter	Measurement
OVmn	Mean velocity
OGmn	Mean gradient
OGp	Peak gradient

Tricuspid

TRICUSPID

This group includes the following parameters:

Parameter	Abbreviation	Measurement
TRIC FVI	Tricuspid flow profile	TFVI Profile
PEAK V-E	Tricuspid Velocity E Wave	TEVp Velocity
PEAL V-A	Tricuspid Velocity A Wave	TAVp Velocity

Once the measurements are complete, the following parameters are calculated automatically:

Parameter	
TEGp	Tricuspid Peak Gradient (E)
TAGp	Tricuspid Peak Gradient (A)
TVmn	Tricuspid Mean Velocity
TGmn	Tricuspid Mean Gradient
E/A	Tricuspid Velocity E Wave/ Tricuspid Velocity A Wave

Tricuspid Regurgitation

TRIC REGURGE

This group includes the following parameter:

Parameter	Abbreviation	Measurement
REG VEL	Tric reg velocity	TRV
		Velocity

Once the measurements are complete, the following parameters are calculated:

Parameter	
TRG	Tricuspid Flow Gradient
RVPs	RV Systolic Pressure

Pulmonary Veins

PULM VEINS

This group includes the following parameters:

Parameter	Abbreviation	Measurement
PV SYS V	PV systolic velocity	PVVs
DIAST VEL	PV diastolic velocity	PVVD
REV ATR V	Reverse arterial velocity	RAV
		Velocity

Once the measurements are complete, the **SVDV** is calculated automatically.

Pulmonary Artery

PULM ARTERY

This group includes the following parameters:

Parameter	Abbreviation	Measurement
PULM FLOW	Pulm flow profile	PFVI
PULM Vp	Pulm peak velocity	PVp
		Profile
		Velocity

Once the measurements are complete, the following parameters are calculated automatically:

Parameter	
PVmn	Mean velocity
PGmn	Mean gradient
PGp	Peak gradient
PAP	Pulmonary artery pressure

Pulmonary Regurgitation

PULM REG

This group includes the following parameters:

Parameter		Abbreviation	Measurement
REG PHT	Pulm reg PHT	PPHT	Slope
PDIAS VEL	Pulm protodiast velocity	PVpd	Velocity
EDIAS VEL	Pulm enddiast velocity	PVed	Velocity

Cardiac Output - LVOT

CO-LVOT

This measurement requires the acquisition of a Doppler trace and a B-Mode image. This group includes the following parameters:

Parameter		Abbreviation	Measurement
LVOT FLOW	LVOT flow profile	OFVI	Profile
R-R INT	R-R interval	R-R	Time
LVOT DIAM	LVOT diameter	LVOD	Distance

The parameters of this group are measured in two phases.

Performing the Measurements

Acquire a Doppler trace and press **FREEZE**.

Press **MEASURE** to activate the menu, select “**CO-LVOT**”

Follow the instructions on the screen to measure the regurgitation profile and the R-R interval.

Press **FREEZE** to return to Real Time and go to **B-MODE**.

Acquire a B-Mode image and press **FREEZE**.

Display the required frame then press **MEASURE** to reactivate the menu.

Follow the instructions on the screen to complete the measurements.

The following parameters are calculated automatically:

Parameter	
HR	Heart rate
OTA	LVOT area
SV	Stroke volume
SI	Stroke index
CO	Cardiac output
CI	Cardiac index

Cardiac Output - Aorta

CO-AORTA

This measurement requires the acquisition of a Doppler trace and a B-Mode image. This group includes the following parameters:

Parameter		Abbreviation	Measurement
AO FLOW	Aortic flow profile	AFVI	Profile
R-R INT	R-R interval	R-R	Time
AO DIAM	AO diameter	AOD	Distance

The parameters of this group are measured in two phases.

Performing the Measurements

Acquire a Doppler trace and press **FREEZE**.

Press **MEASURE** to activate the menu, select the “CO-AORTA” group

Follow the instructions on the screen to measure the regurgitation profile and the R-R interval.

Press **FREEZE** to return to Real Time and go to **B-MODE**.

Acquire a B-Mode image and press **FREEZE**.

Display the required frame then press **MEASURE** to reactivate the menu.

Follow the instructions on the screen to complete the measurements.

Once the measurements are complete, the following parameters are calculated automatically:

Parameter	
HR	Heart rate
AOA	Aorta area
SV	Stroke volume
SI	Stroke index
CO	Cardiac output
CI	Cardiac index

Cardiac Output - Pulmonary

CO-PULM

This measurement requires the acquisition of a Doppler trace and a B-Mode image. This group includes the following parameters:

Parameter		Abbreviation	Measurement
PULM FLOW	Pulmonary flow profile	PFVI	Profile
R-R INT	R-R interval	R-R	Time
PA DIAM	Pulmonary diameter	PAD	Distance

The group parameters are measured in two phases.

Performing the Measurements

Acquire a Doppler trace and press **FREEZE**.

Press **MEASURE** to activate the menu, select the “CO-PULM” group

Follow the instructions on the screen to measure the regurgitation profile and the R-R interval.

Press **FREEZE** to return to Real Time and go to **B-MODE**.

Acquire a B-Mode image and press **FREEZE**.

Display the required frame then **MEASURE** to reactivate the menu.


Follow the instructions on the screen to complete the measurements.

Once the measurements are complete, the following parameters are calculated automatically:

Parameter	
HR	Heart rate
PAA	Pulmonary area
SV	Stroke volume
SI	Stroke index
CO	Cardiac output
CI	Cardiac index

Setup

Under **SETTINGS** tab press the **MENU** button and follow **MEASUREMENTS & ANALYSIS** – Cardiac to open the Cardiac configuration page:



This page allows the operator to:

- Include/exclude measurement pages per operating modes.
- Include/exclude single measurements for each operating mode.
- For the enabled measurements, it is possible to define how the parameters have to be measured:
 - Off disables the parameter calculation from the selected measurement.

MyLab – ADVANCED OPERATIONS

- On enables the parameter calculation in the selected measurement.
- Enable the average computation for the repeated measurements.

Cardiac Report

Pressing the **REPORT** tab on the control panel opens the Cardiac Report.

Cardiac Ultrasound Report

Date : 11-26-2010

Owner's name : Doctor's name :
Animal's name : Performing Physician :
Species : EQUINE Birthdate :
ID : 2010-11-26-05 Age : 0 Years
Accession number : Weight : 0.00 Kg
Diagnosis : Breed :
Gender : Unknown Neutered : Unaltered
BSA :

B-MODE

EF (AREA-LENGTH)

LV Diastolic Area : 31.38 cm ²	Ejection Fraction : -268.11 %
LV Diastolic Axis : 5.32 cm	Stroke Volume : 421.74 cc
LV Diastolic Volume : 157.30 cc	Stroke Index : --
LV Systolic Area : 76.25 cm ²	Heart Rate : -- bpm
LV Systolic Axis : 8.52 cm	Cardiac Output : -- l/min
LV Systolic Volume : 579.05 cc	Cardiac Index : --

LEFT VENTRICLE

IV Septum-Diastole : 8.53 cm	LV Diameter-Systole : -- mm
LV Diameter-Diastole : 5.95 cm	LV Fract Shortening : -- %
Post Wall-Diastole : 6.59 cm	

M-MODE

LEFT VENTRICLE

RV Diameter-Diastole : 1.93 cm	Ejection Fraction : -- %
IV Septum-Diastole : 7.87 cm	LV Fract Shortening : -- %
LV Diameter-Diastole : -- mm	Septum Thickening : -16.10 %
Post Wall-Diastole : 6.07 cm	PW Thickening : 16.48 %
IV Septum-Systole : 6.60 cm	LV Mass : -- g


MEASUREMENTS
COMMENTS
IMAGES
PRINT
PDF

The **COMMENTS** key allows the user to add comments at the end of the report. When pressing the comments key, the comments window appears allowing data to be edited in the empty field by use of the keyboard. Also predefined comment can be added to the report by checking the involved checkbox. These comments can be defined in the **REPORT** menu of the **DISPLAY** menu in the **SYSTEM SETUP** menu


The **MEASUREMENTS** key opens a window which displays all the measurements and calculations performed. Each measurement can be selected to be shown in the report. Different sections can be selected by the corresponding tabs.

The **IMAGES** key makes it possible to add pictures to the report that were taken during the exam. In the submenu **IMAGES** the pictures can be selected by ticking their corresponding boxes.

MyLab – ADVANCED OPERATIONS




To add the pictures to the report press **IMAGES** once more. The result can be seen in the picture below.



Examples of pages with different sections are reported hereinafter.

B-Mode

The B-Mode page displays all the measurements and calculations that have been performed in B-Mode.



In case of repeated measurements, single measurements can be unchecked or checked in order to exclude or include them in the mean value calculation (if the average calculation is enabled in the setup).

M-Mode

The M-Mode page displays all the measurements and calculations that have been performed in M-mode.

The screenshot displays the 'Doppler' mode interface. At the top, there are tabs for 'B-Mode', 'M-Mode', and 'Doppler'. A central table titled 'Left Ventricle' lists measurements and their values across three trials (1, 2, 3). Each value has a corresponding checkbox. Below this table, there are summary statistics for 'LEFT VENTRICLE' and 'M-MODE'. On the right side, there are buttons for 'MEASUREMENTS', 'COMMENTS', and 'PRINT'. On the left, there are fields for patient information such as Date, Owner's name, Animal's name, Species, ID, Accession number, Diagnosis, Gender, and BSA.

Label	Value	1	2	3
Left Ventricle				
RV Diameter-Diastole	1.93 cm	<input checked="" type="checkbox"/> 1.93 cm	<input type="checkbox"/> -- mm	<input type="checkbox"/> -- mm
IV Septum-Diastole	7.87 cm	<input checked="" type="checkbox"/> 6.53 cm	<input checked="" type="checkbox"/> 9.20 cm	<input type="checkbox"/> -- mm
LV Diameter-Diastole	-- mm	<input type="checkbox"/> -- mm	<input type="checkbox"/> -- mm	<input type="checkbox"/> -- mm
Post Wall-Diastole	5.07 cm	<input checked="" type="checkbox"/> 6.07 cm	<input type="checkbox"/> -- mm	<input type="checkbox"/> -- mm
IV Septum-Systole	5.60 cm	<input checked="" type="checkbox"/> 6.50 cm	<input type="checkbox"/> -- mm	<input type="checkbox"/> -- mm
LV Diameter-Systole	5.80 cm	<input checked="" type="checkbox"/> 6.80 cm	<input type="checkbox"/> -- mm	<input type="checkbox"/> -- mm
Post Wall-Systole	7.07 cm	<input checked="" type="checkbox"/> 7.07 cm	<input type="checkbox"/> -- mm	<input type="checkbox"/> -- mm

LEFT VENTRICLE			
IV Septum-Diastole	: 8.10 cm	IV Diameter-Systole	: -- mm
LV Diameter-Diastole	: 5.56 cm	LV Fract Shortening	: -- %
Post Wall-Diastole	: 5.59 cm		

M-MODE			
LEFT VENTRICLE			
RV Diameter-Diastole	: 1.93 cm	Ejection Fraction	: -- %
IV Septum-Diastole	: 7.87 cm	LV Fract Shortening	: -- %
LV Diameter-Diastole	: -- mm	Septum Thickening	: -16.10 %
Post Wall-Diastole	: 5.07 cm	PW Thickening	: 16.48 %
IV Septum-Systole	: 5.60 cm	LV Mass	: -- g

In case of repeated measurements, single measurements can be unchecked or checked in order to exclude or include them in the mean value calculation (if the average calculation is enabled in the setup).

Doppler

The Doppler page displays all the measurements and calculations that have been performed in Doppler mode.

MyLab – ADVANCED OPERATIONS

The screenshot displays the 'Doppler' mode of the software. A central window titled 'Mitral' shows a table of measurements:

Label	Value	1	2	3
Flow Profile	-- m	-- m	-- m	-- m
Peak Vel. E Wave	16.96 cm/s	16.96 cm/s	-- cm/s	-- cm/s
Peak Vel. A Wave	21.73 cm/s	21.73 cm/s	-- cm/s	-- cm/s
PHT	11.24 ms	11.24 ms	-- ms	-- ms
E Wave Acc. Time	63.22 ms	63.22 ms	-- ms	-- ms
E Wave Dec. Time	189.67 ms	189.37 ms	-- ms	-- ms
Isov. Relax. Time	105.37 ms	105.37 ms	-- ms	-- ms

Below this table, there are sections for 'LEFT VENTRICLE' and 'M-MODE LEFT VENTRICLE' with various measurements like 'IV Septum-Diastole' and 'Ejection Fraction'.

In case of repeated measurements, single measurements can be unchecked or checked in order to exclude or include them in the mean value calculation (if the average calculation is enabled in the setup).

Additional

The additional page displays all generic measurements which have been added to the report with the **ADD TO REPORT** key. In addition there is the option to rename the label of the measurement. To rename the label, select a measurement from the list and enter a new name in the “Edit Label” field and confirm with **enter**.

The screenshot shows the 'Additional' mode. An 'Edit Label' dialog box is open, displaying a table:

Label	Value	Label	Value
Distance			
D3	1.08 cm		

The background shows the 'Additional' mode interface with a sidebar for patient information and buttons for MEASUREMENTS, COMMENTS, and IMAGES.

Bibliographic References and Formulas in B-Mode

LV Volume Simpson Formula	Measure unit	Derived parameters
$LVVd \text{ and } LVVs = (\pi/4) * (h/20) * \sum_{1-20} d_h D_h$ <p>h: Long axis d_h: A2C diameter D_h: A4C diameter V: Volume</p> <p>Accuracy: $\pm 15\%$</p>	ml	-
<p>Schiller N.B., ..., Two-Dimensional Echocardiographic Determination of Ventricular Volume, Systolic Function and Mass. In: <i>Summary and Discussion of the 1989 Recommendations of the American Society of Echocardiography</i></p>		

LV Volume (A-L) Formula	Measure unit	Derived parameters
$LVVd \text{ and } LVVs = 8 * A^2 / (3 * \pi * D)$ <p>A: Area D: Long axis</p> <p>Accuracy: $\pm 21\%$</p>	ml	-
<p>Schiller N.B., ..., Two-Dimensional Echocardiographic Determination of Ventricular Volume, Systolic Function and Mass. In: <i>Summary and Discussion of the 1989 Recommendations of the American Society of Echocardiography</i></p>		

Ejection Fraction (Simpson and A-L) Formula	Measure unit	Derived parameters
$EF = (A-B)*100/A$ <p>A: LV diastolic volume B: LV systolic volume</p> <p>Accuracy: $\pm 42\%$</p>	-	-
<p>Feigenbaum H., Echocardiography, 4th Ed., Lea & Febiger, Philadelphia, 1986, pp. 153-155.</p>		

Stroke Volume Formula	Measure unit	Derived parameters
$SV = A - B$ <p>A: LV diastolic volume B: LV systolic volume</p> <p>Accuracy: $\pm 42\%$</p>	ml	-
<p>Weyman A., Principles and Practice of Echocardiography, Lea & Febiger, 1994, p. 605</p>		

MyLab – ADVANCED OPERATIONS

Stroke Index Formula	Measure unit	Derived parameters
SI = A/B	-	-
A: Stroke volume		
B: BSA		
Oh J, Seward J, Tajik A The echo manual-Second edition, Lippincott Williams &Wilkins		

Cardiac Output Formula	Measure unit	Derived parameters
CO = (A – B) * HR	l/min	-
A: LV diastolic volume		
B: LV systolic volume		
Accuracy: ± 45%		
Weyman A., Principles and Practice of Echocardiography, Lea & Febiger, 1994, p. 605		

Cardiac Index Formula	Measure unit	Derived parameters
CI = A/B	-	-
A: Cardiac output		
B: BSA		
Oh J, Seward J, Tajik A The echo manual-Second edition, Lippincott Williams &Wilkins		

Area Fractional Shortening Formula	Measure unit	Derived parameters
FAC = (A-B)*100/A	-	-
A: LV diastolic area		
B: LV systolic area		
Accuracy: ± 16%		

Diameter Fractional Shortening Formula	Measure unit	Derived parameters
FS = (A - B)*100/A	-	-
A: LV diastolic diameter		
B: LV systolic diameter		
Accuracy: ± 10%		
Quinones M.A., Gaasch W.H., Alexander J.K., "Echocardiographic Assessment of Left Ventricular Function with Special Reference to Normal Velocities". In: Circulation, 1974 , 50, p. 42.		

MyLab – ADVANCED OPERATIONS

Ejection Fraction (LV Group)	Measure unit	Derived parameters
Formula		
EF = (A - B)*100/A	-	A: [7 * D ³] / [2.4 + D] D: LV Diameter in diastole B: [7 * D ³] / [2.4 + D] D: LV Diameter in systole

LV Mass	Measure unit	Derived parameters
Formula		
LVM = 0.8 * {1.04[(A + B + C) ³ - A ³]} + 0.6	g	-
A: LV internal diameter – Diastole		
B: Post wall - Diastole		
C : Intraventricular septum - Diastole		
Accuracy: ± 15%		
Lang R, Bierig M, Devereux R,.. Recommendations for chamber quantification: a report from the American Society of Echocardiography’s Guidelines and Standards Committee and the Chamber Quantification Writing Group, developed in conjunction with the European Association of Echocardiography, a branch of the European Society of Cardiology In: J Amer. Soc. Echocardiography, 2005, Vol.18; N.12; pp1440-1463.		

Outflow Tract Area	Measure unit	Derived parameters
Formula		
OTA=Π * (D/2) ²	cm ²	-
D: Outflow tract diameter		

Aortic Valve Opening	Measure unit	Derived parameters
Formula		
AAO=Π * (D/2) ²	cm ²	-
D: Aortic diameter		

RV Volume	Measure unit	Derived parameters
Formula		
RVV = A * D ² /3	ml	-
A: RV Area		
D: RV long axis		
Accuracy: ± 21%		

Pulmonary Artery Area	Measure unit	Derived parameters
Formula		
PAA=Π * (D/2) ²	cm ²	-
D: Pulmonary Artery diameter		

LA Volume Formula	Measure unit	Derived parameters
LAV= (0.85 * A * B)/C	ml	-
A: Left atrium area – 4AC		
B: Left atrium area – 2AC		
C: Left atrium length		
Accuracy: ± 24%		
Oh J, Seward J, Tajik A The echo manual-Second edition, Lippincott Williams &Wilkins		

Bibliographic References and Formulas in M-Mode

LV Ejection Fraction Formula	Measure unit	Derived parameters
EF = (A-B)*100/A	-	A=(7 * D ³)/(2,4 * D) D: LV diastolic diameter (cm)
		B=(7 * D ³)/(2,4 * D) D: LV systolic diameter (cm)
Accuracy: ± 30%		
Teichholz L.E., ... , Problems in Echocardiographic Volume Determinations: Echocardiographic/Angiographic Correlations in the Presence or Absence of Asynergy. In: American Journal of Cardiology, 37, January 1976.1986, pp. 153-155.		

LV Fractional Shortening Formula	Measure unit	Derived parameters
FS = (A - B)*100/A	-	-
A: LV diastolic diameter		
B: LV systolic diameter		
Accuracy: ± 10%		
Feigenbaum H., Echocardiography, 4th Edition, Lea & Febiger, Philadelphia, 1986, pp. 153-155.		

Septum Thickening Formula	Measure unit	Derived parameters
S% = (A - B)*100/B	-	-
A: systolic interventricular septum		
B: diastolic interventricular septum		
Accuracy: ± 10%		
Feigenbaum H., Echocardiography, 4th Ed., Lea & Febiger, Philadelphia, 1986, pp. 153-155.		

MyLab – ADVANCED OPERATIONS

Posterior Wall Thickening Formula	Measure unit	Derived parameters
$PW\% = (A - B) * 100 / B$	-	-
A: systolic posterior wall B: diastolic posterior wall		
Accuracy: $\pm 10\%$		
Feigenbaum H., Echocardiography, 4th Ed., Lea & Febiger, Philadelphia, 1986, pp. 153-155.		

LV Mass Formula	Measure unit	Derived parameters
$LVM = 1.04[(A+B+C)^3 - B^3] - 13.6$	g	-
A: diastolic interventricular septum B: LV diastolic diameter C: diastolic posterior wall		
Accuracy: $\pm 15\%$		
Devereux R.B., Reichek N., ..., Echocardiographic Determination of Left Ventricular Mass in Man - Anatomic Validation of the Method. In: Circulation, n.55, 1977, pp. 613-8		

LV Mass Index Formula	Measure unit	Derived parameters
$LVM I = LVM / BSA$	-	-
LVM: LV Mass BSA: Body Surface Area		
Accuracy: $\pm 15\%$		
Devereux R.B., Reichek N., ..., Echocardiographic Determination of Left Ventricular Mass in Man - Anatomic Validation of the Method. In: Circulation, n.55, 1977, pp. 613-8		

LA/Aorta Diameter Formula	Measure unit	Derived parameters
$LA/AD = A/B$	-	-
A: Left Atrium B: Aortic diameter		
Accuracy: $\pm 10\%$		

Exentricity Index Formula	Measure unit	Derived parameters
$EXC IND = A/B$	-	-

MyLab – ADVANCED OPERATIONS

- A: Aortic diameter
- B: Aortic coaptation line

Accuracy: ±10%

Nanda N.C., Gramiak R.: Evaluation of Bicuspid Valves by Two-Dimensional Echocardiography. In: American J. Cardiol. 1987, 11 p.372

Bibliographic References and Formulas in Doppler

Peak Gradient Formula	Measure unit	Derived parameters
$G_p = 4 \cdot V^2$	mmHg	-
V: Peak velocity		
Accuracy: ± 16%		
Weyman A., Principles and Practice of Echocardiography, Lea & Febiger, 1994, p. 605		

Mean Velocity Formula	Measure unit	Derived parameters
$V_{mn} = FVI/t$	m/s	t=flow duration
Accuracy: ± 11%		

Mean Gradient Formula	Measure unit	Derived parameters
$G_{mn} = 4 * (V_1^2 + V_2^2 + \dots + V_n^2)/n$	mmHg	-
Vi: Instant velocities		
Accuracy: ± 11%		
Weyman A., Principles and Practice of Echocardiography, Lea & Febiger, 1994, p. 605		

PHT Formula	Measure unit	Derived parameters
$PHT = V_{MAX} * (1 - 0.707)/Slope$	ms	-
Accuracy: ± 28%		
Hatle L., Angelsen B., Noninvasive Assessment of Atrioventricular Pressure Half-Time by Doppler Ultrasound. In: Circulation 60, n.5, 1979, pp. 1096-1104		

Mitral Valve Area Formula	Measure unit	Derived parameters
$MVA = 220/PHT$	cm ²	-
Accuracy: ± 28%		
Weyman A., Principles and Practice of Echocardiography, Lea & Febiger, 1994, p. 605		

E Wave/A Wave Formula	Measure unit	Derived parameters
------------------------------	---------------------	---------------------------

MyLab – ADVANCED OPERATIONS

$$E/A=A/B$$

A: E wave peak velocity

B: A wave peak velocity

Accuracy: $\pm 10\%$

Miocardiac Performance Index Formula	Measure unit	Derived parameters
---	---------------------	---------------------------

$$MPI = (A+B)/C$$

A: Isovolumetric Contraction Time

B: Isovolumetric Relax Time

C: Ejection Time

Accuracy: $\pm 6\%$

C.Bruch, A.Schmermund ..“TEI-index in patients with mid-to-moderate congestive heart failure”, In: Eu. H.J. 2000, n.21 pp.1888-1895

dP/dt Ratio Formula	Measure unit	Derived parameters
--------------------------------	---------------------	---------------------------

$$dP/dt = 32/t$$

mmHg/s

-

t: time elapsed between $-1m/s$ to $-3m/s$ velocity values

Accuracy: $\pm 3\%$

Bargiggia GS, Bertucci C,..., Anew method for estimating left ventricular dP/dt by continuous wave Doppler echocardiography. Validation studies at cardiac catheterisation, In: Circulation 1989; 80; 1287-1292

Regurgitation Flow (PISA) Formula	Measure unit	Derived parameters
--	---------------------	---------------------------

$$REG = 6.28 * R^2 * V$$

ml/s

-

R: Radius

V: Aliasing velocity

Accuracy: $\pm 14\%$

Bargiggia G.S., Tronconi L., Sahn D.J., A New Method for Quantitation of Mitral Regurgitation Based on Color Flow Doppler Imaging of Flow Convergence Proximal to Regurgitant Orifice. In: Circulation, 1991, 84: pp. 1481-1489

Effective Regurgitation Orefice (PISA) Formula	Measure unit	Derived parameters
---	---------------------	---------------------------

$$O = 6.28 * R^2 * V1 / V2$$

ml

-

R: Radius

MyLab – ADVANCED OPERATIONS

V1: Aliasing velocity

V2: Regurgitation peak velocity

Accuracy: ± 22%

Oh J, Seward J, Tajik A The echo manual-Second edition, Lippincott Williams &Wilkins

Mitral Regurgitation Volume (PISA) Formula	Measure unit	Derived parameters
$MVOL = 6.28 * R^2 * V / 3.25$	ml	-
R: Radius		
V: Aliasing velocity		
Accuracy: ± 14%		
Rossi A., Dujardin K.S., ..., Rapid Estimation of Regurgitant Volume by the Proximal Isovelocivity Surface Area Method in Mitral Regurgitation: Can Continuous-Wave Doppler Echocardiography Be Omitted? In: Journal of the American Society of Echocardiography. Volume 11, Number 2, pp. 138-148.		

Aortic Regurgitation Volume (PISA) Formula	Measure unit	Derived parameters
$REG VOL = 6.28 * R^2 * V1 * FVI / V2$	ml	-
R: R Radius		
V1: aliasing velocity		
FVI: flow velocity integral		
V2: Regurge peak velocity		
Accuracy: ± 30%		
Shiota T., Jones M., Yamada I., ..., Effective Regurgitant orifice Area by the Color Doppler Flow Convergence Method for Evaluating the Severity of Chronic Aortic Regurgitation. An Animal Study. In: Circulation, 1996; 93; pp. 594-602.		

E' Wave/A' Wave Formula	Measure unit	Derived parameters
$E'/A' = A/B$	-	-
A: Peak velocity E' wave		
B: Peak velocity A' wave		
Accuracy: ± 16%		

E Wave/E' Wave Formula	Measure unit	Derived parameters
$E/E' = A/B$	-	-
A: Peak velocity E wave		
B: Peak velocity E' wave		

MyLab – ADVANCED OPERATIONS

Accuracy: $\pm 16\%$

Intraventricular Mechanical Delay Formula	Measure unit	Derived parameters
--	---------------------	---------------------------

IMD = A-B	ms	-
-----------	----	---

A: Aorta Pre-ejection time

B: Pulmonary Pre-ejection time

Accuracy: $\pm 9\%$

F.Knebel, R.K.Reibeis., "Tissue Doppler Echocardiography and Biventricular Pacing Heart Failure: Patient Selection, Procedural Guidance, Follow up, quantification of Success", IN: Card Ultr 2004, n.2-17

Effective Aortic Valve Area Formula	Measure unit	Derived parameters
--	---------------------	---------------------------

AVA = $A \cdot FVI1 / FVI2$	cm ²	$A = \pi \cdot (D/2)^2$
-----------------------------	-----------------	-------------------------

A: LVOT area

D: LVOT diameter

FVI1: LVOT flow velocity integral

FVI2: aortic tract flow velocity integral

Accuracy: $\pm 28\%$

Huntsman L., Stewart D., ..., Noninvasive Doppler Determination of Cardiac Output in Man, In: Circulation 67, n. 3, March 1983

Maximal Aortic Valve Area Formula	Measure unit	Derived parameters
--	---------------------	---------------------------

AMVA = $A \cdot V1 / V2$	cm ²	$A = \pi \cdot (D/2)^2$
--------------------------	-----------------	-------------------------

A: LVOT area

D: LVOT diameter

V1: Aortic Peak velocity in LVOT

V2: Aortic Peak velocity

Accuracy: $\pm 22\%$

Zaghbi WA, Farmer KL, ..., Accurate non-invasive quantification of stenotic aortic valve area by Doppler echocardiography, In:Circulation 1986; 73; 452-459

Systolic Pressure Formula	Measure unit	Derived parameters
--------------------------------------	---------------------	---------------------------

$P_s = 4 \cdot V^2 + 10$	mmHg	-
--------------------------	------	---

V: regurge velocity

MyLab – ADVANCED OPERATIONS

Accuracy: $\pm 16\%$

Currie P.J., ..., Continuous Wave Doppler Determination of Left Ventricular Pressure: a Simultaneous Doppler Catheterization Study in 127 Patients. In: J. Amer. College Cardiol. 1985, 6, p.750

Systolic Velocity/Diastolic Velocity Formula	Measure unit	Derived parameters
VS/VD = A/B	-	-
A: systolic velocity		
B: diastolic velocity		
Accuracy: $\pm 10\%$		

Heart Rate Formula	Measure unit	Derived parameters
HR = 60/T	bpm	-
T: R-R interval		
Accuracy: $\pm 3\%$		

Stroke Volume Formula	Measure unit	Derived parameters
SV = A * FVI	ml	$A = \Pi * (D/2)^2$
A: Area		D: Diameter
FVI: Flow velocity		
Accuracy: $\pm 19\%$		

Huntsman L., Stewart D., ..., Noninvasive Doppler Determination of Cardiac Output in Man, In: Circulation 67, n. 3, March 1983

Stroke Index Formula	Measure unit	Derived parameters
SI = A/B	-	-
A: Stroke Volume		
B: BSA		
Accuracy: $\pm 19\%$		

Huntsman L., Stewart D., ..., Noninvasive Doppler Determination of Cardiac Output in Man, In: Circulation 67, n. 3, March 1983; Skjaerpe T, Hegrenaes L., ..., Non invasive estimation of valve area in patients with aortic stenosis by Doppler ultrasound and two-dimensional echocardiography, In: Circulation 1985; 72; 810-818

Cardiac Output Formula	Measure unit	Derived parameters
CO = A*FVI*HR	l/min	$A = \Pi * (D/2)^2$
A: Area		D: Diameter

MyLab – ADVANCED OPERATIONS

FVI: Flow velocity integral

HR: heart rate

Accuracy: $\pm 21\%$

Huntsman L., Stewart D., ..., Noninvasive Doppler Determination of Cardiac Output in Man, In: Circulation 67, n. 3, March 1983

Cardiac Index Formula	Measure unit	Derived parameters
----------------------------------	---------------------	---------------------------

CI = A/B

A: Cardiac Output

B: BSA

Accuracy: $\pm 19\%$

Huntsman L., Stewart D., ..., Noninvasive Doppler Determination of Cardiac Output in Man, In: Circulation 67, n. 3, March 1983; Skjaerpe T, Hegrenaes L....., Non invasive estimation of valve area in patients with aortic stenosis by Doppler ultrasound and two-dimensional echocardiography, In: Circulation 1985; 72; 810-818

Qp/Qs Ratio Formula	Measure unit	Derived parameters
--------------------------------	---------------------	---------------------------

Qp/Qs = A/B

A: Pulmonary artery stroke volume

B: LVOT stroke volume

Accuracy: $\pm 42\%$

Sanders S.P., ..., Measurement of Systemic and Pulmonary Blood Flow and Qp/Qs Ratio using Doppler and Two-Dimensional Echocardiography. In: Am. J. Cardiol. 1983, 51, p.952

Formulas and Bibliographic References for Body Surface Area in Canine and Feline Cardio

In the Canine Cardio and Feline Cardio applications, the system automatically calculates the Body Surface Area (BSA), based on the animal weight and the following tables:

Weight to Body Surface Area Conversion Chart – Canine (Dog)									
kg	m ²	kg	m ²	kg	m ²	kg	m ²	kg	m ²
0.5	0.064	10.0	0.469	20.0	0.744	30.0	0.975	40.0	1.181
1.0	0.101	11.0	0.500	21.0	0.759	31.0	0.997	41.0	1.201
2.0	0.160	12.0	0.529	22.0	0.785	32.0	1.018	42.0	1.220
3.0	0.210	13.0	0.553	23.0	0.817	33.0	1.029	43.0	1.240
4.0	0.255	14.0	0.581	24.0	0.840	34.0	1.060	44.0	1.259
5.0	0.295	15.0	0.608	25.0	0.864	35.0	1.081	45.0	1.278
6.0	0.333	16.0	0.641	26.0	0.886	36.0	1.101	46.0	1.297
7.0	0.370	17.0	0.668	27.0	0.909	37.0	1.121	47.0	1.302
8.0	0.404	18.0	0.694	28.0	0.931	38.0	1.142	48.0	1.334
9.0	0.437	19.0	0.719	29.0	0.953	39.0	1.162	49.0	1.352
								50.0	1.371
Weight to Body Surface Area Conversion Chart – Feline (Cat)									
kg	m ²	kg	m ²	kg	m ²	kg	m ²	kg	m ²
0.1	0.022	1.4	0.125	3.6	0.235	5.8	0.323	8.0	0.400
0.2	0.034	1.6	0.137	3.8	0.244	6.0	0.330	8.2	0.407
0.3	0.045	1.8	0.148	4.0	0.252	6.2	0.337	8.4	0.413
0.4	0.054	2.0	0.159	4.2	0.260	6.4	0.345	8.6	0.420
0.5	0.063	2.2	0.169	4.4	0.269	6.6	0.352	8.8	0.426
0.6	0.071	2.4	0.179	4.6	0.277	6.8	0.360	9.0	0.433
0.7	0.079	2.6	0.189	4.8	0.285	7.0	0.366	9.2	0.439
0.8	0.086	2.8	0.199	5.0	0.292	7.2	0.373	9.4	0.445
0.9	0.093	3.0	0.208	5.2	0.300	7.4	0.380	9.6	0.452
1.0	0.100	3.2	0.217	5.4	0.307	7.6	0.387	9.8	0.458
1.2	0.113	3.4	0.226	5.6	0.315	7.8	0.393	10.0	0.464

Reference: Withrow SJ, MacEwen, EG: *Small Animal Clinical Oncology, Second Edition, 1995.* WB Saunders Company, Philadelphia, PA, USA

3 - Abdominal Measurements Package



Equine



Canine



Feline

Abdominal application icon

This chapter describes the Abdominal measurements package (measurements and report environments) and lists the available measurements.

The Abdominal application package applies for the equine, canine and feline species.

Abdominal

ABDO

Note

This application calculation package is only available when the corresponding license is activated. Please refer to the System Configuration manual for all the options/licenses of the **MyLab** and their status.

Abdominal Application

The Abdominal application can be activated selecting the **ABDO** key from the application list.

When the Abdominal application is selected, in freeze mode, press the **MEASURE** key to activate the Abdominal measurements.

Abdominal Measurements

B-Mode Measurements

In each application this package allows distance measurements and volume calculations of up to 14 different organs (groups). Measurements and volume calculations are also provided for certain abnormalities (cysts, masses and calculi) For bilateral organs, a software toggle is available for selecting the right or left organ. The organ groups are detailed further:

KIDNEY

KIDNEY Group

This package allows the following measurements:

Subgroup	Parameter		Measure
SAGITTAL LENGTH	Sagittal Length	Distance	Distance
SAGITTAL HEIGHT	Sagittal Height	Distance	Distance
SAGITTAL AREA	Sagittal Area	Area	Area
TRANSVERSE HEIGHT	Transverse Height	Distance	Distance
TRANSVERSE WIDTH	Transverse Height	Distance	Distance
TRANSVERSE AREA	Transverse Area	Area	Area
VOLUME	Volume	Volume	Volume
CORTEX	Cortex	Distance	Distance
MEDULLA	Medulla	Distance	Distance
CORTEX/ MEDULLA RATIO	Cortex/ Medulla ratio	%Distance	%Distance
KIDNEY CALCULI DIAMETER	KIDNEY CALCULI Diameter	Distance	Distance
KIDNEY MASS1	KIDNEY MASS1	Complex volume	Complex volume

URINARY Group

This package allows the following measurements:

URINARY

Subgroup	Parameter		Measure
Ureter Diameter	Ureter Diameter	Distance	Distance
Urinary Bladder Sagittal Length	Urinary Bladder Sagittal Length	Distance	Distance
Urinary Bladder Sagittal Height	Urinary Bladder Sagittal Height	Distance	Distance
Urinary Bladder Sagittal Area	Urinary Bladder Sagittal Area	Area	Area
Urinary Bladder Wall Thickness	Urinary Bladder Wall Thickness	Distance	Distance
Urinary Bladder Mass	Urinary Bladder Mass	Complex volume	Complex volume

ADRENAL GLANDS Group

This package allows the following measurements:

ADRENAL GLANDS

Subgroup	Parameter		Measure
Sagittal Length	Sagittal Length	Distance	Distance
Sagittal Area	Sagittal Area	Area	Area
ADRENAL GLAND MASS	ADRENAL GLAND MASS	Complex volume	Complex volume

PANCREAS Group

This package allows the following measurements:

PANCREAS

Subgroup	Parameter		Measure
Body thickness	Body thickness	Distance	Distance
Right lobe thickness	Right lobe thickness	Distance	Distance
Left lobe thickness	Left lobe thickness	Distance	Distance

SPLEEN Group

This package allows the following measurements:

SPLEEN

Subgroup	Parameter		Label	Measure
Head thickness	Head thickness	Distance	SPL HEAD TH	Distance
Body thickness	Body thickness	Distance	SP BODY TH	Distance
Tail Thickness	Tail Thickness	Distance	SPL TAIL TH	Distance

GASTROINTESTINAL Group

This package allows the following measurements:

**GASTRO
INTESTINAL**

Subgroup	Parameter		Label	Measure
STOMACH BODY Thickness	STOMACH BODY Thickness	Distance	ST BODY TH	Distance
Duodenum Thickness	Duodenum Thickness	Distance	DUODENUM TH	Distance
Jejunum Thickness	Jejunum Thickness	Distance	JEJUNUM TH	Distance
ILEUM Thickness	ILEUM Thickness	Distance	ILEUM	Distance
COLON Thickness	COLON Thickness	Distance	COLON TH	Distance

LIVER Group

This package allows the following measurements:

LIVER

Subgroup	Parameter		Measure
Longitudinal distance	Longitudinal distance	Distance	Distance
Transverse distance	Transverse distance	Distance	Distance
LIVER MASS	LIVER MASS	Complex volume	Complex volume

GALLBLADDER Group

This package allows the following measurements:

GALLBLADDER

Subgroup	Parameter		Label	Measure
Sagittal Length	Sagittal Length	Distance	SAG L	Distance
Sagittal Height	Sagittal Height	Distance	SAG H	Distance
Sagittal Area	Sagittal Area	Area	SAG A	Area
Transverse Height	Transverse Height	Distance	TRANSV H	Distance
Transverse Width	Transverse Width	Distance	TRANSV W	Distance
Transverse Area	Transverse Area	Area	TRANSV A	Area
Volume	Volume	Volume	VOLUME	Volume
			Only shown in worksheet	
GALLBLADDER WALL Thickness	GALLBLADDER WALL Thickness	Distance	GALLB W TH	Distance
Common Bile Duct	Common Bile Duct	Distance	C BILE D	Distance
GALLBLADDER CALCULI	GALLBLADDER CALCULI	Distance	GLLB CAL1	Distance

PROSTATE Group

This package allows the following measurements:

PROSTATE

Subgroup	Parameter		Label	Measure
Sagittal Length	Sagittal Length	Distance	SAG L	Distance
Sagittal Height	Sagittal Height	Distance	SAG H	Distance
Sagittal Area	Sagittal Area	Area	SAG A	Area
Transverse Height	Transverse Height	Distance	TRANSV H	Distance
Transverse Width	Transverse Width	Distance	TRANSV W	Distance
Transverse Area	Transverse Area	Area	TRANSV A	Area
Volume	Volume	Volume	VOLUME	Volume
			Only shown in worksheet	

UTERUS Group

This package allows the following measurements:

	Subgroup	Parameter		Label	Measure
UTERUS	Diameter	Diameter	Distance	DIAMETER	Distance
	Wall Thickness	Wall Thickness	Distance	WALL	Distance

OVARY Group

This package allows the following measurements:

	Subgroup	Parameter		Label	Measure
OVARY	Sagittal Length	Sagittal Length	Distance	SAGIT L	Distance
	Sagittal Height	Sagittal Height	Distance	SAGIT H	Distance
	Sagittal Area	Sagittal Area	Area	SAGIT A	Area
	Transverse Height	Transverse Height	Distance	TRANSV H	Distance
	Transverse Width	Transverse Width	Distance	TRANSV W	Distance
	Transverse Area	Transverse Area	Area	TRANSV A	Area
	Volume	Volume	Volume	VOLUME	Volume
				Only shown in worksheet	

CAUDAL VENA CAVA Group

This package allows the following measurements:

	Subgroup	Parameter		Label	Measure
CAUDAL VENA CAVA	Transverse Diameter	Transverse Diameter	Distance	TRANSV D	Distance
	Transverse Area	Transverse Area	Area	TRANSV A	Area

PORTAL VEIN Group

This package allows the following measurements:

	Subgroup	Parameter		Label	Measure
PORTAL VEIN	Transverse Diameter	Transverse Diameter	Distance	TRANSV D	Distance
	Transverse Area	Transverse Area	Area	TRANSV A	Area

Doppler Measurements

In each application this package allows velocity and resistance index measurements of up to 4 different organs (groups). The organ groups are detailed further

HEPATIC VEIN Group

This package allows the following measurements:

HEPATIC VEIN	Subgroup	Parameter	Measure	
	SYSTOLIC VELOCITY	Systolic velocity	Velocity	Velocity
DIASTOLIC VELOCITY	Diastolic velocity	Velocity	Velocity	

KIDNEY Group

This package allows the following measurements:

KIDNEY	Subgroup	Parameter	Measure	
	SYSTOLIC VELOCITY	Systolic velocity	Velocity	Velocity
DIASTOLIC VELOCITY	Diastolic velocity	Velocity	Velocity	
RESISTANCE INDEX (R.I.)	RESISTANCE INDEX	Dedicated formula	Resistance Index	

PORTAL VEIN Group

This package allows the following measurements:

PORTAL VEIN	Subgroup	Parameter	Measure	
	SYSTOLIC VELOCITY	Systolic velocity	Velocity	Velocity
DIASTOLIC VELOCITY	Diastolic velocity	Velocity	Velocity	


HEPATIC ARTERY Group

This package allows the following measurements:

HEPATIC ARTERY	Subgroup	Parameter	Measure	
	SYSTOLIC VELOCITY	Systolic velocity	Velocity	Velocity
DIASTOLIC VELOCITY	Diastolic velocity	Velocity	Velocity	
RESISTANCE INDEX (R.I.)	RESISTANCE INDEX	Dedicated formula	Resistance Index	

Setup

The Abdominal measurement package can be configured via the **SETTINGS** tab – **MENU - MEASUREMENTS & ANALYSIS – ABDOMINAL**. The Setup environment consists of the following page:



This page permits the user to include or exclude pages and, inside each page, measurements.

For the enabled measurements, it is possible to disable the sub-parameters or enable them and define how to measure:

- When the parameter is set to On it can be measured only if selected (not via the folder related to the main measurement).
- When the parameter is set to Off, the parameter is excluded from the measurement.

In the same configuration page it is possible to enable the average for the measurements.

Abdominal Report

Pressing the **REPORT** tab on the control panel opens the Abdominal Report.

Abdominal Ultrasound Report

Date : 11-29-2010

Owner's name : Doctor's name :
 Animal's name : Performing Physician :
 Species : EQUINE Birthdate :
 ID : 2010-11-29-00 Age : 0 Years
 Accession number : Weight : 0.00 Kg
 Diagnosis : Breed :
 Gender : Unknown Neutered : Unaltered

B-MODE

LEFT KIDNEY TRANSVERSE

Height : 2.63 cm Area : 5.94 cm²
 Width : 2.07 cm

RIGHT KIDNEY SAGITTAL

Length : 2.02 cm Area : 4.64 cm²
 Height : 2.29 cm

RIGHT KIDNEY TRANSVERSE

Height : 2.64 cm Area : 5.12 cm²
 Width : 1.87 cm

RIGHT KIDNEY VOLUME

Volume : 8.66 cc

RIGHT KIDNEY MASS

Distance 1 : 1.78 cm Distance 3 : 1.45 cm
 Distance 2 : 1.64 cm Volume : 2.23 cc


PANCREAS BODY

MEASUREMENTS
 COMMENTS
 IMAGES
 PRINT
 PDF


The **COMMENTS** key allows the user to add comments at the end of the report. When pressing the comments key, the comments window appears allowing data to be edited in the empty field by use of the keyboard. Also predefined comment can be added to the report by checking the involved checkbox. These comments can be defined in the **REPORT** menu of the **DISPLAY** menu in the **SYSTEM SETUP** menu

The **MEASUREMENTS** key opens a window which displays all the measurements and calculations performed. Each measurement can be selected to be shown in the report. Different sections can be selected by the corresponding tabs.

The **IMAGES** key makes it possible to add pictures to the report that were taken during the exam. In the submenu **IMAGES** the pictures can be selected by ticking their corresponding boxes.



To add the pictures to the report press **IMAGES** once more. The result can be seen in the picture below.



Examples of pages with different sections are reported hereinafter.

B-Mode

The B-Mode page displays all the measurements and calculations that have been performed in B-Mode.

In case of repeated measurements, single measurements can be unchecked or checked in order to exclude or include them in the mean value calculation (if the average calculation is enabled in the setup).

The screenshot displays the B-Mode software interface. At the top, there are tabs for 'B-Mode' and 'Doppler'. A central table lists measurements for the 'Left Kidney Transverse' and 'Right Kidney Sagittal' views. Each measurement has a 'Value' column and three columns for individual measurements (1, 2, 3), each with a checkbox and a unit. Below the table, there are sections for 'RIGHT KIDNEY TRANSVERSE', 'RIGHT KIDNEY VOLUME', 'RIGHT KIDNEY MASS', and 'PANCREAS BODY'. On the right side, there are buttons for 'MEASUREMENTS', 'COMMENTS', 'IMAGES', 'PRINT', and 'PDF'. On the left side, there are input fields for 'Date', 'Owner's name', 'Animal's name', 'Species', 'ID', 'Accession num', 'Diagnosis', and 'Gender'.

	Label	Value	1	2	3
Left Kidney Transverse					
	Height	2.63 cm	<input checked="" type="checkbox"/> 2.08 cm	<input checked="" type="checkbox"/> 3.19 cm	<input type="checkbox"/> ~mm
	Width	2.07 cm	<input checked="" type="checkbox"/> 2.08 cm	<input checked="" type="checkbox"/> 2.06 cm	<input type="checkbox"/> ~mm
	Area	5.94 cm ²	<input checked="" type="checkbox"/> 5.48 cm ²	<input checked="" type="checkbox"/> 6.43 cm ²	<input type="checkbox"/> ~mm ²
Right Kidney Sagittal					
	Length	1.77 cm	<input checked="" type="checkbox"/> 1.77 cm	<input checked="" type="checkbox"/> 2.09 cm	<input checked="" type="checkbox"/> 2.19 cm
	Height	1.90 cm	<input checked="" type="checkbox"/> 1.90 cm	<input checked="" type="checkbox"/> 2.79 cm	<input checked="" type="checkbox"/> 2.19 cm
	Area	4.64 cm ²	<input checked="" type="checkbox"/> 5.53 cm ²	<input checked="" type="checkbox"/> 3.47 cm ²	<input checked="" type="checkbox"/> 4.92 cm ²

RIGHT KIDNEY TRANSVERSE

Height : 2.64 cm Area : 5.12 cm²
 Width : 1.87 cm

RIGHT KIDNEY VOLUME

Volume : 8.66 cc


RIGHT KIDNEY MASS

Distance 1 : 1.78 cm Distance 3 : 1.45 cm
 Distance 2 : 1.64 cm Volume : 2.23 cc

PANCREAS BODY

Doppler


The Doppler page displays all the measurements and calculations that have been performed in Doppler mode.



In case of repeated measurements, single measurements can be unchecked or checked in order to exclude or include them in the mean value calculation (if the average calculation is enabled in the setup).

Additional

The additional page displays all generic measurements which have been added to the report with the **ADD TO REPORT** key. In addition there is the option to rename the label of the measurement. To rename the label, select a measurement from the list and enter a new name in the “Edit Label” field and confirm with **enter**.



4 - Equine MSK / Tendon Package

This chapter describes the Equine MSK / Tendon application package (measurements and report environments) and lists the available measurements.



Equine species icon

Equine Tendon



Note

This application calculation package is only available when the corresponding license is activated. Please refer to the System Configuration manual for all the options/licenses of the **MyLab** and their status.

Equine MSK/Tendon Application

The Equine MSK/Tendon application can be activated selecting the **MSK TENDON** key from the application list, within the Equine species selected.

When the Equine MSK/Tendon application is selected, in freeze mode, press the **MEASURE** key to activate the Tendon measurements.

Equine Tendon Measurements

B-Mode Measurements

Pressing the **MEASURE** button makes the list of the equine tendon measurements appear on the right side of the screen.

The structure of the measurement package is based on a “core measurement“

MyLab – ADVANCED OPERATIONS

This consists of a repetitive pattern of measurements of distance, area and distance or area ratio and distance. Two alternative measurement methods can be used:

- distance, where the position of the measurement is indicated in distance to a specific reference point;
- zone, where the measurement position is indicated traditionally, in the classical manner of dividing the limbs into zones.

Measurement items list

A selection of the limb, structure and leg zone is performed prior to starting the measurement.

This is done by pressing the **MEASURE** button, then the dedicated **LIGAMENTS/ TENDON** software key.

A dialog window appears, as indicated below:

Measurement > Ligament / Tendon

Fore Hind

Left Right

Anatomic Structure Selection

Region	Structure
<input checked="" type="checkbox"/> Carpus	<input checked="" type="checkbox"/> Acc Lig DDF
<input type="checkbox"/> Fetlock	<input type="checkbox"/> Acc Lig SDF
<input type="checkbox"/> Palmar	<input type="checkbox"/> CDE tendon
<input type="checkbox"/> Pastern	<input type="checkbox"/> CDE tendon sheath
<input type="checkbox"/> Shoulder	<input type="checkbox"/> Carpal Sheath
<input type="checkbox"/> Spine	<input type="checkbox"/> Dorsal Carpal Sheath
	<input type="checkbox"/> DDFT
	<input type="checkbox"/> ECO tendon
	<input type="checkbox"/> ECO tendon sheath
	<input type="checkbox"/> ECR tendon

Ok Cancel

The user can select the specific limb, the region and the structure upon which the measurement should be performed. After that, the measurements can be performed sequentially or selectively.

A complete list of regions and structures is given below:

Palmar Plantar MC/ MT	Parameter	Selectable in	Label	Measure
	SDF (superficial digital flexor) tendon	Front, Hind	SDFT	As selected in Method
DDF (deep digital flexor) tendon	Front, Hind	DDFT	As selected in Method	
Inferior Check Ligament	Front, Hind	ICL	As selected in Method	
Accessory Ligament SDF	Front, Hind	AccLigSDF	As selected in Method	
Accessory Ligament DDF	Front, Hind	AccLigDDF	As selected in Method	
SL (Suspensory Ligament) Body	Front, Hind	SLBody	As selected in Method	
SL (Suspensory Ligament) Lateral Branch	Front, Hind	SLLatBr	As selected in Method	
SL (Suspensory Ligament) Medial Branch	Front, Hind	SLMedBr	As selected in Method	
Superior Check Ligament	Front, Hind	SupChkLig	As selected in Method	

Fetlock	Parameter	Selectable in	Label	Measure
	SDFT – lateral	Front, Hind	LatBrSDFT	As selected in Method
SDFT – medial	Front, Hind	MedBrSDFT	As selected in Method	
DDFT	Front, Hind	DDFT	As selected in Method	
Palmar/plantar Annular ligament	Front, Hind	PAL	As selected in Method	
SL lateral extensor branch	Front, Hind	LatExBrSL	As selected in Method	
SL medial extensor branch	Front, Hind	MdExtBrSL	As selected in Method	
Lateral Collateral ligament Superficial portion	Front, Hind	LCLigSup	As selected in Method	
Lateral Collateral ligament Deep portion	Front, Hind	LCLigDeep	As selected in Method	
Medial Collateral ligament Superficial portion	Front, Hind	MCLigSup	As selected in Method	
Medial Collateral ligament Deep portion	Front, Hind	MCLigDeep	As selected in Method	
Proximal Digital Annular ligament	Front, Hind	PrDgAnLig	As selected in Method	
Straight sesamoidean ligament	Front, Hind	StrSesLig	As selected in Method	
Oblique (middle) Distal Sesamoidean ligament-Lateral	Front, Hind	ODSLLat	As selected in Method	
Oblique (middle) Distal Sesamoidean ligament-Medial	Front, Hind	ODSLMed	As selected in Method	
Cruciate Distal Sesamoidean ligaments	Front, Hind	CDSL	As selected in Method	

MyLab – ADVANCED OPERATIONS

Digital sheath –lateral	Front, Hind	DigShLat	As selected in Method
Digital sheath –medial	Front, Hind	DigShMed	As selected in Method
Dorsal Synovial Pad	Front, Hind	DrsSynPad	As selected in Method
Long Digital Extensor	Front, Hind	LngDigExt	As selected in Method
Common Digital Extensor	Front, Hind	ComDigExt	As selected in Method
Lateral Digital Extensor	Front, Hind	LatDigExt	As selected in Method

Pastern

Parameter	Selectable in	Label	Measure
SDFT	Front, Hind	SDFT	As selected in Method
SDF lateral branch	Front, Hind	LatBrSDFT	As selected in Method
SDF medial branch	Front, Hind	MedBrSDFT	As selected in Method
DDFT	Front, Hind	DDFT	As selected in Method
DDF lateral lobe	Front, Hind	LatLDDFT	As selected in Method
DDF medial lobe	Front, Hind	MedLDDFT	As selected in Method
Distal Digital Annular ligament	Front, Hind	DistDAL	As selected in Method
Proximal Digital Annular ligament	Front, Hind	PrDgAnLig	As selected in Method
Superficial Abaxial Palmar ligament	Front, Hind	SpAbPmLig	As selected in Method
Deep Abaxial Palmar ligament medial	Front, Hind	DAbPLMed	As selected in Method
Deep Abaxial Palmar ligament lateral	Front, Hind	DAbPLLat	As selected in Method
Deep Axial Palmar Ligament medial	Front, Hind	DAxPLMed	As selected in Method
Deep Axial Palmar Ligament lateral	Front, Hind	DAxPLLat	As selected in Method

Shoulder

Parameter	Selectable in	Label	Measure
Bicipital Tendon median lobe	Front	BicTMedL	As selected in Method
Bicipital Tendon lateral lobe	Front	BicTLatL	As selected in Method
Bicipital bursa	Front	BicBursa	As selected in Method
Transverse Humeral Ligament	Front	TRSVHumL	As selected in Method

Carpus

Parameter	Selectable in	Label	Measure
Extensor Carpi Radialis (ECR) tendon	Front	ECRT	As selected in Method
Extensor Carpi Radialis (ECR) tendon sheath	Front	ECRTSh	As selected in Method
Common Digital Extensor (CDE) tendon	Front	CDET	As selected in Method

MyLab – ADVANCED OPERATIONS

Common Digital Extensor (CDE) tendon sheath	Front	CDEtSh	As selected in Method
Lateral Digital Extensor (LDE) tendon	Front	LDET	As selected in Method
Lateral Digital Extensor (LDE) tendon sheath	Front	LDEtSh	As selected in Method
Extensor Carpi Oblique (ECO) tendon	Front	ECOT	As selected in Method
Extensor Carpi Oblique (ECO) tendon sheath	Front	ECOTSh	As selected in Method
Ulnaris Lateralis tendon	Front	ULT	As selected in Method
Ulnaris Lateralis tendon sheath	Front	ULTSh	As selected in Method
Carpal Sheath	Front	CarpSh	As selected in Method
SDFT	Front	SDFT	As selected in Method
DDFT	Front	DDFT	As selected in Method
Dorsal Carpal Sheath	Front	DCSh	As selected in Method
Acc Lig DDF	Front	AccDigDDF	As selected in Method
ICL	Front	ICL	As selected in Method
Acc Lig SDF	Front	AccDigSDF	As selected in Method
Superior Check Ligament	Front	SupChLig	As selected in Method
Medial Collateral ligament	Front	MedCollig	As selected in Method
Lateral Collateral ligament	Front	LatCollig	As selected in Method

Tarsus

Parameter	Selectable in	Label	Measure
Gastrocnemius Tendon	Hind	GastTend	As selected in Method
DDFT	Hind	DDFT	As selected in Method
SDFT	Hind	SDFT	As selected in Method
Long Plantar ligament	Hind	LgPLLig	As selected in Method
Peroneus Tertius tendon	Hind	PerTerT	As selected in Method
Long Digital Extensor tendon	Hind	LongDigET	As selected in Method
Tibialis cranialis tendon	Hind	TibCrTend	As selected in Method
Cunean tendon (medial branch of tibialis cranialis)	Hind	CunTendon	As selected in Method
Tarsal Sheath	Hind	TarsalSh	As selected in Method
Calcaneal bursa (between gastrocnemius and SDF)	Hind	CalcBrs	As selected in Method
Gastrocnemius bursa (between gastrocnemius and tuber)	Hind	GastBrs	As selected in Method
Superficial bursa (between skin and SDF)	Hind	SupBrs	As selected in Method
Cunean tendon bursa	Hind	CunTndBrs	As selected in Method
Medial Collateral ligament	Hind	MedCollig	As selected in Method
Lateral Collateral ligament	Hind	LatCollig	As selected in Method
Short Medial Collateral ligament	Hind	ShMdClLig	As selected in Method

MyLab – ADVANCED OPERATIONS

Long Medial Collateral ligament	Hind	LgMdCLig	As selected in Method
Short Lateral Collateral ligament	Hind	ShLtCLig	As selected in Method
Long Lateral Collateral ligament	Hind	LgLtCLig	As selected in Method

Spine

Parameter	Selectable in	Label	Measure
Atlantal Bursa (between nuchal ligament and atlas)	Front	AtBursa	As selected in Method
Supraspinous Bursa (between supraspinous ligament and underlying proximal thoracic vertebrae)	Front	SupBursa	As selected in Method
Supraspinous Ligament	Front	SupLig	As selected in Method
Interspinous Ligament	Front	IntLig	As selected in Method
Articular Facet	Front	ArtFacet	As selected in Method


Stifle

Parameter	Selectable in	Label	Measure
Collateral ligament- lateral	Hind	LatColLig	As selected in Method
Collateral ligament-medial	Hind	MedColLig	As selected in Method
Meniscus-lateral	Hind	MenLat	As selected in Method
Meniscus-medial	Hind	MenMed	As selected in Method
Cranial cruciate ligament	Hind	CrCrucLig	As selected in Method
Caudal cruciate ligament	Hind	CaCrucLig	As selected in Method
Middle patellar ligament	Hind	MidPatLig	As selected in Method
Lateral patellar ligament	Hind	LatPatLig	As selected in Method
Medial patellar ligament	Hind	MedPatLig	As selected in Method
Lateral femoropatellar ligament	Hind	LFemPatLg	As selected in Method
Medial femoropatellar ligament	Hind	MFemPatLg	As selected in Method

The specific measurement can then be performed.

The specific selected limb, zone and region will be indicated on the left side of the screen, as reference, alongside the measurements.


An example is indicated in the picture below:



Setup

The equine tendon measurement package can be configured via the **SETTINGS - MENU – MEASUREMENTS & ANALYSIS – TENDON**.

The following page appears:



In the “Measure Method” list, the user can choose between “Distance Method” or “Zone Method”

DISTANCE METHOD

LABEL
Proximal
Mid
Distal
1 cm
2 cm
4 cm
8 cm
12 cm
16 cm
20 cm
24 cm
28 cm

ZONE METHOD

LABEL
Region 1
Region 1A
Region 1B
Region 2
Region 2A
Region 2B
Region 3
Region 3A
Region 3B
Region 4
Region 4A
Region 4B
Region P1a
Region P1b
Region P1c
Region P1d

“Overall measurement type” will select which type of actual measurements will be performed.

For each selection, the following measurements can be performed:

- for “Distance” – 3 distance measurements
- for “Tr-Area” – 2 area measurements (trace method) and 1 distance measurement

MyLab – ADVANCED OPERATIONS

- for “Ax-Area” – 2 area measurements (ellipse method) and 1 distance measurement
- for “A-Ratio” – 1 area ratio measurement (trace method) and 1 distance measurement
- for “E-Ratio” – 1 area ratio measurement (ellipse method) and 1 distance measurement

In the “Measurement Menu”, the user can select which of the predefined labels will be shown in the measurement menu (activated by pressing the MEASURE button).

By doing this, the number of typical measurements can be tailored to individual preferences.

MSK / Tendon Report

Pressing the **REPORT** tab on the control panel opens the MSK / Tendon Report.

MSK / Tendon Ultrasound Report		
Date	: 11-29-2010	
Owner's name	:	
Animal's name	:	
Species	:	EQUINE
ID	:	2010-11-29-01
Accession number	:	
Diagnosis	:	
Gender	:	Unknown
Doctor's name	:	
Performing Physician	:	
Birthdate	:	
Age	:	0 Years
Weight	:	0.00 Kg
Breed	:	
Neutered	:	Unaltered
FORE LEFT LEG		
CARPUS ACCLIGDDF		
REGION 1		
Distance 1	: 3.92 mm	Distance 2 : 5.40 mm
Distance 3	: 5.80 mm	
REGION 1B		
Distance 1	: 4.12 mm	Distance 2 : 4.13 mm
Distance 3	: 4.74 mm	
REGION 2B		
Distance 1	: 3.94 mm	Distance 2 : 5.59 mm
Distance 3	: 5.47 mm	
HIND LEFT LEG		
FETLOCK COMMONDIGEXTENSOR		
REGION 1		
Distance 1	: 3.19 mm	Distance 2 : 4.01 mm
Distance 3	: 3.66 mm	

The **COMMENTS** key allows the user to add comments at the end of the report. When pressing the comments key, the comments window appears allowing data to be edited in the empty field by use of the keyboard. Also predefined comment can be added to the report by checking the involved checkbox. These comments can be defined in the **REPORT** menu of the **DISPLAY** menu in the **SYSTEM SETUP** menu


The **OBSERVATIONS** key gives the opportunity for quick and easy reporting of general characteristics observed during the examination.

The **MEASUREMENTS** key opens a window which displays all the measurements and calculations performed. Each measurement can be selected to be shown in the report. Different sections, representing different legs, can be selected by the


MyLab – ADVANCED OPERATIONS

corresponding tabs. The additional page displays all generic measurements which have been added to the report with the **ADD TO REPORT** key. In addition there is the option to rename the label of the measurement. To rename the label, select a measurement from the list and enter a new name in the “Edit Label” field and confirm with **enter**.

The **IMAGES** key makes it possible to add pictures to the report that were taken during the exam. In the submenu **IMAGES** the pictures can be selected by ticking their corresponding boxes.



To add the pictures to the report press **IMAGES** once more. The result can be seen in the picture below.



5 - Reproductive Calculation Packages

In this chapter the application package for Reproductive Calculation is described.

Note

This application calculation package is only available for farm animals (porcine, ovine and llama) when the corresponding license is activated. Please refer to the System Configuration manual for all the options/licenses of the **MyLab** and their status.

Reproductive icon




REPRO

Reproductive package includes calculation packages for large animals: equine and bovine, small animals: canine and feline and farm animals: porcine, ovine and llamas.

When Reproductive application is selected by the **REPRO** application key, pressing the **MEASURE** key opens the list of measurements on the right of the screen.

Once activated the application, the system displays the following menu:



Reproductive Calculations in B-Mode

In B-mode, the user can calculate the gestational age in the following applications: Equine, Bovine, Canine and Feline. Optionally, also Farm Animals: Swine, Ovine and Llamas are provided.

Refer to “System Configuration” section for customization of veterinary reproductive calculations

Gestational Age in Reproductive Canine Application

In Reproductive Canine application, the following parameters can be measured:

Repro Canine Application Icon



Parameter		Label	Measure
GS	Gestational Sac diameter	GS	Distance
CRL	Crown Rump Length	CRL	Distance
HD	Head Diameter	HD	Distance
BD	Body Diameter	BD	Distance

Every time a parameter is measured, the gestational age is automatically calculated.

All measurements will appear in the report.

Up to 4 fetuses can be measured, either in single or compared mode.

Based on the calculated gestational age, the machine will calculate also the DBP (Days Before Parturition) value, which will appear on the first page of the report.

Repro Feline
Application Icon



Gestational Age in Reproductive Feline Application

In Reproductive Feline application, the following parameters can be measured:

Parameter		Label	Measure
HD	Head Diameter	HD	Distance
BD	Body Diameter	BD	Distance

Every time a parameter is measured (with the exception of the amniotic index), the gestational age is automatically calculated.

All measurements will appear in the report.

Up to 4 fetuses can be measured, either in single or compared mode.

Based on the calculated gestational age, the machine will calculate also the DBP (Days Before Parturition) value, which will appear on the first page of the report.

Formulas and Bibliographic References

Equations to calculate the fetal age in canine, feline, and the number of days before parturition are indicated below:

CANINE

GSD FORMULA CANINE- Nyland, Mattoon

(Based on days post-luteinizing hormone (LH) surge in the female dog; valid for gestational age less than 40 days)

$$\text{Gestational Age} = (6 \times \text{GSD}) + 20$$

GSD (in cm) = gestational sac diameter

Precision : ± 3 days

CRL FORMULA CANINE- Nyland, Mattoon

(Based on days post-luteinizing hormone (LH) surge in the female dog ; valid for gestational age less than 40 days)

$$\text{Gestational Age} = (3 \times \text{CRL}) + 27$$

CRL (in cm) = crown – rump length

Precision : ± 3 days

HD FORMULA CANINE- Nyland, Mattoon

(Based on days post-luteinizing hormone (LH) surge in the female dog ; valid for gestational age more than 40 days)

$$\text{Gestational Age} = (15 \times \text{HD}) + 20$$

HD (in cm) = head diameter

Precision : ± 3 days

BD FORMULA CANINE- Nyland, Mattoon

(Based on days post-luteinizing hormone (LH) surge in the female dog ; valid for gestational age more than 40 days)

$$\text{Gestational Age} = (7 \times \text{BD}) + 29$$

MyLab – ADVANCED OPERATIONS

BD (in cm) = body diameter

Precision : ± 3 days

DBP FORMULA CANINE- Nyland, Mattoon

(Based on 65 +/- 1 days post breeding in the dog)

DBP = 65 – GA

DBP (in days) = days before parturition

GA (in days) = gestational age

FELINE

HD FORMULA FELINE- Nyland, Mattoon

(Based on days post breeding in the cat; valid for gestational age more than 40 days)

Gestational Age = (25 x HD) + 3

HD (in cm) = head diameter

Precision: ± 2 days

BD FORMULA - Nyland, Mattoon

(Based on days post breeding in the cat; valid for gestational age more than 40 days)

Gestational Age = (11 X BD) + 21

BD (in cm) = body diameter

Precision: ± 3 days

DBP FORMULA FELINE- Nyland, Mattoon

(Based on 61 +/- 1 days post breeding in the cat)

DBP = 61 – GA

DBP (in days) = days before parturition

GA (in days) = gestational age

Bibliography for reproductive canine and feline applications measurements and calculation

The gestational age and number of days before parturition have been estimated based on the following bibliographic references:

Parameter	Bibliography
GSD CANINE	Nyland, Thomas G; Mattoon, John S – <i>Veterinary Diagnostic Ultrasound</i> – 1995 – ISBN 0-7216-2745-5
CRL CANINE	England GCW, Allen WE, Porter DJ – <i>Studies on canine pregnancy using B-mode ultrasound: development of the conceptus and determination of the gestational age</i> – Journal of Small Animal Practice 1990; 31:324-329
HD CANINE	Yeager AE, Mohammed HO, Meyers-Wallen V, et. al – <i>Ultrasonographic appearance of the uterus, placenta, fetus, and fetal membranes throughout accurately timed pregnancy in beagles</i> – American Journal of Veterinary Research – 1992; 53:342 – 351
BD FELINE	Beck KA, Baldwin CJ, Bosu WTK – <i>Ultrasound prediction of parturition in the queen</i> – Veterinary Radiology – 1990:31:32-35)

Gestational Age in Reproductive Equine Application

In Reproductive Equine application, the following parameters can be measured:



Parameter		Label	Measure
GSD	Gestational Sac Diameter	GS	Distance
ERD	Equine Rump Diameter	ERD	Distance
ESD	Equine Skull Diameter	ESD	Distance
EED	Equine Eye ball Diameter	EED	Distance

Every time a parameter is measured, the gestational age is automatically calculated.

All measurements will appear in the report.

Formulas and Bibliographic References

Equations to calculate the fetal age in equine and the number of days before parturition are indicated below:

EQUINE

GSD FORMULA EQUINE - Pipers

$$\text{Gestational Age} = 6.666666667 * \text{GSD} + 2.733333333$$

GSD (in cm) = Gestational Sack Diameter

Range : 0 weeks 3days // 6 weeks 0days

RUMP DIAMETER FORMULA EQUINE - Kähn

$$\text{Gestational Age} = 190.0 - \frac{\sqrt{144400 - [4 * (250 * \text{ERD} + 14082.5)]}}{2}$$

ERD (in mm) = Equine Rump Diameter

Range = 6 weeks 0days // 27 weeks 0days

SKULL DIAMETER FORMULA EQUINE - Kähn

$$\text{Gestational Age} = \sqrt{\frac{\text{ESD}}{0.0013}} - 7753.8462$$

ESD (in mm) = Equine Skull Diameter

Range = 11 weeks 2days // 30 weeks 0days

EYE BALL DIAMETER FORMULA EQUINE - Kähn

$$\text{Gestational Age} = (\text{EED} / 0.14) - 5.5$$

EED (in mm) = Equine Eye ball Diameter

Range : 11 weeks 2days // 47 weeks 0days

Bibliography for reproductive equine application measurements and calculation

The gestational age has been estimated based on the following bibliographic references:

Parameter	Bibliography
GSD	Pipers FS, Zent W, Holder R, et al. - Ultrasonography as an adjunct to pregnancy assessment in the mare, Journal of the American Veterinary Medical Association 1984;184: 328-334
EQUINE	
ERD	Kähn W, Leidl W - Deutsche Tierärztliche Wochenschrift: 505-515, 9/1987
EQUINE	
ESD	Kähn W, Leidl W - Deutsche Tierärztliche Wochenschrift: 505-515, 9/1987
EQUINE	
EED	Kähn W, Leidl W - Deutsche Tierärztliche Wochenschrift: 505-515, 9/1987
EQUINE	

Gestational Age in Reproductive Bovine Application

In Reproductive Equine application, the following parameters can be measured:

Repro Bovine
Application Icon



Parameter		Label	Measure
CRL	Crown Rump Length	CRL	Distance
SD	Bovine Skull Diameter	SD	Distance
TD	Bovine Trunk Diameter	TD	Distance
UD	Bovine Uterine Diameter	UD	Distance

Every time a parameter is measured, the gestational age is automatically calculated.

All measurements will appear in the report.

Formulas and Bibliographic References

Equations to calculate the fetal age in bovine and the number of days before parturition are indicated below:

BOVINE

CRL FORMULA BOVINE - Kähn

$$\text{Gestational Age} = \sqrt{(50.251256 * CRL) + 540.70352}$$

CRL (in mm) = Crown-Rump Length

Range = 4 weeks 2 days // 11 weeks 4days

SKULL DIAMETER FORMULA BOVINE – Kähn

$$\text{Gestational Age} = 26.5965583174 + \frac{SD}{0.523}$$

SD (in mm) = Skull Diameter

Range = 8 weeks 4 days// 30 weeks 0 days

TRUNK DIAMETER FORMULA BOVINE -White

$$\text{Gestational Age} = 41.152263 * \ln(TD) + 35.946502$$

TD (in cm) = Trunk Diameter

Range = 9 weeks 1day // 19 weeks 4 days

UTERINE DIAMETER FORMULA BOVINE - White

$$\text{Gestational Age} = 46.511628 * \ln(UD) + 0.5674419$$

UD (in cm) = Uterine Diameter

Range = 4 w 2 // 18 w 4

Bibliography for reproductive bovine application measurements and calculation

The gestational age has been estimated based on the following bibliographic references:

Parameter	Bibliography
CRL BOVINE	Kähn, W. - Sonographic Fetometry in the Bovine.: 1105-1121, 1989
SD BOVINE	Kähn,W; Leidl, W. - Deutsche Tierärztliche Wochenschrift.: 505-515,9/1987
TD BOVINE	White IR, Russell AJF, Wright, Whyte IK. - The Veterinary Record.:5-8, July 6, 1985
UD BOVINE	

Gestational Age in Reproductive Porcine Application

In Reproductive Equine application, the following parameters can be measured:

Repro Swine
Application Icon



Parameter		Label	Measure
HLA	Heart Long Axis	HLA	Distance
SLA	Stomach Long Axis	SLA	Distance

Every time a parameter is measured, the gestational age is automatically calculated.

All measurements will appear in the report.

Formulas and Bibliographic References

Equations to calculate the fetal age in porcine and the number of days before parturition are indicated below:

HEART LONG AXIS FORMULA PORCINE - Fraunholz

$$\text{Gestational Age} = 10.0344827586 + \frac{\text{HLA}}{0.29}$$

HLA (in mm) = Heart Long Axis distance

Range = 3 weeks 0 days // 16 weeks 0 days

STOMACH LONG AXIS FORMULA PORCINE – Fraunholz

$$\text{Gestational Age} = 16.0263157894 + \frac{\text{SLA}}{0.38}$$

SLA (in mm) = Stomach Long Axis distance

Range = 3 weeks 0 days // 16 weeks 0 days

Bibliography for reproductive porcine application measurements and calculation

The gestational age has been estimated based on the following bibliographic references:

Parameter	Bibliography
HLA	
SWINE	Fraunholz, J.; Kähn, W.; Leidl, W. - Mh. Vet.-Med. 44(1989)
SLA	: 425-430
SWINE	

Repro Ovine
Application Icon



Gestational Age in Reproductive Ovine Application

In Reproductive Equine application, the following parameters can be measured:

Parameter	Label	Measure
CRL	Crown Rump Length	CRL Distance

Every time a parameter is measured, the gestational age is automatically calculated.

All measurements will appear in the report.

Formulas and Bibliographic References

Equations to calculate the fetal age in ovine and the number of days before parturition are indicated below:

CRL TABEL OVINE - Van den Boogard

Diameter (mm)	Weeks
13.3	4
28.3	5
45.0	6
65.0	7
85.0	8
105.0	9
123.3	10

Bibliography for reproductive ovine application measurements and calculation

The gestational age has been estimated based on the following bibliographic references:

Parameter	Bibliography
CRL Ovine	van den Boogaard AEJM, Smeets FAM, Slangen FMA, Muijtens A.- Real time ultrasound pregnancy testing in ewes and estimation of fetal age - European Journal of Obstetrics & Gynecology and Reproductive Biology, Vol.22, 4,1986, 266-267



Gestational Age in Reproductive Llamas Application

In Reproductive Equine application, the following parameters can be measured:

Parameter	Label	Measure
BPD Biparietal Diameter	BPD	Distance

Every time a parameter is measured, the gestational age is automatically calculated.

All measurements will appear in the report.

Formulas and Bibliographic References

Equations to calculate the fetal age in llama and the number of days before parturition are indicated below:

BIPARIETAL DIAMETER FORMULA LLAMAS - Haibel

$$\text{Gestational Age} = 18.8 + (3.79 * \text{BPD})$$

BPD (in mm) = Biparietal Diameter

Range = 9 weeks 3 days // 33 weeks 4 days

Bibliography for reproductive ovine application measurements and calculation

The gestational age has been estimated based on the following bibliographic references:

Parameter	Bibliography
BPD Llamas	Haibel, G.K. ; Fung, E.D. - Real-time biparietal diameter measurement for the prediction of gestational age in llamas, Theriogenology. 1988 Dec ;30 (6):1053-7 17087893

Setup

The obstetric measurement list can be configured via the [SETTINGS – MENU - MEASUREMENTS AND ANALYSIS – REPRO](#).

The Repro Setup section includes two sub-pages to configure the biometric measurements and the calculations.

Gestational Age



Up to 6 measurements can be configured in the system.

Every measurement can be configured with the following parameters:

Item	Description
Enable	Enables the measurement. If the measurement is enabled it will be included in the list when the Measure key is pressed.
Average	Enables or disables the average calculation on the measurement. If Average is disabled, only one measurement will be saved in the report. In case of repeated measurements, only the last one will be saved. If Average is enabled, up to three measurements can be saved in the report and used for the average calculation.
Measure	This pull-down menu allows the application measurements available in the system to be selected.
Table	This pull-down menu allows the Fetal Age table to be selected.
Method	This pull-down menu allows the way the measurement is performed to be selected (Area by Ellipse, Area by Trace, etc).


Editing User Tables

Edit User Tables

When the Gestational Age page of Setup is displayed, the **EDIT USER TABLES** key is available on the screen.

Pressing **EDIT USER TABLES** displays the following customization environment:

MyLab – ADVANCED OPERATIONS



This environment allows the user to define new tables for gestational age for new measurements or for existing parameters.

A maximum of 10 custom tables can be generated (selected by Edit User Table pull-down menu).


Item	Description
Select Parameter	When Select Parameter is checked, the new custom table will be associated to an existing parameter (like BPD, TAD, etc) selected via the pull-down menu beside. In this case, values measured for these parameters will be used for calculations including the parameters.
New Parameter	When New Parameter is selected, the new custom table will be associated to a new parameter made by the user and not included in the OB parameter's list.
Measure	When New Parameter is selected, the user is allowed to edit in this field a description for the new parameter.
Table	In this field the user can edit the name of the new custom table.
Method	This pull-down menu is active only if New Parameter has been selected and allows the user to select how the new parameter has to be measured (by distance, length, area, etc). If Select Parameter has been checked, this menu is not available for the user and displays the measurement method for the selected parameter.

In the lower-left part of the screen, the user can specify the values in mm and relative GA in days with positive and negative deviations.

According to the input values, the graph is updated in the lower-right part of the window.

NOTE: the maximum number of input values that can be entered in a new table is 80.

General settings



In this page it is possible to select:

1. Enable Instant Ref display to have the gestational age displayed when a measurement is performed.
2. Enable the graphs to put in the report:
No Graphs – no graphs will be entered in the report;
All Graphs – all available graphs will be entered in the report;

Reproductive Report

Pressing the **REPORT** tab on the control panel opens the Repro Report.

Repro Ultrasound Report

Date : 11-30-2010

Owner's name : Doctor's name :
 Animal's name : Performing Physician :
 Species : EQUINE Birthdate :
 ID : 2010-11-30-02 Age : 0 Years
 Accession number : Weight : 0.00 Kg
 Diagnosis : Breed :
 Gender : Female Neutered : Unaltered

FETUS A

AVERAGE ULTRA SOUND AGE **DAYS BEFORE PARTURATION**

GA : 4w1d EDD : 09-27-2011 DBP : 301 +/- 10

MEASUREMENTS

GS : 1.59 cm (1.25 cm, 1.80 cm, 1.73 cm)	GA : 1w6d	FA Table : Pipers
ERD : 2.47 mm (3.26 mm, 2.42 mm, 1.73 mm)	GA : 6w2d	FA Table : Kaehn
ESD : 1.45 cm (1.62 cm, 1.64 cm, 1.08 cm)	GA : --w-d	FA Table : Kaehn
EED : 1.07 cm (0.89 mm, 1.06 cm, 1.10 cm)	GA : --w-d	FA Table : Kaehn

GRAPHS

GS Pipers Fetus A

50.0
42.3
33.9

AGA (w)

Control Panel: FETUS 0, MEASUREMENTS, GRAPHS, COMMENTS, IMAGES, PRINT, PDF

With the **FETUS** key, the report for different fetuses can be selected and configured.


The **MEASUREMENTS** key opens a window which displays all the measurements and calculations performed. Each measurement can be selected to be shown in the report. The additional page displays all generic measurements which have been added to the report with the **ADD TO REPORT** key. In addition there is the option to rename the label of the measurement. To rename the label, select a measurement from the list and enter a new name in the “Edit Label” field and confirm with **enter**.

The **GRAPHS** key gives the opportunity to configure the graphs displayed in the report.


The **COMMENTS** key allows the user to add comments at the end of the report. When pressing the comments key, the comments window appears allowing data to be edited in the empty field by use of the keyboard. Also predefined comment can be added to the report by checking the involved checkbox. These comments can be defined in the **REPORT** menu of the **DISPLAY** menu in the **SYSTEM SETUP** menu

The **IMAGES** key makes it possible to add pictures to the report that were taken during the exam. In the submenu **IMAGES** the pictures can be selected by ticking their corresponding boxes.

MyLab – ADVANCED OPERATIONS



To add the pictures to the report press **IMAGES** once more. The result can be seen in the picture below.




Graphs Display

In the **GRAPHS** window both Curve and Bar graphs can be selected to be displayed.

Bar Structure

All performed measurements are displayed in one graph. Gestational weeks and days are on the horizontal axes (one axis per measured item). In the same screen, the referenced author/table is also shown.

The GA (US) is the gestational age estimated by the performed measurements and is displayed by a vertical line. The gestational age from each single measurement is displayed by a cross. A range (present according to the selected table) is displayed by a colored rectangle around the cross.



Curve Structure

All performed measurements are displayed in one graph. Gestational weeks and days are on the horizontal axis. The measured value is shown on the vertical axis. The referenced author/table is also shown.

The AGA (US) is the gestational age estimated by the performed measurements and is displayed by a vertical line. A range (present according to the selected table) is displayed by a lines going along the normal value according to the specific reference.

6 - Animal Science Package

In this chapter the application package for Animal Science Calculation is described.

Note

This application calculation package is only available for large and farm animals (bovine, swine, ovine).

Animal Science icon



The Animal Science application package contains a measurement package for measuring the hide-, fat- and loin thickness within cattle, sheep and swine and a percentage intramuscular fat calculation for cattle (Yield Grade).

The Animal Science package can be extended with the following optional license:

- Animal Science IMF license (article number: **302012100**), to activate the QUIP index and/or QUIP index Lean in Bovine species. The Animal Science IMF license also contains the Animal Science license mentioned above.

When the Animal Science application is selected by selecting the **ANIMAL SCIENCE** application key, the **MEASURE** key (in **FREEZE** mode) opens the list of measurements on the right of the screen.

Once activated the application, depending on the activated license and the selected species, different calculation packages will be available:


MyLab – ADVANCED OPERATIONS

- **Bovine**


Refer to “System Configuration” section for customization of veterinary animals science calculations



- **Porcine**



- **Ovine**



Animal Science Calculations in B-Mode

In B-mode, depending upon the settings of the measurements (see below SETUP), the user can calculate Yield Grade, QUIP index and/or QUIP index Lean in Bovine species and Yield Grade in Porcine and Ovine species.

Yield Grade in Bovine Animal Science Application

In Bovine Animal science application, the following parameters can be measured:

*Animal Science Application
Icon*



Parameter		Label	Measure
HT	Hide thickness	HT	Distance
FT	Fat thickness	FT	Distance
LT	Loin thickness	LT	Distance
LA	Loin Area	LA	Area
YG	Yield grade	YG	

Selecting the item **Hide thickness** activates the measurement of the hide layer. Touch the screen on the position where the measurement has to start, a cross will appear in the image. Touch the screen on the end position of the measurement, the result will be displayed on the left side of the screen.

MyLab – ADVANCED OPERATIONS

Selecting the item **Fat thickness** activates the measurement of the fat layer. Touch the screen on the position where the measurement has to start, a cross will appear in the image. Touch the screen on the end position of the measurement, the result will be displayed on the left side of the screen.

Selecting the item **Loin thickness** activates the measurement of the Loin thickness. Touch the screen on the position where the measurement has to start, a cross will appear in the image. Touch the screen on the end position of the measurement, the result will be displayed on the left side of the screen.

Selecting the item **Loin area** activates the Area and circumference measurement of the loin-eye. Touch the screen on the position where the measurement has to start and don't release the screen, draw a contour around the loin-eye and release the screen. The result will be displayed on the left side of the screen.

Note

This measurement can only be done in a transverse image of the Loin Eye.

To get the calculation of the **Yield Grade**, insert the weight of the animal in Patient & Archive page and measure the fat thickness and Loin area.

All performed measurements will appear in the report.

Yield Grade in Porcine Animal Science Application

In Porcine Animal science application, the following parameters can be measured:

Parameter		Label	Measure
FT	Fat thickness	FT	Distance
LT	Loin thickness	LT	Distance
LA	Loin Area	LA	Area

Selecting the item **Fat thickness** activates the measurement of the fat layer. Touch the screen on the position where the measurement has to start, a cross will appear in the image. Touch the screen on the end position of the measurement, the result will be displayed on the left side of the screen.

Selecting the item **Loin thickness** activates the measurement of the Loin thickness. Touch the screen on the position where the measurement has to start, a cross will appear in the image. Touch the screen on the end position of the measurement, the result will be displayed on the left side of the screen.

Selecting the item **Loin area** activates the Area and circumference measurement of the loin-eye. Touch the screen on the position where the measurement has to start and don't release the screen, draw a contour around the loin-eye and release the screen. The result will be displayed on the left side of the screen.

Animal Science Application
Icon



Note

This measurement can only be done in a transverse image of the Loin Eye.

Based on the measured fat thickness, **Loin area** and weight (previously insert in the Patient & Archive page), the system will calculate Percent Lean (cross) which will appear on the first page of the report.

Based on the measured fat thickness, **Loin thickness** and weight (previously insert in the Patient & Archive page) , the system will calculate Percent Lean (lateral) which will appear on the first page of the report.

All performed measurements will appear in the report.

Yield Grade in Ovine Animal Science Application

In Bovine Animal science application, the following parameters can be measured:

*Animal Science Application
Icon*



Parameter		Label	Measure
FT	Fat thickness	FT	Distance
LT	Loin thickness	LT	Distance
LA	Loin Area	LA	Area
YG	Yield grade	YG	

Selecting the item **Fat thickness** activates the measurement of the fat layer. Touch the screen on the position were the measurement has to start, a cross will appear in the image. Touch the screen on the end position of the measurement, the result will be displayed on the left side of the screen.

Selecting the item **Loin thickness** activates the measurement of the Loin thickness. Touch the screen on the position were the measurement has to start, a cross will appear in the image. Touch the screen on the end position of the measurement, the result will be displayed on the left side of the screen.

Selecting the item **Loin area** activates the Area and circumference measurement of the loin-eye. Touch the screen on the position were the measurement has to start and don't release the screen, draw a contour around the loin-eye and release the screen. The result will be displayed on the left side of the screen.

Note

This measurement can only be done in a transverse image of the Loin Eye.

To get the calculation of the **Yield Grade** insert the weight of the animal in Patient & Archive page and measure the fat thickness and Loin area.

All performed measurements will appear in the report.

Bibliographic References

- International Study Guide, Real-time Ultrasound Beef cattle applications; Jerry D. Gresham PhD.
- International Study Guide, Real-time Ultrasound Swine applications, Jerry D. Gresham PhD.

QUIP index in Bovine Animal Science Application

For the calculation of the Intramuscular Fat, two methods are available: one is according to Dr. Gresham’ equation(Quip Index), one is according to Dr.Aass’s equation (Quip Index Lean). (refer to bibliographic references).

In Animal science Bovine application, the following parameters can be measured:

- Quip Index

Parameter		Label	Measure
QUIP1	Quip Index 1	QI1	
QUIP2	Quip Index 2	QI2	
BOX H	Box Height	BH	cm
IMF	Intramuscular Fat	IMF	
AVIMF	Average Intarmuscular Fat	AVIMF	

- Quip Inedx Lean

Parameter		Label	Measure
QUIP1	Quip Index 1	QI1	
QUIP2	Quip Index 2	QI2	
BOX H	Box Height	BH	cm
IMF L	Intramuscular Fat Lean	IMF	
AVIMF	Average Intarmuscular Fat	AVIMF	

The Quip index is a single measurement based on the average gray level inside the box which is used to calculate the Intramuscular Fat.

Once selected **QUIP INDEX/QUIP INDEX LEAN**, touch the screen where the ROI (rectangle) has to be placed. By touching and dragging the corners of the box it is possible to change the dimension of the ROI (rectangle).

All performed measurements will appear in the report.

*Animal science
Bovine Application
Icon*




Bibliography for animals science bovine and porcine applications measurements and calculation

- International Study Guide, Real-time Ultrasound Beef cattle applications; J.D. Gresham.
- Ultrasound prediction of intramuscular fat content in lean cattle. L.Aass, C.G. Fristedt, J.D. Gresham

Setup

The animal science measurement list can be configured via the **SETTINGS – MENU - MEASUREMENTS AND ANALYSIS – ANIMAL SCIENCE**.

Bovine Yield Grade, Bovine quip Index, Bovine Quip Index Lean can be enabled or disabled.



Every measurement can be configured with the following parameters:

Item	Description
Enable	Enables the measurement. If the measurement is enabled it will be included in the list when the Measure key is pressed.
Average	Enables or disables the average calculation on the measurement. If Average is disabled, only one measurement will be saved in the report. In case of repeated measurements, only the last one will be saved. If Average is enabled, up to three measurements can be saved in the report and used for the average calculation.

Animal Science Report

Pressing the **REPORT** tab on the control panel opens the Animal Science Report.

In the report the measurements of up to three images will be displayed (between brackets) and the average value of them will be displayed as first number.

Animal Science Ultrasound Report

Date : 10-Sep-2013

Owner's name :	Doctor's name :
Animal's name :	Performing Physician :
Species : BOVINE	Birthdate :
ID :	Age :
Accession number :	Weight : 500.00 Kg
Diagnosis :	Breed :
Gender : Unknown	Neutered : Unknown

B-MODE

YIELD GRADE

Hide thickness :	3.18 mm	Loin area :	1.25 cm ²
Fat thickness :	2.78 mm	Yield grade :	6.15
Loin thickness :	2.75 mm		

QUIP INDEX

Quip index 1 :	0.08	(0.08; 0.08)
Quip index 2 :	0.07	(0.07; 0.07)
Box height :	2.72 cm	(3.55 cm; 1.89 cm)
Intramuscular Fat :	1.51	(1.46; 1.55)

QUIP INDEX LEAN


Quip index 1 :	0.19	(0.19; 0.19)
Quip index 2 :	0.19	(0.19; 0.19)
Box height :	1.62 cm	(1.31 cm; 1.99 cm)
Intramuscular Fat Lean :	1.79	(1.62; 1.96)

MEASUREMENTS
COMMENTS
IMAGES
PDF

Refer to "Archive" section for Export options


In parallel with the report, it is possible to create a file with all the measured and calculated values in a structured format. Please refer to the ARCHIVE section.

The **MEASUREMENTS** key opens a window which displays all the measurements and calculations performed. Each measurement can be selected to be shown in the report.




The additional page displays all generic measurements which have been added to the report with the **ADD TO REPORT** key. In addition there is the option to rename the label of the measurement. To rename the label, select a measurement from the list and enter a new name in the “Edit Label” field and confirm with **enter**.

The **COMMENTS** key allows the user to add comments at the end of the report. When pressing the comments key, the comments window appears allowing data to be edited in the empty field by use of the keyboard. Also predefined comment can be added to the report by checking the involved checkbox. These comments can be defined in the **REPORT** menu of the **DISPLAY** menu in the **SYSTEM SETUP** menu




MyLab – ADVANCED OPERATIONS

The **IMAGES** key makes it possible to add pictures to the report that were taken during the exam. In the submenu **IMAGES** the pictures can be selected by selecting their corresponding boxes.



To add the pictures to the report press **IMAGES** once more. The result can be seen in the picture below.



Animal Science measurements, explanation

Bovine measurements

- Hide thickness
- Fat thickness
- Loin thickness
- Loin area
- Quip index
- Quip index Lean

Porcine measurements

- Fat thickness
- Loin thickness
- Loin area

Ovine measurements

- Fat thickness
- Loin thickness
- Loin area

Introduction

In order to get the best results it is important to have a good contact between probe and skin or between probe, standoff and skin. To obtain good acoustical contact, it is sometimes advised to clip the hair before scanning. If the scanning site has not been clipped, the image quality may be improved by insuring the removal of loose hair and dirt by use of a brush or metal comb. The best coupling material for scanning is vegetable oil and not scanning gel.

Note


Never use mineral oil, This oil can damage the probe surface and cables.

For cross-sectional scanning a stand-off has to be attached to the ASP probe to meet the curvature of the animal.

Note

Esaote Europe can provide a stand-off for swine scanning and a stand-off for beef scanning.

Apply some scanning gel on the stand-off before attaching the stand-off to the ASP probe. This gel ensures acoustical contact between the probe and standoff.



A sketch of the position of the probe on the animal is displayed above

Cleaning of the stand-off for swine and beef

The stand-off should be cleaned with warm water and soap after each use and stored in its container to prevent the damage and drying out of the flexible material.

BEEF MEASUREMENTS

The imaging of the backfat and the longissimus muscle (loin eye) is performed approximately 5 cm lateral from the spinous processes of the spine and centered over the 12th rib. Ultrasound gel or vegetable oil can be used as a couplant. With some bulls shaving of the skin is necessary to get a better contact. The probe is placed directly on the surface of the animal or a standoff path is placed between the probe and the surface of the animal.

The probe can be placed in a longitudinal or transverse plane. In the longitudinal plane only distance measurements are made and the intramuscular fat measurement in cattle. In a transverse plane a distance measurement of the Hide-, and fatlayer, the loin thickness and a circumference measurement of the Longissimus muscle are performed.

In the transverse and longitudinal image one can distinguish three distinct curved lines that will be parallel. The outer line will be the outer surface of the skin or animal hide. The middle line will be the interface between of hide and outer surface of the fat layer. The inner line will be the reflection of the interface between the inner surface of the fatlayer and the outer surface of the longissimus dorsi (loin eye) muscle.

In the transverse image the loin eye muscle should be very visible to include both the lateral and medial surfaces and the intercostals visible at the bottom

MyLab – ADVANCED OPERATIONS

of the loin eye muscle.

Total fat depth is determined by measuring the distance from the outer layer of the skin to the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle.

The muscle diameter is determined by measuring the distance between the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle to the ventral surface of the longissimus muscle. The Yield Grade is the ratio between the weight and the measurements of the fat layer and the muscle.

Hide thickness

The hide thickness can be measured at the 12th and 13th rib, in a transverse or longitudinal plane.

Select the item Hide thickness and touch the screen at the outer edge of the outer line. Touch the screen again at the outer edge of the middle line. The result will be displayed on the left side of the image.

Fat thickness

The fat thickness can be measured at the 12th and 13th rib, either in a transverse or longitudinal plane. In the transverse plane the fat thickness is measured by selecting a site to approximate of the length of the loin eye muscle.

Select the item Fat thickness and touch the screen at the outer edge of the middle line. Touch the screen again at the outer edge of the inner line. The result will be displayed on the left side of the image.

Loin thickness

The Loin thickness can be measured at the 12th and 13th rib, in a transverse or longitudinal plane.

Transverse plane

Select the item loin thickness and touch the screen at the outer edge of loin eye muscle. Touch the screen again at the outer edge of the bottom border of the loine eye muscle. The result will be displayed on the left side of the image.

Longitudinal plane

Select the item loin thickness and touch the screen at the outer edge of loin eye muscle. Touch the screen again at a midpoint between the ribs. The result will be displayed on the left side of the image.

Loin area

The Loin area can be measured at the 12th and 13th rib, in a transverse plane. The loin area is a circumference measurement.

Select the item loin area and Touch the screen at the outer edge of loin eye Muscle and don't release the screen. Draw a contour along the outer edge of the loin eye, back to the original starting point. The result will be displayed on the left side of the image.

Quip index/ Quip index Lean

The Quip index measurement has to be done in a longitudinal scan of the loin eye at the 12th and 13th rib. The probe is placed in a plane parallel to the mid-line of

the animal at a point to approximate the midpoint of the loin eye muscle lateral of the spinal column.

Select the item Quip index/Quip index Lean, this activates the intra muscular fat measurement of the loin-eye. A rectangle will appear in the image. Touch the screen and drag the bottom line of the rectangle down and position this line just above the ribs. By touching and dragging the corners of the box it is possible to change the dimension of the rectangle. The result will be displayed on the left side of the screen.

SWINE MEASUREMENTS

The imaging of the backfat and the longissimus muscle (loin eye) is performed approximately 5 cm lateral from the spinous processes of the spine and centered over the 10th rib. Ultrasound gel or vegetable oil is used as a couplant.

The probe is placed directly on the surface of the animal or a standoff path is placed between the probe and the surface of the animal. The probe can be placed in a longitudinal or transverse plane. In the longitudinal plane only distance measurements are made and in a transverse plane a distance measurement of the fat layer and a circumference measurement of the Longissimus muscle.

In the image one can distinguish from top to bottom the fat, which consists of three layers, the muscle and the bowel. The first and second layer of fat are always visible on the ultrasound image. the third layer is not always present, depending on the leanness of the hog. The third layer will then present itself as a thick white band in the image.

Total Fat depth is determined by measuring the distance from the outer layer of the skin to the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle.

The muscle diameter is determined by measuring the distance between the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle to the ventral surface of the longissimus muscle. The % lean is the ratio between the measurements of the fat layer and the muscle.

Fat thickness

Select the item Fat thickness and touch the screen at the outer layer of the skin. Touch the screen again at the outer edge of the longissimus muscle. The result will be displayed on the left side of the image.

Loin thickness

Select the item loin thickness and touch the screen at the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle. Touch the screen again at the ventral surface of the longissimus muscle. The result will be displayed on the left side of the image.

Loin area

Select the item loin area and touch the screen at the outer edge of loin eye muscle and don't release the screen. Draw a contour along the outer edge of the loin eye back to the original starting point. The result will be displayed on the left side of the image.

SHEEP MEASUREMENTS

The imaging of the backfat and the longissimus muscle (loin eye) is performed approximately 5 to 10 cm lateral from the spinous processes of the spine and centered over the 12th rib. Ultrasound gel or vegetable oil can be used as a couplant. With some bulls shaving of the skin is necessary to get a better contact. The probe is placed directly on the surface of the animal or a standoff path is placed between the probe and the surface of the animal.

The probe can be placed in a longitudinal or transverse plane. In the longitudinal plane only distance measurements are made and in a transverse plane a distance measurement of the fatlayer and a circumference measurement of the Longissimus muscle are performed. In the image one can distinguish from top to bottom the fatlayer, the muscle and the bowel.

Total fat depth is determined by measuring the distance from the outer layer of the skin to the the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle.

The muscle diameter is determined by measuring the distance between the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle to the ventral surface of the longissimus muscle. The Yield Grade is the ratio between the weight and the measurements of the fatlayer and the muscle.

Fat thickness

Select the item Fat thickness and touch the screen at the outer layer of the skin. Touch the screen again at the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle. The result will be displayed on the left side of the image.

Loin thickness

Select the item loin thickness and touch the screen at the interface of the bottom layer of fat and the dorsal surface of the longissimus muscle. Touch the screen again at the ventral surface of the longissimus muscle. The result will be displayed on the left side of the image.

Loin area

Select the item loin area and touch the screen at the outer edge of loin eye muscle and don't release the screen. Draw a contour along the outer edge of the loin eye back to the original starting point. The result will be displayed on the left side of the image.

7- Measurements Accuracy

This chapter addresses measurements accuracy.

Table A reports each measurement accuracy as a function of the scales (Column **ACCURACY**) and the worst case values (Column %).

TABLE A			
Mode	Calculation	Accuracy	%
2D	Distance(mm)	$\pm[1.5\% \times \text{Depth}(\text{mm}) + 0.1]\text{mm}$	± 5
	Perimeter (mm)	$\pm[6\% \times \text{Depth}(\text{mm}) + 1]\text{mm}$	± 5
	Area(mm ²)	$\pm[1.5\% \times (D_1 + D_2) \text{Depth}(\text{mm}) + 0.025\% \times \text{Depth}(\text{mm})^2 + 1]\text{mm}^2$	± 8
M-Mode split and dual screen	Distance(mm)	$\pm[1.6\% \times \text{Depth}(\text{mm}) + 0.1]\text{mm}$	± 5
	Time(s)	$\pm[1\% \times \text{Time}(\text{s}) + 0.005]\text{s}$	± 3
Doppler full screen	Inst.velocity(m/s)	$\pm[2\% \times \text{VR}(\text{m/s}) + 0.01]\text{m/s}$	± 6
	Time(s)	$\pm[1\% \times \text{Time}(\text{s}) + 0.005]\text{s}$	± 3
Doppler split and dual screen	Inst.velocity(m/s)	$\pm[2.5\% \times \text{VR}(\text{m/s}) + 0.01]\text{m/s}$	± 8
	Time	$\pm[1\% \times \text{Time}(\text{s}) + 0.005]\text{s}$	± 3

Where VR is the Doppler velocity range.

On QUAD formats the measurement accuracy could be reduced up to the half in comparison with the same measure done on the same image displayed in full screen.

Note

If angle correction is used, a 0.1% of computation error must be added to the accuracy of the Doppler measurements.

Worst case values are calculated with the following assumptions:

- Measurement values equal to one third of the analysis depth (ex.: with a depth of 18 cm., a distance measurement of 6 cm.).
- Ultrasound speed to be constant at 1540 m/s.

Derived Data

Derived data accuracy can be calculated through the law of error propagation; worst case accuracy, based on the above mentioned assumptions, are reported together with the formulas in the following chapters.

To minimize the measurement uncertainty:

- optimize image quality.
- whenever possible, use the zoom function for maximum resolution.
- optimize the probe alignment with the Doppler flow.

ARCHIVING SECTION

This section explains how to use the archiving services and is organized as follows:

- Chapter 1: Digital Archiving
This chapter describes the digital archiving characteristics and how exams are archived.
- Chapter 2: Review of Archived Exams
This chapter explains how to review archived exams.
- Chapter 3: Archive Menu
This chapter explains how to use the menus of the Archive icons.
- Chapter 4: MyLabDesk
This chapter explains what MyLabDesk is how to use it.
- Appendix A: MyLabDesk PC Requirements, Calculation Packages and Advanced Tools
This chapter details the suggested characteristics for the MyLabDesk PC and the available calculations packages and advanced tools.

Table of Contents

1 - DIGITAL ARCHIVING	1-1
Characteristics	1-1
Archive Icons	1-1
Archiving Exams	1-2
Current exam	1-3
Import patient data	1-4
Archive	1-4
Archiving Media	1-8
Writable CDs (CD R)	1-8
Rewritable CDs (CD RW)	1-8
Double Layers DVDs	1-8
Writable DVDs	1-8
USB Keys	1-8
Wireless connectivity for networks	1-8
2 - REVIEW OF ARCHIVED EXAMS	2-1
Access to the Archive	2-1
How to Select an Exam	2-3
How to Review an Archived Exam	2-4
3 - ARCHIVE MENU	3-1
Hard Disk	3-1
Burner	3-2
Operations	3-2
Erase Device	3-3
USB Medium	3-4
Export logfile to USB	3-4
Operations	3-4
Erase USB Device	3-5
Safely remove USB Device	3-6
Network	3-6
Operations	3-7
DICOM Functions	3-7
Operations	3-8
4 - MYLABDESK	4-1
MyLabDesk Description	4-1
MyLabDesk Installation	4-2
Language	4-2
How to use MyLabDesk	4-3
Short Cut Keys	4-3
Complex Measurements	4-3
Setup Menu Key	4-3

DICOM Configuration	4-4
General Configuration	4-4
Image Review Configuration	4-4
QUIT and e Keys	4-4
Navigation	4-4
APPENDIX A - PC REQUIREMENTS	A-1
PC Requirements	A-1

1 - Digital Archiving

This chapter describes the digital archiving characteristics and how to archive exams.

Characteristics








The system is equipped with an internal hard disk (local archive) onto which exams can be archived. Data can also be stored on the external supports, listed hereinafter.

Data can be archived both in native format and in DICOM¹ format (for systems equipped with a Dicom licence) and exported as BMP, PNG, JPEG and AVI files (codec Microsoft[®] MPEG-4 V2 or MS-Video1 compatible). Exported data cannot be reviewed by the system.

Refer to chapter 2 of the System Configuration Manual for the optimal format choice

Archive Icons

The icons identifying the archiving media are displayed on the right side of the header bar when activated.

Hard Disk	Burner	USB medium	Wireless connectivity	Wireless connectivity	Network	DICOM
						



A marked icon indicates that there are problems in the management of that specific archiving medium.



A crossed icon indicates that a specific archive medium is inactive or not connected.

While data is being saved, the icon corresponding to the destination medium is flashing.

¹ DICOM is a trademark of NEMA (National Electrical Manufacturers Association); Microsoft is a trademark of Microsoft Corporation.

Note


Do not switch the system off or remove the archiving medium while saving; this could cause damages to data or to the hard disk.

Before removing the archiving medium check the operation status: the medium can be removed only when the operation shows a “Completed” status.

Archiving Exams

During the exam, still images and clips (in 2D and CFM formats) are saved into the system’s hard disk. Still images can be saved with full (BMP format) or compressed resolution (PNG and JPEG formats); clips are compressed into MPEG format, with a minimum loss of information. For more information on the codecs used for video compression please check the System configuration part of this manual


The **PATIENT & ARCHIVE** tab displays the following window:



The three sections of the Patient & Archive screen are explained below:

Current exam

An exam with all the data can be archived and stored on the system after pressing the **CLOSE EXAM** key. Each exam will be stored with a unique ID. Pressing the **CLOSE EXAM** button will clear out all data in the patient exam fields and a new exam can be initialized. Various archive options can be chosen and those options will be visible when the **CLOSE EXAM** key is pressed.



The **RESET** key resets all the data that is loaded or put in the editable fields.



Pressing **PREVIOUS PATIENT** key reloads the name and data of the previous patient, but will not re-open the previous exam.




Pressing the **EXPORTS** key offers the following options:



The **EXPORT** key allows the images or sequences to be exported into an external format (in BMP, PNG or JPEG for single frames and AVI for clips).



The **EXPORT.CSV** key allows the calculated values from one or more exams to be exported into a structured format.



The **COPY** key allows the user to copy the exams in native format on an external support.

ARCHIVING



The **DICOM** key allows archiving of data in DICOM format.

Note

The **DICOM** key is only available when the DICOM option is enabled in the system.



Pressing the **EXPORTS** key allows the user to return to the previous menu.

Note

Data can be made anonymous for Export and DICOM by the Archive & Export page of the System Setup in the menu.



Import patient data

Reimports patient data from the database into the patient fields of the current exam.

Note

The **FROM ARCHIVE** key only appears when an exam from the list is selected and highlighted AND when not more than two exams are selected by Query.




The **FROM WORKLIST** key is used to retrieve information about patient data from the DICOM server.



Archive

The **QUERY** key activates the search criteria for exams. A result will be visible if the criteria matches the data input. Pressing the **QUERY** key offers the following options:



Pressing the **TODAY** key will display a list of exams performed today.




Pressing the **YESTERDAY** key will display a list of exams performed yesterday.



Pressing the **BY FIELDS** key will display exams as a result of the search performed in the database according to data entered in the fields (Patient Name, ID, Accession Number, Diagnosis, Referring and Performing Physician, Birthdate, Age, Weight and Height).





Pressing the **ALL** key will display all exams in the archive.



Pressing the **QUERY** key allows the user to return to the previous menu.



The **SELECT** key activates the possibility for multiple selection of exams. Pressing the **SELECT** key offers the following options:



Pressing the **ALL** key will select all exams in the archive.



Pressing the **MULTI** key allows multiple selection of exams by selecting a set of exams one by one in the archive.




Pressing the **BLOCK** key allows multiple selection of exams by selecting a first and last exam to highlight a set in the archive.



Pressing the **SINGLE** key will select only one exam in the archive.



Pressing the **SELECT** key allows the user to return to the previous menu.



Select the exam(s) and press the **DELETE** key to delete the selected archived exam(s).




The **REVIEW** key enables reviewing of the stored images and sequences from an archived exam.

Note

The **REVIEW** key only appears when an exam from the list is selected and highlighted AND when not more than two exams are selected by Query.



Pressing the **EXPORTS** key offers the following options:



The **EXPORT** key allows the images or sequences from an archived exam (selected in the list) to be exported into an external format (in BMP, PNG or JPEG for single frames and AVI for clips).



The **COPY** key allows the user to copy the selected exams in native format on an external support.



The **DICOM** key allows archiving of data in DICOM format.

Note

The **DICOM** key is only available when the DICOM option is enabled in the system.



Pressing the **EXPORTS** key allows the user to return to the previous menu.

Note

Data can be made anonymous for Export and DICOM by the Archive & Export page of the System Setup in the menu

Note

When the system is switched on, the system archives the last exam performed, if the machine was switched off without first closing the exam underway.

Note

Make sure to perform regular back-ups of the exams stored in the database. Back-ups can be made on external media.

The following media can be selected for archiving and exporting operations:

Medium	Archiving	Exportation
Hard disk (DB)	Yes	No
CD (R and RW) (peripheral)	Yes	Yes
DVD (+R, -R, singlelayer) (peripheral)	Yes	Yes
USB medium	Yes	Yes
Network Directory	Yes	Yes
Dicom Storage SCU Server	Yes	No

Refer to chapter 2 of the System Configuration manual for more information about the export and archiving features

Selection can be made using the drop-down menus.

When the exam is archived on CD or DVD in DICOM format, the Biopacs Lite² viewer is automatically stored on CD or DVD, allowing the user to review the exams on any PC.

CAUTION

² Biopacs Lite is a DICOM viewer developed by Esaote

The system is equipped with several USB ports. Insert only one USB medium in either of the ports during data archiving or exporting.

The exported exams are organized in folders: every exam is contained in a specific folder together with its images, clips and reports.


If no option is selected, all stored data is deleted.

The patient data can be inserted any time during the examination. When the patient data is edited, the whole exam (with saved images and clips) is associated and named under that patient.

Archiving Media

Writable CDs (CD R)

Clean disks must be used. If the CD contains data, the system will not allow it to be written on and displays the following message:



Insert an empty disk and press OK to continue. Multisession is not supported.

Rewritable CDs (CD RW)

Rewritable CDs can be used to archive data as long as they are empty. Multisession is not supported. Read the following chapters for instructions on how to delete data from rewritable CDs.

Double Layers DVDs

Clean disks must be used. If the DVD contains data, the system will not allow it to be written. Multisession is not supported. Insert an empty disk and press OK to continue.

Writable DVDs

Clean disks must be used. If the DVD contains data, the system will not allow it to be written. Multisession is not supported.

USB Keys

USB keys are managed in multi-session: data can be added to those already containing data.

Wireless connectivity for networks

The system can only remotely control printers for printing in B&W and color. Data and images of exams can also be transferred to a network.

2 - Review of Archived Exams

This chapter explains how to review archived exams.

Access to the Archive

Images can be reloaded for each patient and a specific exam can be reviewed.

The **PATIENT & ARCHIVE** tab displays the following window:

The screenshot shows the MyLab software interface. At the top, it displays 'esaote MyLab', the date 'Jun-25-2011', and the time '09:59:08am'. The interface is divided into several sections:

- Patient Information:** Fields for Animal name, Owner name, Species (OVINE), ID (2011-06-29-19), Accession number, Diagnosis, Doctor's name, Performing Physician, Birthdate (mm-dd-yyyy), Age, Weight (kg), Diastolic Pressure (mmHg), Systolic (mmHg), Breed, and Gender (M, F, U). There are also radio buttons for Neutered (Yes, No, Unknown).
- Current exam:** A section with buttons for CLOSE EXAM, RESET, PREVIOUS PATIENT, EXPORTS, and LOG OFF.
- Import patient data:** Buttons for FROM ARCHIVE and FROM WORKLIST.
- Archive:** Buttons for QUERY (TODAY), SELECT (SINGLE), DELETE, REVIEW, and EXPORTS.
- Exam List Table:** A table with columns: Animal name, Owner name, ID, Exam Date, Exam Time, and Application. It contains four rows of data.

Animal name	Owner name	ID	Exam Date	Exam Time	Application
		2011-06-29-17	Jun-29-2011	08:53am	REPRD
		2011-06-29-15	Jun-29-2011	08:51am	REPRD
		2011-06-29-04	Jun-29-2011	08:21am	KLPRD
		2011-06-29-04	Jun-29-2011	08:13am	REPRD

The bottom part of the window lists all the exams performed and saved in the system.

Searching exams

A search can be performed through the list in order to display a specific exam or a limited list of exams matching defined criteria keys.

Below is the section of the **PATIENT & ARCHIVE** window enclosing patient and exam data:

*Searching by
Patient/exam data*

To define a search by patient/exam data:

- Select a field in the patient/exam area.
- The alphanumeric keyboard will then appear on screen.
- Enter the patient name or ID or other parameters to use as search key.
- Press the **QUERY** key and then choose the **BY FIELDS** key..

The exams matching with the search data key will be displayed in the bottom area.

Animal name	Owner name	ID	Exam Date	Exam Time	Application
test		2011-11-30-04	30-Nov-2011	14:51	ABDOMINAL

Archive

FROM WORKLIST

QUERY (BY FIELDS) SELECT (ALL)

The selection of the exam list can be performed by the **TODAY** and **YESTERDAY** keys displaying the list of exams performed the same day or the day before.

Pressing the **ALL** key will display all exams in the archive.



The **RESET** key resets the complete search procedure; the set of search criteria in the editable fields



After selecting an exam in the list, pressing the **REVIEW** key enables the reviewing of the stored images and sequences and deletion of any files that are not required.

How to Select an Exam

Archived exams are listed in alphabetic order (if **QUERY (ALL)** is selected). Each line in the exam list encloses the following data:

- Animal name
- Owner name
- Birthdate
- Exam number (ID)
- Exam date
- Application

The exams to be reviewed can be selected by touching the line in the exam list. By pressing the button **SELECT** more buttons will appear:

- **ALL** will select all the exams in the list
- **MULTI** will make it possible to select multiple exams
- **BLOCK** will make it possible to select multiple exams in blocks (first en last in the list that has been selected and everything in between)
- **SINGLE**: selection of one exam is possible


After selecting it is possible to:

- Review it
- Copy it/them in proprietary format on external memory media
- Export it/them on external memory media
- Delete it/them

How to Review an Archived Exam

Once the exam to be reviewed is selected press **REVIEW**.

The system automatically shows the first image of the exam on the screen and the list of thumbnails corresponding to the other saved images.



By selecting the thumbnails, the correspondent images or clips will be opened and displayed on screen.

*See the
“Calculations”
section for taking
measurements.*

Clips are played in motion: The **STOP** key de-activates the kinetic presentation and allows the sequence to be scrolled image-by-image, using the memory-scrolling bar. The **BEGIN** and **END** keys position the scroll memory cursor at the begin or end of the selected image or sequence.

Reloaded images can also be printed.

How to Delete an Image

To delete a stored image or sequence, select the image or sequence and press the **DELETE** key.

Pressing the **BACK** button will return to the Patient & Archive environment.


3 - Archive Menu

This chapter explains how to use the archive icon dialogue windows. The archive icons are located on the top right corner of the screen.

Hard Disk



By pressing on the icon automatically displays the space still available in memory. A popup menu will appear like the following:



The same icon graphically indicates the free space available in the internal hard disk as shown in the table below:

More than 20% of free space	Less than 20% of free space

Note

When the free space is lower than 20%, make a copy of the archive and then delete all copied exams to free space on the hard disk.

WARNING

Do not switch the unit off while executing this procedure. The hard disk could be permanently damaged.

Burner



Pressing the burner icon displays the following menu:



Operations

This option displays the operations status. The dialogue window displays the list of exams (in the Details column), the operation status (if completed, underway or failed), the date and time of the operation and the type of operation (DICOM or Export). This window can be displayed by pressing the **OPERATIONS** key, the following will be visible:

Details	Status	Date	Type	Destination details
max	Done	15 Oct 2009 10:16	Export	
kiwi	Done	15 Oct 2009 10:15	Export	

Buttons: Show all, Details, Ok




When one or more operations have been unsuccessful, the icon is marked with an exclamation mark. The dialogue window shows which operations have failed.

Note


Whenever there is a failed operation status, the **DETAIL** key will be activated. This function provides more information according to the type of failure.

Erase Device

This option is used to delete data stored on rewritable CDs. Insert the CD in the burner, select **ERASE USB DEVICE** to begin the procedure. The following message appears:



Press the **YES** to continue the procedure. While erasing the files from the CD the following message will appear:



WARNING

Do not switch the unit off while executing this procedure. The hard disk could be permanently damaged.

USB Medium




When a USB device is connected to the system a USB icon is displayed on the top right part of the screen. Touching the icon will open the following dialog window:



Export logfile to USB

This option allows the user to save the logfiles onto a USB key. To save the logfiles, insert a USB key into one of the two connectors and activate the procedure.

While exporting the logfiles onto a USB key the following message will appear:



WARNING

Do not switch the unit off while executing this procedure. The hard disk could be permanently damaged.

Operations

This option displays the operations status. The dialogue window displays the list of exams (in the Details column), the operation status (if completed, underway or failed), the date and time of the operation and the type of operation (DICOM or Export). This window can be displayed by pressing the **OPERATIONS** key, the following will be visible:

Details	Status	Date	Type	Destination details
max	Done	15-Oct-2009 10:16	Export	
klwl	Done	15-Oct-2009 10:15	Export	


Buttons: Show all, Details, Ok

Note


Whenever there is a failed operation status, the **DETAIL** key will be activated. This function provides more information according to the type of failure.

Erase USB Device

This option is used to delete stored data. Insert the USB, select **ERASE USB DEVICE** to begin the procedure. The following message appears:



Press **YES** to continue the procedure. While erasing the files from the USB key the following message will appear:



WARNING

Do not switch the unit off while executing this procedure. The hard disk could be permanently damaged.

Safely remove USB Device

This option is used to safely remove the USB device from the system.

WARNING


Do not unplug USB media during the system operation. Please make use of the safety removal procedure or unplug the USB device when the system is off.

Network



Network Icon

Pressing the network icon will open the following dialog window:



Operations

This option allows the user to select the desired network directory and to check its status, its date, the type of the network and the destination details. Pressing its **OPERATIONS** will give the following:


Details	Status	Date	Type	Destination details
max	Done	15-Oct-2009 10:16	Export	
kiwi	Done	15 Oct 2009 10:15	Export	

Note

Whenever there is a failed operation status, the **DETAIL** key will be activated. This function provides more information according to the type of failure.

DICOM Functions

Pressing the DICOM icon will open the following dialog window:



See the “System Configuration” section for DICOM configuration

The same instructions as those given for the network also apply to this option.

Operations

The option allows to check the operation status. The window displays the patients list, the operation status (completed, in progress or failed), the date and time of the operation. This window can be displayed by pressing the **OPERATIONS** key:

Details	Status	Date	Type	Destination details
max	Done	15-Oct-2009 10:16	Export	
kiwi	Done	15 Oct 2009 10:15	Export	

Buttons: Show all, Details, Ok

Note

Whenever there is a failed operation status, the **DETAIL** key will be activated. This function provides more information according to the type of failure.

4 - MyLabDesk

This chapter explains how to install and use the **MyLabDesk** utility.

MyLabDesk Description

MyLabDesk, once installed on a PC¹, reproduces the **MyLab** working environment: its working procedures are equivalent to what is available on **MyLab**.



MyLabDesk offers the **MyLab** major features for exam management: exams can be archived, copied and imported in native format; patient data can be modified; measurements, annotations, body marks can be activated; images and reports can be printed and advanced tools such as Stress Echo can be activated.

Note

The PC mouse works as a cursor in **MyLabDesk**.

¹ Appendix A lists the PC characteristics and the available advanced calculations and tools.

MyLabDesk Installation

MyLabDesk Set Up is organized in two folders: the “Archive” folder, containing copied exams, and the “MyLabDeskSetUp” folder with the installation files.

Note

MyLabDesk installation is reserved to users with an Administrator profile.

Refer to the Appendix A for the PC characteristics. Insert the medium containing the set up in the PC, select the “MyLabDeskSetUp” folder out of the File Management utility, copy the folder into a local disk and run the setup.exe file. The installation is guided by a wizard: follow the given instructions to successfully complete the installation.


Note

Make sure to select Vet during the installation, otherwise the archived Vet images cannot be loaded.

Once the installation is over, the desktop will include the **MyLabDesk** icon

To copy the exam archive into the PC, follow the following procedure:

- Run **MyLabDesk** by double clicking on its icon.
- Select the icon of the source medium (CD/DVD, USB, Network).



- The system lists all exams available in the medium: using the cursor select the exams to be locally imported.
- Press the **COPY** key and select “Local Archive”
- Press OK to start.

Language

MyLabDesk is available in the following languages: Italian, English, French, German and Spanish. The default language is English.

Refer to the next paragraphs for further information.

How to use MyLabDesk

With the exception of the **SETUPMENU**, **e** and **QUIT** keys, **MyLabDesk** keys work as the corresponding **MyLab** commands: **ANNOT** activates the annotation menu, **IMAGE** saves single frames etc.

Short Cut Keys

In some cases the cursor is controlled by the active function (e.g. when tracing distances, area or profiles in the measurements environment, or selecting words in the Annotations function.). In these situations the short cut keys can be used to activate software and control panel keys.

Below each control panel key there's the indication of the correspondent short cut (ALT + *l*). Each software key has two short cuts (Fn and SHIFT + Fn) while each software button has only one (Fn).



Right click to interrupt the active operation and free the cursor.

Complex Measurements

Some measurements of the cardiac calculations package require the selection of a different view or a different modality. Before starting complex measurements, make sure that the available images allow the user to complete them.

Measurements are system guided: operating instructions are given on the lower part of the screen.

Procedure

- Select the desired image.
- Press the **MEASURE** key and select the group.
- Follow the instructions to perform the first set of measurements.
- When requested, select the next image.
- Press **MEASURE** again to proceed with measurements.
- Repeat the procedure to complete measurements.

Setup Menu Key

This key displays the **MyLabDesk** configuration menu. The menu allows the user to configure:

- Reports
- Application Measurements



- Generic Measurements
- Glossary
- Peripherals
- Network Directory
- DICOM Configuration
- Export Settings
- Tools Preset
- General Configuration
- Image Review Configuration

DICOM Configuration

The menu includes three folders: General, Quality and Report. The first folder allows the user to assign the AE Title and set the forwarding modalities of stress echo views. Images characteristics and report forwarding modalities are defined in the other two folders.

General Configuration

The menu allows the user to set language, date and time format, height and weight format and center ID. The same menu allows the user to select the environment, human or vet.

Image Review Configuration

Brightness and contrast of the images displayed on screen can be adjusted through this option.

Refer to the “System Configuration” section for further details on the single options.

QUIT and *e* Keys



e displays the **MyLabDesk** version and **QUIT** closes the program.

Navigation

All **MyLabDesk** features are available both on local exams and on the ones archived on external media. Select the pertaining icon to access the remote archives.

The **COPY** key copies into a local folder exams archived in native format on external media.



Appendix A - PC Requirements

PC Requirements

The suggested PC configuration includes the following requirements:

- Windows® XP (Home or Professional), Windows Vista™, Windows® 7 or Windows® 8
- CPU Pentium IV 1.5 GHz MHz or higher
- HD: at least 2 GB of free space
- RAM: minimum 512 MB, suggested 2 GB
- Video board minimum 128 MB supporting 32 bit true color
- Video resolution: minimum 1024 x 768 (128MB video board suggested), suggested 1280 x 1024 (512MB video board suggested)

With a 1280 x 768 video resolution, the control panel is minimized. To maximize it and be able to access its commands, place the cursor on the icon and press **ENTER**.



Control Panel Icon



Minimizing icon

The minimizing icon reduces **MyLabDesk** to icon

Note

MyLabDesk installation can only be performed by users with an Administrator profile.

To check the active user profile, access the “User Account” utility in the PC Control Panel.

SYSTEM CONFIGURATION SECTION

This section explains the options available in the System Configuration menu. The section is organized as follows:

- Chapter 1: System Menu
This chapter describes the setup menu of the **MyLab** and the configuration options available in the system.
- Chapter 2: System Setup
The chapter includes display options for the user interface and the calibration of the touch screen. Modifications can be made regarding the measurement units and the generic aspects, as well as the archive and export services.
- Chapter 3: Peripherals
This chapter explains how to set the peripheral units.
- Chapter 4: Measurement & Analysis
This chapter explains how to customize the calculation packages of each application.
- Chapter 5: Network Configuration
This page encloses settings for system configuration in the network and configuration of network directory and wireless.
- Chapter 6: Tools
This option allows modifying and saving both the gray map selected in the active application and the settings of the Protocol Editor.
- Chapter 7: Security
This chapter explains how to set the list of users allowed to access the system.
- Chapter 8: Miscellaneous

This chapter encloses information about the system configuration. Setup of the body marks and annotations is also available in this setup environment.

- Chapter 9: Export Settings

This chapter explains which formats can be set to export the exams.

- Chapter 10: Peripherals Controls

This chapter explains the controls on the exterior of the **MyLab**. It also takes into account the buttons of peripheral devices.

Table of Contents

1 - SYSTEM MENU	1-1
System Menu	1-1
2 - SYSTEM SETUP	2-1
System Setup Menu	2-1
Display	2-2
Measurements Units	2-7
Generic Measurements	2-7
Touchscreen Calibration	2-9
Calibrate system rotation	2-10
Archive & Export	2-11
3 - PERIPHERALS	3-1
DICOM	3-1
General	3-1
Worklist	3-2
Store	3-3
DICOM Printer Configuration	3-5
DICOM Printer Management	3-7
Report	3-8
Quality	3-8
Digital Printer	3-9
Image	3-9
Digital Printer Management	3-11
Report	3-11
Wireless Connectivity	3-13
4 - MEASUREMENTS & ANALYSIS	1
Activating the Configuration Menu	1
5 - NETWORK CONFIGURATION	5-1
	5-1
Network General	5-2
Network Drive	5-3
Network Drive Configuration	5-3
Wireless	5-5
Shared Directory	5-7
6 - TOOLS	6-1
Protocol Editor	6-1
Protocol Editor Configuration	6-2
Starting the protocol	6-5
7 - SECURITY	1
Users Accounts	1

Security Access to the System	1
Security Configuration	3
Security settings	3
Change password	5
8 - MISCELLANEOUS	8-1
System Administration	8-1
System Configuration	8-1
System Options	8-2
Backup / Restore	8-3
Body Marks	8-7
Generic	8-7
Groups Setup	8-8
Default Pattern	8-9
Annotation Setup	8-9
Generic	8-9
Libraries Setup	8-10
Default	8-12
9 - EXPORT SETTINGS	1
Images and Clips Formats	1
Clips Formats	1
Image Formats	1
Report	1
10 - PERIPHERAL CONTROLS	1
MyLab Handle-Keys	1
Linear Probe	2
Footswitch	2

1 - System Menu

This chapter describes the setup menu of the **MyLab** and the configuration options available in the system.

System Menu

Under the **SETTINGS** tab the **MENU** key is available and it gives access to the system menu.

A dark green rectangular button with the word "MENU" in white capital letters.

The setup menu is organized by a tree structure composed of main options roots listed on the left of the screen; selecting a root shows the sub-menu pages enclosed and available.

Selecting one of the options opens a configuration page enclosing all the settings for customization.

The new parameters will be operative immediately.

**Setting
Parameters**

- Choose the parameters to be set.
- Select an option from the window menu and press **OK** to confirm.
- If input is needed on a field, press on the field and use the alphanumeric keyboard that appears.
- Pressing the **CANCEL** key closes the setup window (without confirming) and goes back to the main configuration page.


Note

After changing one of the settings, wait 5 seconds before switching the system off in order for the new parameters to take effect.


2 - System Setup

This chapter lists the options that are available in **SYSTEM SETUP** and explains how to configure the parameters.

System Setup Menu

A dark green rounded rectangular button with the text "System setup" in white.

When pressing the **SYSTEM SETUP** key the following screen is visible:



The options available under System Setup are:


- Display
- Measurements Units
- Generic Measurements
- Calibrate Touchscreen

- Calibrate System Rotation
- Archive & Export

Display

The display environment is organized with four tabs: **GENERAL**, **SCREEN**, **ORIENT** and **MYLIBRARY**, **PHYSICIANS** and **REPORT**.

General




General

Parameter	Description
Hospital	Editable field for the hospital name.
Actual date	This field sets the date.
Actual Time	This field sets the time.
Date Format	This setting allows the user to define the format of the date as it is displayed on the screen.
Time Format	This setting allows the user to define the format of the time as it is displayed on the screen.
Transducer Timeout	This setting allows the user to define the time for the screensaver to be activated.
Language	This setting allows the user to define the language of the user interface.
Keyboard	This setting allows the user to define the layout of keyboard.


Screen

Within this menu the appearance of information on the screen can be customized.



Parameter	Description
Settings Autohide	The Settings controls can be set to autohide or be always visible on the screen.
Display Real Time Image Parameters	Enables display of image parameters on the screen.
Display Framerate	Enables display of framerate on the screen.
Display Transducer Name	Enables display of transducer name on the screen.
Display Date & Time	Enables display of date and time on the screen.
Display Logo	Enables display of Esaote logo on the screen.
Display Coordinates	Enables display of x and y coordinates.
Key Click Sound	Allows the user to define the sound when touching the controls on the screen.

Orient




The display orientation can be adjusted according to the user's wishes. An activated 'follow system orientation' box will make sure when turning the system, the screen will turn accordingly. Orientations whose checkboxes are disabled will not be available on the screen when turning the system.

MyLibrary


MyLibrary can be enabled by clicking the box 'Yes' showing the MyLibrary tab in the control tabs on top of the screen.

Note

The Library is only available when the corresponding license is activated.




Physicians



With this menu users have the possibility to add names of physicians that can be recalled by pull-down menu in the **PATIENT & ARCHIVE** menu.


Report



In the **REPORT** menu several properties of the report can be personalized.

The first menu is the **HEADER** menu. On this screen the user can modify the name of the Department executing the examination. Another option on this screen is to add a logo to the report. The **LOAD** button opens a screen where the user can select a logo image file, for instance from an USB memory stick. The **CLEAR** button removes the added logo file. The size of the logo can be adjusted with the **SIZE** button or it can be set to auto size.


The second menu is the **COMMENTS** menu.



On this screen the user can define comments which will be available as a checkbox in the **COMMENTS** field of the report. The user can select the application, select a comment number and edit the text.

Measurements Units

This option allows the user to set the preferred measurement units.



Parameter	Description
Patient Height/Weight	Defines the units of height and weight of the patient.
Velocity	Sets the unit for velocity measurements.
Distance	Sets the unit for distance measurements.
Volume	Sets the unit for volume measurements.
Area	Sets the unit for area measurements.

Distance and Area have available an Automatic option. This option sets the unit automatically to mm and cm or mm² and cm² according to the following criteria:

- If the distance is less than 1 cm, then mm is adopted as unit; otherwise cm is used.


The area unit is set to the same unit as used for the distance.

Generic Measurements

This option allows the user to program the generic measurements available for each application. The menu is organized in two tabs: **2D MEASUREMENTS** (listing generic measurements available in 2D mode), **TAM MEASUREMENTS** (for M-Mode and Doppler).

2D Measurements

Under the **2D MEASUREMENTS** tab the following view is shown:



The 2D Measurements configuration page allows measurements to be included or excluded from this list in 2D mode.

The Autorepeat option enables the auto-repeat function for the generic measurements. This means that, when a measurement is completed, the same type of measurement is automatically enabled to be performed again.

The Default Measurements area allows to set the starting measurement if the **CALIPER** key is pressed in B-mode and Color Mode.

Note

The Default Measurements list includes also the option 'None'. When 'None' is selected, pressing the **CALIPER** key displays the list of measurements without enabling any measurement automatically.

TAM Measurements


The TAM Measurements configuration page selects measurements to be included or excluded in M-Mode and Doppler-mode.

Configuration procedure

- Enables or disables measurements from the measurement list in TAM modes.
- Set the default measure. This measure will be active as soon as the **CALIPER** key is pressed. When choosing the “none” option, no

measure will be automatically enabled when the **CALIPER** key is pressed.


The Autorepeat option enables the auto-repeat function. This means that, when a measurement is completed, the same type of measurement is automatically enabled to be performed again.



The Trace Type area configures whether the Doppler spectrum has to be traced manually (following the spectrum profile on one cardiac cycle) or automatically (by setting the beginning and the end on the spectrum and let the system trace the profile).

Touchscreen Calibration

This feature allows the user to re-calibrate the touchscreen.




This operation is very important and delicate. Please perform the operation with the touch screen stylus and follow the instructions as given on the screen by the system.

Calibrate system rotation

In this menu the tilting sensor can be calibrated. Follow the instructions as displayed on the screen.




When the tasks on the screen are correctly executed, green signs will appear as shown in the picture below:



Archive & Export


This option allows the user to enable or disable export and archiving features in the system.



- ‘Actions enabled’ allows to automatically execute archive actions as specified in the Archive & Export menu of the System Setup. When enabled the warning “Exam not Closed” page will not appear when pressing the **CLOSE EXAM** button and the specified archive actions will be executed.
- Local Archive activates the capability to archive exams in the local database of the system.
- Copy To enables the system to copy the exam in proprietary format into an external memory medium.
The option Copy MyLab Desk Setup, when enabled, copies the executable file into the same memory medium for the installation of MyLab Desk on an external PC.
- Report To activates the capability to export the report into an external memory medium.
- Dicom To enables the export in DICOM format of the exam.
- Export Images/ Clips To activates the capability to export images and clips into a memory medium in commercial formats (BMP, AVI, etc).
- Image Export Format allows to choose the image format (BMP, PNG or JPEG).
- Clips Export Format allows to choose the clip format (Codec MS-PEG4 V2, Codec MS-VIDEO1).
- Anonymize Export/Dicom allows to remove the patient data from the images and clips that are exported.

3 - Peripherals

Peripherals of **MyLab** can be configured accessing the Peripherals setup environment. When pressing the **PERIPHERALS** key the following screen is visible:




DICOM

Selecting this option allows the configuration of **MyLab** in DICOM networks and DICOM functions. DICOM option is available if the DICOM license is enabled.

The configuration menu is organized in six tabs: **GENERAL**, **WORKLIST**, **STORE**, **PRINT**, **REPORT** and **QUALITY**. Pressing on the tabs will access the sub-menu.

General

The option sets the **MyLab** AE Title. This field encloses the Application Entity Title of the scanner. The factory setting is "MYLAB". Also the user has the option to add a header with patient data to DICOM exports.




OK and CANCEL keys

In the sub-menu, the **OK** key exits from the menu saving the settings while the **CANCEL** key without saving.

Worklist

The DICOM **WORKLIST** functions as an exchange server for patient data. To retrieve this data, activate the worklist by clicking the ON button.



Procedure

Add the DICOM server of the facility or edit servers that are already in the list. By pressing the button **ADD** and completing the form below.

The form of adding the WORKLIST server is displayed below.

The screenshot shows a web-based configuration window titled "Peripherals > DICOM > Add Worklist Server". The form contains the following fields and controls:

- Worklist Server Name:** A text input field.
- AE Title:** A text input field.
- Hostname / IP Address:** A text input field.
- Port Number:** A numeric input field with a value of "104" and minus/plus navigation buttons.
- Buttons:** On the right side, there are four buttons: "Ok", "Cancel", "Ping", and "Test Connection".
- Content Area:** A large, empty light-green rectangular area, likely intended for a list of servers or a preview.

Store

This function enables or disables the possibility to export images to DICOM servers via network. The DICOM class is used only when the Enable option is 'On'.

The screenshot shows the "Store" configuration tab within the "Peripherals > DICOM" configuration window. The interface includes:

- Navigation Tabs:** "General", "Worklist", "Store" (active), "Print", "Report", and "Quality".
- Enable Control:** A toggle switch labeled "Enable" with "On" and "Off" options. The "Off" option is currently selected.
- Servers List:** A vertical list box titled "Servers:" which is currently empty.
- Action Buttons:** "Add", "Delete", and "Edit" buttons are located at the bottom of the main configuration area.
- Confirmation Buttons:** "Ok" and "Cancel" buttons are located on the right side of the window.

The DICOM class is active only when the Enable option is 'On'.


The configuration menu of these DICOM classes allows the user to add (ADD option), delete (DELETE option) one Storage server, or to modify (EDIT option) its parameters.

To add a new Storage server press **ADD** button and compile the following form.

Procedure

- Press the **ADD** button.

The Storage server configuration window is displayed:



Parameter	Description
Store Server Name	Sets the name of the Storage server.
AE Title	Sets the Application Entity Title of the scanner.
Hostname / IP Address	Network Path or IP address of the shared folder.
Port Number	Sets the port number of the Storage server.
PING	This key is used to check whether a particular DICOM server is reachable across the IP network.
Test Connection	This key is used to check the network connection to the store server.


Note

To use DICOM functions, a static IP address is recommended.

Please refer to the local network administrator for network data of the DICOM modalities and of **MyLab**.

DICOM Printer Configuration

The option allows the user to set the DICOM printer model and its printing configuration.



Preset	Action
Printers	This field lists all the printers configured in the system.
Enable	Enables or disables the possibility to print images by DICOM printers via network.
Layout	Selects the layout of the printout.
Film Orientation	Sets whether printing in portrait mode or landscape mode.
Film Size	Selects the film to print (if more than one is available for the selected printer).
Medium Type	Selects the type of medium to print on (clear film, blue film, paper).
Color Capability	Sets whether printing in color mode or grayscale mode.
Print Key selection	Sets the button on the keyboard to use for printing.
Number of copies	Sets the number of copies to be made.


The configuration menu of the DICOM printer allows the user to add (**ADD** key), delete (**DELETE** key) one DICOM printer, or to modify (**EDIT** key) its parameters.

To add a new Printer server press **ADD** button and compile the following form.

Procedure

- Press the **ADD** button.

The Printer server configuration window is displayed:



Parameter	Description
Printer Name	Name of the Printer server.
AE Title	Sets the Application Entity Title of the scanner.
Hostname / IP Address	Network Path or IP address of the shared folder.
Port Number	Sets the port number of the Printer server.
Printer Model	This pull-down menu lists the printer models that can be connected to DICOM.
PING	This key is used to check whether a particular DICOM server is reachable across the IP network.
Test Connection	This key is used to check the network connection to the dicom printer.

DICOM Printers

Refer to the site www.esaote.com for the detailed list of the compatible DICOM printers.

Print Profiles

Each DICOM printer can have different print layouts. The Layout-menu lists the set printer profiles. Scroll down the list to select the desired printer profile.

The print layout depends on the selected DICOM printer. The menu allows the user to set the print layout, the orientation, the size, medium type (sheet, film...), the color scale and the number of copies.

The option of the print key selection enables the print key for the DICOM printer to be visible on the screen.

Note

The DICOM Printer function requires the DICOM license.



DICOM Print Icon


The DICOM print icon is displayed on the right upper part of the screen. Two numbers are displayed on the icon: the left number indicates the number of images sent to be printed and the right number indicates the number of images that can be printed in one sheet. In the example icon shown in the left: one image has been sent to the printer; the maximum number of images per sheet is four.



*DICOM Printer
Icon*

DICOM Printer Management

The user can access the DICOM printer menu. Press on the icon and the system displays the following menu:

**Force Print**


To print before formatting is complete, press the **FORCE PRINT** key to start printing.

Reset Page

The option cancels all the images sent to be printed: the image counter is automatically reset.

Report

In this menu the choice can be made to include a report when exchanging data via DICOM.




The following selections for the DICOM Report format are possible:

- Biopacs report. Biopacs is an Esaote product for ultrasound exam management.
- DICOM image format for other DICOM servers.
- On the same menu the user can set not to export the exam

Quality

This option allows the user to set three different compression levels both for clips and for single images (image quality). The following values can be set: high quality (minimum compression), medium and low quality (maximum compression).

The set quality is used for any DICOM archiving operation (on server or on any other medium).



Preset	Action
Clip Quality	Sets the quality (by setting the compression level) of the clips to be exported in DICOM. The higher the quality, the higher the memory required (and transfer time).
Image Quality	Sets the quality (by setting the compression level) of the images to be exported in DICOM. The higher the quality, the higher the size of the files.

Digital Printer


*For DICOM
Printer configuration
see further in this
section.*

MyLab manages both laser and inkjet USB printers. Contact ESAOTE personnel for suggested models and assistance with proper configuration of the printer or visit ESAOTE website (www.esaote.com).

To access the configuration menu, choose **DIGITAL PRINTER**. Two tabs are shown in the following overview: **IMAGE** and **REPORT**.

Image

Pressing the **IMAGE** key will give the following screen:




Field	Action
Printer	This pull-down menu lists the digital printers that can be connected to the MyLab .
Printer Type	Enables the print key for the digital printer to be visible on the screen.
Print Key Selection	Sets the number of pictures that should be printed on one page. The images are sized properly accordingly.
Images Per Sheet	Sets the printer to print in color or BW.
BW - Color	Sets the orientation of the page (portrait or landscape).
Orientation	Sets the printer to print the parts of the pictures outside the ultrasound image in white.
White Background	Includes Physician Name and/or Hospital Name in the report page to print out.
Include Header	Includes Patient Name and/or Date in the report page to print out.

Format

The following print format icons are available, depending on the number of pictures set that should be printed on one page. The settings of images per sheet is described above.

Print format icons



The format icon gradually becomes blue as images are sent to the printer.



The figure displays the print format set on 3x2 and three images (colored in blue) are sent to the printer. Printing takes place when the icon is completely filled.


The printing function is also available off line, i.e. when the printer is not physically connected to the system. In this case, printouts are temporarily stored on the hard disk. Printing is automatically activated once the printer is connected.



Digital Printer icon

Digital Printer Management

The user can open the menu to manage the digital printer. Press on the icon and the system displays the following menu:



Force Print

To print before formatting is complete, press the **FORCE PRINT** key to start printing.

Reset Page

The option cancels all images sent to be printed: the image counter is automatically reset.


Cancel all printjobs

This option cancels all printjobs instantly.

Report

The Report page refers to the configuration of the printer that is used to print reports.

The **REPORT** tab will give access to the configuration screen of the report.




Field	Action
PRINTER	Sets the printer model to be selected among the installed ones
BW - COLOR	Black/White or Color

CAUTION

Do not switch the system off until the printing stage has been completed

Wireless Connectivity

This is a protocol for exchanging data over short distances from fixed and mobile devices. Keyboards or mouse devices can be compliant with this protocol of the **MyLab** and can be connected to the system.



The list of devices recognized via the wireless protocol by **MyLab** is shown in the window.

The list can be refreshed by **REFRESH** button.

Given the list of devices, it is possible to enable the connection (by **CONNECT** button) or disable it (**DISCONNECT** button).

4 - Measurements & Analysis

This chapter explains how to customize the calculation packages of each application, available when the **MEASURE** key is pressed. This section also lists the report editors for the application packages.

Activating the Configuration Menu



To access the configuration menu, press **MENU** (under **SETTINGS** tab) and choose **MEASUREMENT & ANALYSIS**. When pressing the **MEASUREMENT & ANALYSIS** key the following screen is visible listing the applications whose packages can be configured:




For the configuration and customization of the application packages please refer to the Calculations section of this manual.

5 - Network Configuration

This menu encloses four submenus:

- Network general
- Network drive
- Wireless
- Shared directory¹

To access the configuration menu, press **MENU** (**SETTINGS** tab) and choose **NETWORK CONFIGURATION**. The following screen is then available:



¹ This option is available for the **MyLabSatVet**

Network General

This section explains the configuration settings that can be found for **LAN** and **WIFI**.

The General Network menu encloses two tabs: **LAN** and **WIFI**. Selecting the **LAN** tab will show the following window:

Note

The setup page under the **WIFI** tab is similar to the **LAN** tab.

Parameter	Description
Computer Name	Sets the name of the scanner in the network.
IP Address Subnet Mask Default Gateway DNS Server	These settings allow to configure the system in the network.
Enable DHCP	If DHCP is enabled, the system will acquire network configuration settings from the network server (if the server is enabled to provide this service). If DHCP set off, network configuration settings have to be specified by providing IP Address, Subnet Mask and Default Gateway).

Note

Please refer to the local network administrator for network data.


Network Drive

This section explains how to configure the network drives that can be used as archiving media.

Network Drive Configuration

The Network Drive option allows the user to select a shared folder in the network as destination driver for image or report export.

Once selected the option, the system displays the following window:



The Network Drive window lists the network drivers configured and available as destinations in the system.

In this environment it is possible to:

- Add a new network driver
- Delete an existing network driver
- Update the data of an existing network driver (Edit)
- Reconnect to a previous selected network

Note

The network has to be configured before setting network directory: see the “Archive” section for further details.

Please refer to the local network administrator for network data.

To add a new shared folder press **ADD** button and compile the following form.

Procedure

- Press the **ADD** button.

The network directory configuration window is displayed:

Parameter	Description
Description	Name of the shared folder.
Drive	Name of the drive associated to the folder.
Hostname/ IP Address	Network Path or IP address of the shared folder.
User Name	Username to access the shared folder through network (according to the network configuration).
Password	Password to access the shared folder through network (according to the network configuration).

Note


Please refer to the local network administrator for network data.

To update an existing shared folder:

Procedure


- Select the folder and press the **EDIT** button

Update the driver information in form as shown below.



Wireless

Pressing the **WIRELESS** key in the **NETWORK CONFIGURATION** menu will activate the following screen:



Parameter	Description
Name	Name of the network.
State	If a network is available, unavailable or new.
Signal quality	Shows the signal strength of the network.
Mode	Shows whether the connection of the network is automatic or manual.

The screen shows a table with all the available networks listed. The signal quality depends of the position of the system in the area. The **REFRESH** key can be used to restore signal strength of the available networks. The **BACK** key returns to the menu. Any actions performed on this page will NOT be canceled

The user can add a Network by pressing the **ADD** button. Pressing the **EDIT** key allows to modify an existing network.

To add a new network press **ADD** button:

Procedure

- Press the **ADD** button.


Fill in the following form

To update an existing shared folder:

Procedure

- Select the folder in the list
- Press the **EDIT** button

Modify the data in the following form



Parameter	Description
SSID	The user can enter the name of a Wireless network.
Connection Mode	The network can be set to automatic or manual connection.
Authentication	Two methods of authentication can be used: Open System authentication and Shared Key authentication, or other security protocols.
Encryption	The encryption type depends on the authentication that is chosen.

Note

Please contact the network administrator before continuing with procedure.

Shared Directory²

When the **SHARED DIRECTORY** key is pressed in the Network Configuration menu, the following screen appears:

² This option is available for the **MyLabSatVet**

System Setup > Network share

Share permission

User: MyLabUser

Password: QggHnx3E Allowed Input: A-Z, a-z, 0-9

Enabled: Yes No

Buttons: Ok, Cancel, Delete Share

Information can be shared between the **MyLabSatVet** and the **MyLabTwiceVet** with the creation of a shared directory. In order to do so, both systems need to be connected to the network.

This connection can be configured in the window Network Share:

Procedure

- Enable the connection by clicking YES
- Press **OK** in order to activate the shared directory

All data (clips, images, reports) that is saved on the **MyLabSatVet** will automatically go to this shared directory. The **MyLabTwiceVet** can access this directory via the Network Drive by entering the similar password as in the **MyLabSatVet**.

The **DELETE SHARE** key can be used to erase the data in the shared directory.

6 - Tools

This menu includes the Protocol Editor.

To access the configuration menu, press **MENU** (under **SETTINGS** tab) and choose **TOOLS**. The following screen is then available:




Protocol Editor

This section explains how to create and configure protocols in the system.

A protocol is a sequence of messages and/or actions that can be defined by a user to support procedures.

To access the Protocol Editor menu, press **SETTINGS** tab and choose **TOOLS**. Press the **PROTOCOL EDITOR** key to activate the setup environment.



Protocol Editor Configuration

The screen is organized in several editor sets.

Activation of Protocol




This section defines whether the protocol has to be associated to the application or to a specific preset in the application.

In these cases the relative selections have to be made inside the checkboxes and the pull-down menus.

- Depending on the selection made by the user, the specific selections apply. E.g. if “Only Species”, “Only Application” and “Only Preset” are all enabled, the protocol will be active for those specific species, application and preset, or if “Only Species” is enabled, the protocol will be active for that specific species, but for all applications and presets included.

At the bottom of the page, the checkbox “Auto-Start at Preset Selection” is present. This setting activates the protocol at the selection of the preset at the beginning of the procedure.

Protocol Configuration

BACK and **NEXT** buttons select the step of the protocol to define (by text and actions).

The text field allows the user to edit text to be displayed on the screen when the specific step is reached during the procedure.

Shown image field provides information over the figure to be displayed at the specific step.

The section “Action before showing text” lists all the possible actions that can be enabled during the procedure when the specific protocol step is reached.

Note

The first item cannot be allocated with an action. The first item can only be allocated with a text message, an image or a sound.

Configuration Procedure

To start editing, the user must first input text. Pressing on the text field will make the keyboard appear. Press the **NEXT** key to continue with the following item.

From this point, the user can choose to allocate text **AND/OR** an action **AND/OR** an image **AND/OR** a sound.

Action

Choose an action by highlighting it first and then press the **NEXT** key to continue.

Images

An image can be uploaded via USB by using the **LOAD IMAGE** key. Connect a USB to the system. The system will show a list of images available in the USB medium.

A preview is shown on the right of the screen. Choose an image and press the **OK** key. An example of the image is shown on the left side of the Protocol Editor screen. Text can also be added beneath the image; press on the text field and the keyboard will appear. Press the **NEXT** key to continue.

The user can scroll through all the edited items by using the **BACK** and **NEXT** key.

Advanced Functions



IMPORT and **EXPORT** buttons permit to import / export a protocol from / to another **MyLab** unit. To export a protocol, connect a USB to the system and press the **EXPORT** key. If preferred, the user can change the name by keyboard. Press **OK** to confirm.

To import protocols, connect a USB to the system and press the **IMPORT** key. Choose a protocol from the list and press **OK**.

Note

Only one protocol at time can be imported. When the protocol is assigned to the application and/or preset, the protocol will be saved. Similar allocations of protocols to applications and/or presets will be overwritten.

CLEAR button erases the protocol in the specific unit.



LOAD TEXT and **LOAD IMAGE** buttons permit to load text or an image from USB media. The text or the image will be added to protocol at the specific step.

At the end of the protocol definition, press **OK** and **CANCEL** buttons to close the protocol editor confirming or not confirming the modifications.

Starting the protocol

The protocol sequence will start as defined in the configuration (at the selection of the preset, etc) and be visible at the bottom of the screen in real time.



To move to the next step of the protocol press the **NEXT** key.

If images are edited with text, they will be displayed on the left side of the screen.

Pressing the **CLOSE** key ends the protocol.

7 - Security

Access to the system and to the archive can be reserved only to authorized users. In this case access to the system can be managed via password.

The access under password allows a secure management of the archive: its data can be reviewed and modified only by authorized personnel.

This chapter explains the archive security management offered by **MyLab** and how to define the list of the authorized users.

Users Accounts

Two different accounts are available: administrator and user.

The system administrator can decide whether to activate the security or not. When enabled, he can add, delete users and define their profiles. More administrators can be defined.

Note


The default administrator user name and password are: ADMINISTRATOR and MYLAB. Change this account if the security management is activated.

Security Access to the System

When security is enabled, a password is required to access the system. When starting up, the system requires user name and password.

Note

The system is case sensitive; please consider this when entering passwords.



When the **EMERGENCY** key is active, exams can be performed (Emergency Access 'On') without entering any user name and password. The Emergency Access allows performing exams and review images saved during the current exam, but will not allow access to the Archive (**PATIENT & ARCHIVE** key).

Security Configuration



The **SECURITY** key includes two option keys: **SECURITY SETTINGS** and **CHANGE PASSWORD**. The last folder is displayed only when the security is enabled.




Security settings

Only administrators can access this option. The configuration menu has two tabs: **GENERAL** and **USERS**.

General Tab

The following settings are available:



Field	Action
Minimum Password Length	Sets the minimum number of characters for the password (maximum 20).
Password Expiration Time	Sets the time (in days) after which the password expires.
Login Disabling Inactivity Time	Sets the inactivity time (in days) after which the account automatically expires.
Access Control	Enables/disables the security access by password.
Emergency Access	Enables / disables the emergency access.
Auto Last User Login	Enables / disables automatic access with the last user authenticated in the system.
Locking Timeout	Sets the time (in minutes) after which the system will lock, when current user didn't use the system.

Saving the configuration

The settings can be saved on USB medium when the **EXPORT** key is pressed.

Users Tab

User accounts can be added (**NEW** key), or deleted (**DELETE** key).

Note

It will be required to enter a new password at once when a new user is added.

A new user account is identified by "ID", "First name", "Middle name" and "Last name". The user can be set as administrator when the relevant box is checked.

The Status field allows to enable a certain user and defines whether a user acquires administrator rights.


The **CLEAR PASSWORD** button allows to change the password for existing users.

Users can be assigned to specific groups: Operator, Referring Physician or Diagnostic Physician. This function can be activated to display a list of users. If no allocation is preferred, choose the option 'None'.

The privileges section allows the user to be able to access certain modules, such as Printing, Archiving, Exporting and access to the Settings page.

Change password

This option is available for all authorized users and allows to change the password. Old and new passwords have to be inserted.



The screenshot shows a 'Change Password' dialog box with a teal header. It contains three text input fields labeled 'Old Password', 'New Password', and 'Retype Password'. To the right of these fields are two buttons: 'Ok' and 'Cancel'. Below the dialog is a virtual keyboard with a light green background. The keyboard includes keys for numbers 1-0, letters q-z, symbols like @, #, \$, %, ^, &, *, (,), -, +, {, }, [,], :, ", ', /, and navigation keys like Backspace, Tab, Enter, Shift, and Clear.

8 - Miscellaneous

A green rectangular button with rounded corners and a slight shadow, containing the word "Miscellaneous" in white text.

The Miscellaneous section encloses the following configuration pages:


- System Administration
- Body Marks
- Annotation Setup



System Administration


System Configuration

This option provides information about the system Hardware, Software and Firmware configuration.






System Options

This tab summarizes the status of the system licenses. The user can activate licenses by selecting them and entering the relative license code.



License codes are generated according to the Hardware ID of the system: different codes cannot activate the same license.

If a demo license has been activated, this shows the time left.

The Options page lists all the options/licenses of the **MyLab** and their status (enabled  / demo license  / disabled ). When an option is enabled, it is indicated by a green check mark.

In order to enter a new license, perform the following steps:

Procedure

- Select the option to be enabled in the Options list.
- Edit the license code in the New License field.
- Press **SUBMIT**.

If the number is correct, the status will change into a positive check and the license becomes instantly functional.

Please consider that the system is case sensitive.


Backup / Restore

This chapter explains how to backup and restore the user presets. The Backup/Restore function allows the user to save the settings and customizations on a USB memory drive and reload them on **MyLab** systems.

Backup

This option is used to back up all user settings on a USB medium. In the menu, press on the **MISCELLANEOUS** key, choose **SYSTEM ADMINISTRATION** and then the **BACKUP/RESTORE** tab.

The list of groups and types of settings that can be saved are listed in the page and can be selectively activated for the backup operation.




The settings that can be backed-up and restored are:

Preset	Action
System Settings	Backs up/restores all configurations of the “System Setup” menu options.
Application Settings	Backs up/restores all configurations of the “User Presets” menu options.
Body Marks	Backs up/restores all configurations of the “Body Marks” menu options.
Annotation	Backs up/restores all configurations of the “Annotation Setup” menu options.
M&A (Measurement & Analysis)	Backs up/restores all configurations of the “Measurement Units”, “Generic Measurements”, “Report layout” and the “Archive and Export Settings” menu options.
Patient Archive	Backs up/restores all configurations of the “Patient database” menu options.
Users	Backs up/restores all configurations of the “General Security Settings” menu options.
Protocols	Backs up/restores all configurations of the “Protocols” menu options.

Procedure

Connect the USB medium to the system and press the **BACKUP** key to start backup. Please wait till when the procedure is finished.

When the procedure is finished the following screen will be shown:




Press the **OK** key to continue.

CAUTION

The **INSTALL DEFAULT SETTINGS** key will restore all settings to factory settings with loss of all user modifications.

Restore


The backed up presets can be restored into **MyLab** units.



It is possible to restore specific user presets: the system lists the user presets that can be restored.

Connect the USB medium to the system and press the **RESTORE** key to start restoring. Please wait until the procedure is finished.

When the procedure is finished the following screen will be shown:



Press the **OK** key to continue.

CAUTION

The restored presets completely overwrite those on the system: the previous presets are then lost.


CAUTION

The **INSTALL DEFAULT SETTINGS** key will restore all settings to factory settings with loss of all user modifications.

Body Marks


The section for the body marks consists of three tabs: **GENERIC**, **GROUP SETUP** and **DEFAULT PATTERN**.

Generic



Parameter	Description
Erase when the Probe is Changed	When this setting is active, changing the active probe removes the body mark from the screen.
Auto Display	If this setting is enabled, the default body mark is displayed on the screen as soon as the image is frozen.

Groups Setup




Parameter	Description
Library Selection	Selects the body marks library.
New Library	Adds a new body marks group. In order to create a new group, edit the name in the editable field (Library Selection) and press New Library. The Pattern List on the left will be empty. It is possible now to select a pattern from the All Patterns list on the right of the screen and copy it into the Pattern List (of the new group).
Delete Library	Deletes the group selected via the pull-down menu.
Patterns List	Lists the body marks belonging to the selected group.
All Patterns	List all the body marks from all the groups in the system.
Factory	Restores the factory body marks groups.
Delete	Removes a body mark from the list.
Copy	Allows to copy a body mark from the All Patterns list on the right of the screen into the Pattern List.

A body mark can be added inside the Pattern List by selecting a cell and then pressing **COPY**.

Selecting a body mark in the Pattern List and pressing **DELETE** removes the body mark from the group.

Selecting **OK** confirms the modifications.

Default Pattern



Parameter	Description
Application	Selects the application.
Body Marks Group	Sets the group to associate to the application.
Default pattern 1/2	Sets the default body mark for the application.

Annotation Setup


The glossary is composed of several user configurable libraries. It is possible to assign a specific library to each application.

This menu is organized in three tabs: **GENERIC**, **LIBRARIES SETUP** and **DEFAULT**.

Generic

This folder allows choosing the font size. It is also possible to modify the home position of the annotation text as well as the free text.

To automatically delete the text when returning to Real Time, select the option Hide with live image 'On'.




Parameter	Description
Annotation Home Position	Sets the coordinates for the text annotation home position.
Text home position	Sets the coordinates for the home position for free text editing.
Font size	Sets the size of the text.
Hide with Live Image	When enabled, removes the annotations from the screen when the system is unfrozen.

Libraries Setup

The user can assign both a list of single words and a list of sentences to a library.

Library by Word

The configuration menu looks like the following figure



Parameter	Description
Library Selection	Selects the annotation library.
New Library	Adds a new library to the list. In order to create a new library, edit the name of the library in the editable field (Library Selection) and press New Library. The Annotation List on the left will be empty. It is possible now to select a word from the All Words list on the right of the screen and Copy it into the Annotations List (of the new library)
Delete Library	Deletes the selected library.
Annotation List	Lists the words enclosed in the library selected in the Library Selection List.
All words	Lists all the words from all the libraries.
Factory	Restores the factory libraries.
Delete	Removes a word from the list.
Copy	Allows to copy a word from the All Words list on the right of the screen into the Annotations List.

Refer to the "Software Keys" section for further details on the glossary

The Library Setup can be activated in the annotation mode By Word or By Sentence.

The annotation list in the mid part of the menu consists of a list of words to be assigned. On the right, there is a list of all the words available, inserted in the different libraries.

To insert a new word or change a word in the Annotations list, press on a word to highlight and the press the keyboard icon. **MyLab** will allow entering text through the alphanumeric keyboard. A blank box can be used as a space inside the list.

Press on the desired item to highlight

New library definition


Words available from the list on the right can be added to libraries. Highlight the desired item and press **COPY** to insert the word into the list. To delete a word, select a word until it is highlighted and press the **DELETE** key.

To create a new library, press the desired field and enter the name of the library with the alphanumeric keyboard. Press **NEW LIBRARY** to create new words following the previous instructions. **DELETE LIBRARY** key deletes the selected library.

In each menu, the **FACTORY** key allows assignment of factory settings.

Library by Sentences

The configuration menu looks like the figure below:




Refer to the "Software Keys" section for further details on glossary

The menu shows an active library with a table composed of four columns. Every row of the table is a sentence, composed of the four words listed in the four columns. At the right of the table the system displays the list of the words present in all libraries.

The procedure to create and modify a library is the same of the glossary by word previously described.

Default

The configuration menu looks like the figure below:



Parameter	Description
Species	Selects the specie.
Annotation Library	Selects the library that is associated with the application.
Automatic vocabulary	Enables the matching word suggestion system during the free text editing.

MyLab allows assignment of a preferred library to each application. The pull-down menu shows the available libraries. To make a selection just press the desired library. Both libraries by word and by sentences can be assigned to each application.

*The **ANNOT** key will automatically activate text insertion during the exam*

When **Automatic Vocabulary** is selected, **MyLab** automatically activates, during the annotation phase, the search of words through the application library. If the characters entered occur in the active library, the system will automatically propose that word. To confirm the proposed word, press **ENTER** or proceed with text insertion.

Note

The automatic Vocabulary option works with an external keyboard only.

9 - Export Settings

This chapter explains which formats can be set to export the exams.

Images and Clips Formats

The option allows to set the compression format of single images and clips and the export format of reports. The defined formats will be used each time images, clips, and reports are exported.

Clips Formats

The following formats are available:

- Codec MS-PEG4 V2, ensuring the best image quality.
- Codec MS-VIDEO1, ensuring a higher compatibility with other software programs for clips management.

Image Formats

The following formats are available:

- No compression (BMP format), ensuring the best image quality.
- Medium compression (PNG format), ensuring a good image quality.
- Compressed (JPEG format), with low image quality.

Report

The report can be exported in:

- XML format.
- Editable PDF file format.

10 - Peripheral Controls

This chapter explains the controls on the exterior of the **MyLab** and peripheral devices.

MyLab Handle-Keys

The system is equipped with four programmable buttons on the battery handle.



The buttons can be programmed by holding down a button on the handle and press any imaging related key on the screen simultaneously. The programmed key will get an icon. In this icon the corresponding button of the handle is marked.



Linear Probe

Several probes are equipped with three programmable buttons at the side. The probes can be configured to be used with image related buttons on the screen.

It can be used e.g. to manage the **FREEZE** button or to control the real-time modes.



The buttons can be programmed with image related softkeys such as **FREEZE**, **FREQUENCY**, **GAIN**, all **SETTINGS** and more. Programming is performed by holding down a button on the probe and press any control panel key on the screen simultaneously. The programmed key is indicated with a mark that is linked to the corresponding button on the probe.



Footswitch

Footswitch

A footswitch (available as an optional accessory) can be connected to **MyLab**. The footswitch can be configured to be used with image related buttons on the screen.



It can be used e.g. to manage the **FREEZE** button or to control the real-time modes. The footswitch should be connected to connector **C16** (see Chapter 3 of the Getting Started) on the side panel of the Desk Stand (or Trolley).

The user can select which control keys shall be linked to the left and right footswitches. Controls on the user interface can be associated to the keys of the footswitch.

The buttons can be programmed by holding down one of the footswitches and press any control key on the screen simultaneously. The programmed key is indicated with a mark that is linked to the corresponding footswitch.



MYLIBRARY SECTION

This section encloses information about the MyLibrary interface. The section is organized as follows:

- Chapter 1: MyLibrary

This chapter describes the concept of the MyLibrary and important warnings regarding the use of MyLibrary.

Table of Contents


1 - MYLIBRARY	1-1
Warnings	1-1
MyLibrary activation	1-2
ACCEPTANCE SCREEN	1-3
MyLibrary Organization	1-4
Warnings	1-5
Warnings	1-6

1 - MyLibrary

The MyLab product incorporates the concept of the MyLibrary. This MyLibrary function gives an 'on-board' viewing possibility showing veterinary procedures according to a reference doctor. The information is set up following a commonly used approach meaning: an image of the relevant anatomical structure, a picture showing the suggested probe positioning on the specific anatomical site, the clinical ultrasound image normally resulting from the scan in the specific anatomical site and some explanatory text. Not all components have to be available in the MyLibrary. MyLibrary is intended to be used for training, review or reminding. All presented information inside MyLibrary is made under the full authority of the reference doctor who is mentioned in the credits page. A basic level of ultrasound scanning is needed before using the information given in MyLibrary.

MyLibrary is a software environment with a limited numbers of image examples that may help the user to perform the examination because, if used properly, it may help the user to recall the correct interpretation of the images and recognize the main anatomical structures. MyLibrary must be used in accordance with the following warnings in order to avoid possible harm to the patient due to the wrong interpretation of the images or inappropriate use of the ultrasound equipment or the inappropriate performance of the procedure.

WARNINGS

 Observe the following warnings for maximum safety

Warnings

- Only a limited numbers of cases/images are represented as examples without considering all possible anatomical variance or pathological findings between different individuals.
- The procedures showed are based on the ultrasound procedures as defined by the MyLibrary reference specialist and local procedures can be different from this one. The MyLibrary reference specialist is indicated in the credit window of the MyLibrary.
- Misinterpretation of the MyLibrary images can give a risk for wrong probe anatomical location.
- Esaote underline the importance for the user to reach through adequate training courses proper skills in the use of ultrasound equipments, ultrasound image interpretation and procedures.
- A basic level of medical ultrasound scanning, ultrasound image interpretation and a basic level of ultrasound probe manipulation are needed to perform ultrasound guided procedures.

- MyLibrary cannot substitute the proper training of the user in the management of the ultrasound equipment, in the interpretation of the ultrasound images and in the performance of procedures.
- MyLibrary is a reference tool providing information and suggestions, not aimed to over-rule or modify in any case the local standard operating procedure. Esaote is not responsible for the results and/or consequences of any performed procedure.
- It is very important to select the proper anatomical site to make sure the right example images are shown.

It is recommended to the user to consider this warning appropriately in order to correctly interpret the ultrasound images and perform the examination and the procedures considering the specific conditions of the single patient under examination.

This chapter explains how to use the MyLibrary function.

MyLibrary activation

MyLibrary is an on-board reference tool integrated in the system that provides application aid and tips in the use of the system during procedures.

To activate MyLibrary, press the **MYLIBRARY** tab from the section tabs. The acceptance page will be displayed as shown below.

The acceptance page emphasizes to the user the importance of using MyLibrary in accordance with the instruction of this User Manual and MyLibrary's inherent limitations and provides to the user the following information:

MYLIBRARY IS A SOFTWARE ENVIRONMENT WITH A LIMITED NUMBERS OF IMAGE EXAMPLES THAT MAY HELP THE USER TO PERFORM THE EXAMINATION BECAUSE, IF USED PROPERLY, IT MAY HELP THE USER TO RECALL THE CORRECT INTERPRETATION OF THE IMAGES AND RECOGNIZE THE MAIN ANATOMICAL STRUCTURES. MYLIBRARY MUST BE USED IN ACCORDANCE WITH THE FOLLOWING WARNINGS IN ORDER TO AVOID POSSIBLE HARM TO THE PATIENT DUE TO THE WRONG INTERPRETATION OF THE IMAGES OR INAPPROPRIATE USE OF THE ULTRASOUND EQUIPMENT OR THE INAPPROPRIATE PERFORMANCE OF THE PROCEDURE.


- ONLY A LIMITED NUMBERS OF CASES/IMAGES ARE REPRESENTED AS EXAMPLES WITHOUT CONSIDERING ALL POSSIBLE ANATOMICAL VARIANCE OR PATHOLOGICAL FINDINGS BETWEEN DIFFERENT INDIVIDUALS

- ESAOTE UNDERLINE THE IMPORTANCE FOR THE USER TO REACH THROUGH ADEQUATE TRAINING COURSES PROPER SKILLS IN THE USE OF ULTRASOUND EQUIPMENTS, ULTRASOUND IMAGE INTERPRETATION AND PROCEDURES.
- MYLIBRARY CANNOT SUBSTITUTE THE PROPER TRAINING OF THE USER IN THE MANAGEMENT OF THE ULTRASOUND EQUIPMENT, IN THE INTERPRETATION OF THE ULTRASOUND IMAGES AND IN THE PERFORMANCE OF PROCEDURES.


IN ADDITION TO THESE WARNINGS, THE USER SHOULD READ AND FOLLOW THE INSTRUCTIONS AND WARNINGS IN THE USER MANUAL IN ORDER TO SAFELY USE MYLIBRARY AND IT IS UP TO THE USER TO CORRECTLY INTERPRET THE ULTRASOUND IMAGES AND PERFORM THE EXAMINATION AND THE PROCEDURES CONSIDERING THE SPECIFIC CONDITIONS OF THE SINGLE PATIENT UNDER EXAMINATION.

To proceed further and activate MyLibrary the user has to press the **ACCEPT** tab shown in the lower part of the acceptance page. The user should press this tab only if he/she has reached the proper level of understandings through reading this User Manual about MyLibrary tool functions and limitations and he/she is able to correctly interpret the information provided in MyLibrary.

ACCEPTANCE SCREEN



To activate MyLibrary press **ACCEPT** and the following screen will be displayed.



On the left side of the screen a list of the anatomical sites is displayed.

Note

The MyLibrary tab can be enabled in the system setup menu. Refer to chapter 2 of the system configuration manual.

MyLibrary Organization

MyLibrary consists of a database with sample images and (sub) application information.

On the left side of the screen, there is a list of anatomical sites, which the user can select. Selecting an anatomical site opens a sub-list of applications that can be selected and displayed in the rest of the page.

The page describing the selected application can include:


- An image of the anatomical structure.
- A picture with an example of scanning in the specific application, showing the suggested probe positioning on the specific anatomical site.
- The clinical ultrasound image example resulting from the scan in the specific anatomical site.

WARNINGS

Observe the following warnings for maximum safety

Warnings

- Only a limited numbers of cases/images are represented as examples without considering all possible anatomical variance or pathological findings between different individuals.
- The procedures showed are based on the ultrasound procedures as defined by the MyLibrary reference specialist and local procedures can be different from this one. The MyLibrary reference specialist is indicated in the credit window of the MyLibrary.
- Misinterpretation of the MyLibrary images can give a risk for wrong probe anatomical location.
- Esaote underline the importance for the user to reach through adequate training courses proper skills in the use of ultrasound equipments, ultrasound image interpretation and procedures.
- A basic level of medical ultrasound scanning, ultrasound image interpretation and a basic level of ultrasound probe. manipulation are needed to perform ultrasound guided.procedures.
- MyLibrary cannot substitute the proper training of the user in the management of the ultrasound equipment, in the interpretation of the ultrasound images and in the performance of procedures.
- MyLibrary is an educational tool providing information and suggestions, not aimed to over-rule or modify in any case the local standard operating procedure. Esaote is not responsible for the results and/or consequences of any performed procedure.
- It is very important to select the proper anatomical site to make sure the right example images are shown.





Overlay icon

Pressing the overlay icon on the right side of the screen will activate an overlay on the ultrasound image describing the relevant anatomical structures. The same anatomical structures are also described in the legenda above the ultrasound image.

The buttons **TRANS** and **SAG** (positioned on the left bottom part of MyLibrary environment) select the page of MyLibrary to display (transversal or sagittal, if applicable for the selected section).

Live Preview¹

The **LIVE PREVIEW** key activates the Real Time image from the scanner in the top right part of MyLibrary environment. The user can compare this way the Real Time image scan with the example ultrasound image in the Library Database. The **FREEZE** key freezes the Real Time image in the live preview window.

WARNINGS



Observe the following warnings for maximum safety

Warnings

- The Live Preview function allows visualizing the real time image inside MyLibrary to compare it with the available image examples. This viewing environment is absolutely not suited to perform any interpretation of the ultrasound images and not suited for the performance of Anaesthesiology and Musculoskeletal procedures.
- In Live Preview there is no possibility to change settings to optimize the ultrasound image. If a change of settings is needed it is mandatory to switch to the user or settings pages

Credits



Credits icon

The Credits Icon is placed on the top right part of the MyLibrary environment and is linked to the Credit Window. Pressing the icon will open the Credit Window showing the MyLibrary version and information about the development of the tool including medical references with whom MyLibrary is developed.

¹ Live Preview is not available for the USA market.



The Credit Window reports also the following warning:

MYLIBRARY IS A REFERENCE TOOL PROVIDING INFORMATION AND SUGGESTIONS, NOT AIMED TO OVER-RULE OR MODIFY IN ANY CASE THE LOCAL STANDARD OPERATING PROCEDURE. THE PROCEDURES SHOWN ARE BASED ON THE ULTRASOUND PROCEDURES AS DEFINED BY THE MYLIBRARY REFERENCE SPECIALIST AND LOCAL PROCEDURES CAN BE DIFFERENT FROM THIS ONE. ESAOTE IS NOT RESPONSIBLE FOR THE RESULTS AND/OR CONSEQUENCES OF ANY PERFORMED PROCEDURE. THE USER SHOULD READ AND FOLLOW THE WARNINGS IN THE USER MANUAL IN ORDER TO CORRECTLY USE THE MYLIBRARY TOOL.

SPECIALTY PROBES AND NEEDLE GUIDE SECTION

*Refer to Getting
Started manual for
available probes.*

This section provides information on how to interact with the Needle Guide.

The section is organized as follows:

- Chapter 1: General Information on Needle Guides

This chapter gives a list of the needle guide kits that can be used and provides general information about the correct way to use the Needle Guide Software.

- Chapter 2: Using the Needle Guide

This chapter explains how to verify the needle guide working and to use it.

Read the “Transducers and Consumables” Manual carefully for detailed information on how to properly and safely handle the probes and their Needle Guide Adaptors.

Table of Contents

1 -	GENERAL INFORMATION ON NEEDLE GUIDES	1-1
	Needle Guide Kits	1-1
2 -	USING THE NEEDLE GUIDE	2-1
	Needle Guide Selection	2-1
	Checking the Guide's Working	2-1
	Using the Needle Guide	2-2
	After the Examination	2-3

1 - General Information on Needle Guides

For **MyLab** there are a series of optional adaptors for the biopsy needle guide available. These are equipped with specific couplings that attach to the probe. **MyLab** can be used to display a needle guideline¹ throughout the ultrasound image. This chapter gives general information about the correct way to use the needle guides.

The entire needle guide program is managed through the **BIOPSY** key. The biopsy procedure is active in B-Mode, in CFM and if the active probe is compatible with an attachment kit.



The user should read and completely understand all information on probes and biopsy kits, as detailed in the “Transducers and Consumables” manual, prior to use of this software.

Needle Guide Kits

The table below lists the available adaptors:

Biopsy adaptor	Probe
WBSL33X	SL3332
CBSL33X	SL3332
IKL3323	SL3323
IKC3421	SC3421
CBSC23X	SC3123

Note

All the listed kits can be used only in biopsy needle guide procedures.



Refer to the “Transducers and Consumables” manual for the needle guides to be used with Esaote probes. Needle guide kits which are not tested by

¹ The IKL3323 and IKC3421 are Needle guides with infinite angle capabilities, no needle guideline is available throughout the ultrasound image for these Needle Guides.

WARNING

Esaote could be non compatible with the Esaote probes and thereby could compromise the patient's safety.

2 - Using the Needle Guide

Refer to Getting Started manual for available probes

This chapter explains how to verify that the needle guide is correctly working and how to use it with **MyLab**.

Needle Guide Selection

After pressing the **BIOPSY** key, the **BIOPSY** menu will appear. In this menu the biopsy kit and biopsy angle can be selected. After selecting OK the screen will display two dotted lines circumscribing the working region of the needle guide. The needle insertion angle can be fixed or variable, according to the probe and guide being used.

According to the probe and guide being used, select when necessary, the guide and then the needle insertion angle.

Probe	Insertion angle
SL3332	35°, 25°*, 40°*, 50°*
SL3323 ¹	42°÷81°*
SC3421 ²	42°÷57°*
SC3123	20°÷35°*

* Available only with disposable kit

Checking the Guide's Working

- Assemble the needle guide on the probe, following the instructions provided by the manufacturer or the “Transducers and Consumables” manual.



Carefully read the “Transducers and Consumables” manual: all features relating to system safety, warnings and cautions also apply to needle guides.

^{1 2} The IKL3323 and IKC3421 are Needle guides with infinite angle capabilities, no needle guideline is available throughout the ultrasound image for these Needle Guides.

- Immerse the probe to the allowed limit (see the “Transducers and Consumables” manual) in a water tank.
- Press the **BIOPSY** key and select the guide and the needle insertion angle, then select **OK** to continue.
- **MyLab** will display two dotted lines circumscribing the working region of the needle guide

WARNING

Before proceeding, be sure the kit has been correctly assembled and the needle has been inserted into the guide corresponding to the angle selected during the calibration procedure. Needle insertion into a guide with a different angle may be dangerous to the patient’s safety!

- Both in B-Mode and Freeze, be sure that the needle is correctly displayed within the working area.
- If necessary, evaluate the needle length required by the biopsy procedure.

WARNING

If the needle is displayed outside the working area, do **NOT** use the guide and immediately contact ESAOTE Service.

Using the Needle Guide

Biopsy procedures can be activated on 2D or 2D-CFM image formats.

WARNING

Biopsy procedures must be performed only on real time image. Never move the needle when the image is frozen.

Procedure



- Assemble the needle guide on the probe, following the instructions provided by the manufacturer or the “Transducers and Consumables” manual.

Carefully read the “Transducers and Consumables” manual: all features relating to system safety, warnings and cautions also apply to needle guides.

- Press the **BIOPSY** key and select the guide and the needle insertion angle, then select **OK** to continue.
- **MyLab** will display two dotted lines circumscribing the working region of the needle guide

WARNING

The line displayed on the screen provides an indication of the needle direction, according to the selected guide. Always watch the ultrasound image while inserting the needle into the patient’s body, and be sure that the needle always stays within the displayed area.

Note

At certain scanning depths the needle insertion point or the needle itself could not be displayed.

WARNING

When scanning vascularized structures, display the needle guide working area keeping the CFM mode active so that vessels can be detected and avoided when inserting the needle. Once identified the optimal zone for biopsy, turn CFM off to gain the maximum needle visibility.

Once the needle guide line is displayed, the system temporarily disables all modes, except 2D or 2D-CFM image format. Select **BIOPSY** key to erase the needle guide line and to access other Real Time modes.

**After the Examination**

When the biopsy procedure has been completed, remove the needle and the guide from the probe. Clean the items following the instructions provided in the “Transducers and Consumables” manual and by the manufacturer and, when applicable, dispose of the items according to the local regulations.